

# FOCUSED LIP & EYE AT-HOME HOW-TO

Get personal with those more delicate areas by firming and smoothing the skin around the eyes and lips to create picture-perfect pout and brow-wow.

## SUPPLIES NEEDED:

- + NuFACE Prep-N-Glow® Facial Towelette or Oil-free Cleanser
- + NuFACE® Ionized Super Booster
- + NuFACE® Clean Sweep Applicator Brush
- + NuFACE® Microcurrent Activator
- + NuFACE Trinity® Device
- + NuFACE Trinity® Effective Lip & Eye Attachment

## MICROCURRENT CONTRAINDICATIONS:


- + Epilepsy or subject to seizures
- + Recent surgery or other surgical procedure
- + Pacemaker or electrical implanted device
- + Active cancer
- + Pregnancy
- + Under the age of 18

## STEPS TO FOLLOW:

- 1 CLEANSE**  
Remove makeup, cleanse with an oil-free cleanser or use the **smooth side** of the **NuFACE Facial Towelette**.
- 2 BOOST**  
Apply 3-5 drops for full face coverage of appropriate **NuFACE Booster** onto fingertips and massage into clean dry skin until fully absorbed.
- 3 ACTIVATE 1 minute**  
Apply a **mask-like layer** of appropriate **NuFACE Activator** using **NuFACE Applicator Brush** in sections as you treat.
- 4 LIFT 5 minutes**  
Using the **NuFACE Trinity Device with Effective Lip & Eye Attachment**, perform **glides/holds** around the eyes, lips and nasolabial folds.

### TRINITY EFFECTIVE LIP & EYE ATTACHMENT

**GLIDES** - Glide spheres in the direction of the arrows until the beep sounds. Complete steps 1-2-3 sequence 3 times for each area. **HOLDS** - Lift and hold each position for 1-2 beeps. No repeating.




**HOLD 1** Place both wands horizontally under your lower lip avoid the midline, lift and hold.

**HOLD 2** Place both wands towards the outer corner of your mouth, lift and hold.



**HOLD 1** Place both wands above your upper lip avoid the midline, lift and hold.

**HOLD 2** Place both wands towards the outer corner of your mouth, lift and hold.



**HOLD 1** Place both wands horizontally along nasolabial fold near the outer corner of your mouth, lift and hold.


**HOLD 2** Move up along the nasolabial fold, stopping just left of your nose, lift and hold.



**GLIDE 1** Place both wands horizontally, begin at the outer corner of the eye. Glide along the orbital bone towards the bridge of the nose.

**GLIDE 2** Begin at the outer corner of the eye. Position the wands closer to the under eye. Glide under your eye towards the bridge of the nose.


**GLIDE 3** Begin at the outer corner of the eye. Position the wands closer to the under eye. Glide under your eye towards the bridge of the nose.



**GLIDE 1** Place both wands horizontally, begin at the outer corner of the eye. Glide following the crease of the eyelid towards the bridge of the nose.

**GLIDE 2** Begin at the outer corner of the eye. Glide above the crease of the eyelid towards the bridge of the nose.


**GLIDE 3** Begin at the outer corner of the eye. Glide under the eyebrow towards the bridge of the nose.



**HOLD 1** Place both wands horizontally under the thick of your eyebrow, lift and hold.

**HOLD 2** Place both wands under the arch of your eyebrow, lift and hold.

**HOLD 3** Place both wands under the tail of your eyebrow, lift and hold.



**GLIDE 1** Center the thick of the brow between the two wands and glide towards the center of the forehead. Repeat this glide two more times.

### HOW-TO VIDEO QR CODE

