

ESSENTIAL FIRST STEP:

**MICROCURRENT
ACTIVATOR**

Always prep your skin with NuFACE® Microcurrent Activators before treatment. Our IonPlex™ Ionized Skincare formula is specifically designed to deliver microcurrent to the skin at every level for optimal results.

 **Hydrating Aqua Gel**

Instantly delivers a thirst-quenching burst of hydration and a skin-cushioning bounce.

N Ū F A C E®



Follow @MyNuFACE for Aesthetician Tips, Demos, and Live Q+A.
Share your routine or results with #MyNuFACE

**TRINITY
MINI**

**N Ū F A C E®
QUICK START**

BEFORE YOU POWER UP

CHARGE DEVICE FOR 12 HOURS

Before your first treatment, charge your new NuFACE® Microcurrent Skincare device and read the enclosed User Manual.

SEE HOW EASY TO USE!

Get more out of your new device by visiting useNuFACE.com to watch step-by-step tutorials and learn more.

GETTING STARTED

4 SIMPLE STEPS

1) PREP

Cleanse skin using an oil-free cleanser. This will wash away any excess oils that could otherwise block or prevent an effective microcurrent treatment.

2) ACTIVATE

Apply a mask-like layer of Microcurrent Activator. Use enough formula to ensure the device glides easily along the skin.

3) LIFT

Push the ON/OFF button on your device and wait for the beep. Start on one side of your face, then complete the same sequence on the other side. Keep both spheres in contact with your skin during treatment.

4) FINISH

When complete, massage in remaining Microcurrent Activator for added skincare benefits.

Visit useNuFACE.com to see how easy to use!

5-MINUTE FACIAL-LIFT™



FOREHEAD

Place spheres slightly above brow and glide toward hairline until you hear the beep. **Complete 1-2-3 sequence 3 times.**

CHEEKS

Start by placing spheres on nasolabial fold and glide toward ear until you hear the beep. **Complete 1-2-3 sequence 3 times.**

NECK

Place spheres slightly off-center, avoiding the midline of the neck. Using light pressure, slowly glide toward the back of your neck until you hear the beep. **Complete 1-2-3 sequence 3 times.**

Optimal results are achieved with regular use. For the first 60 days, we recommend using your device at least 5 days a week, even for just 5 minutes a day.