NĪFACE®

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DÉMARRAGE RAPIDE  
INICIO RÁPIDO  
SCHNELLANLEITUNG  
빠른 시작

## ACTIVATE YOUR RESULTS

### ESSENTIAL FIRST STEP: MICROCURRENT ACTIVATOR

Always prep your skin with NuFACE® Microcurrent Activators before treatment. Our IonPlex™ Ionized Skincare formula is specifically designed to deliver microcurrent to the skin at every level for optimal results.

 **Hydrating Aqua Gel**  
Instantly delivers a thirst-quenching burst of hydration and a skin-cushioning bounce.

## 5-MINUTE FACIAL-LIFT®

### FOREHEAD

Place spheres above brow and glide up until you hear the beep. **Complete 1-2-3 sequence 3 times.**

### CHEEKS

Place spheres on nasolabial fold and glide toward the ear until you hear the beep. **Complete 1-2-3 sequence 3 times.**

### NECK

Avoid midline of the neck. Slowly glide toward the back of your neck until you hear the beep. **Complete 1-2-3 sequence 3 times.**

## BEFORE YOU POWER UP

### CHARGE DEVICE FOR 8 HOURS

Before your first treatment, charge your new NuFACE Microcurrent Skincare device and read the enclosed User Manual.

## GETTING STARTED

### 4 SIMPLE STEPS

#### 1] PREP

Cleanse skin using an oil-free cleanser. This will wash away any excess oils that could otherwise block or prevent an effective microcurrent treatment.

#### 2] ACTIVATE

Apply a mask-like layer of Microcurrent Activator. Use enough formula to ensure the device glides easily along the skin.

#### 3] LIFT

Push the ON/OFF button on your device and wait for the beep. Start on one side of your face, then complete the same sequence on the other side. Keep both spheres in contact with your skin during treatment.

#### 4] FINISH

When complete, massage in remaining Microcurrent Activator for added skincare benefits.

*Optimal results are achieved with regular use.*

For the first 60days, we recommend using your device at least 5 days a week, even for just 5 minutes a day.

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