


vive®

WIRELESS FLOOR SENSOR MAT

User Manual

LVA1074FLR



vivehealth.com

USER GUIDE

The Vive floor sensor mat is designed to alert caregivers of falls or unwanted movement by sounding an alert whenever pressure is applied on sensor mat. The sensor mat can be placed on the floor to provide security and peace of mind that no movement goes unnoticed.

SETTING UP THE SENSOR MAT

Step 1:

Write the start date and expiration date on the mat with a permanent marker. The usable life of this pad is 365 days.


Step 2:

Step 2: Connect the transmitter to the sensor mat.

Step 3:


Place the sensor mat on the desired location. Make sure the transmitter is mounted where the caregiver can access it conveniently.

Step 4: Proceed to programming your transmitter to your wireless pager if not already complete.




PROGRAMMING THE TRANSMITTER TO A WIRELESS PAGER

Step 1: Open the battery door of the pager and insert two AA batteries as shown. Close the battery door.




Step 2: On the pager, press and hold both UP and DOWN arrows for 5 seconds. You will see the number "1" is blinking in red, indicating that the pager is now in programming mode.



Step 3: There are 10 channels on the pager for programming up to 10 individual wireless devices. Select a number for programming by pressing the UP or DOWN arrow on the front of the pager.

Step 4: While the pager is in pairing mode, a red light should flash over the device number (#1-10) that you are attempting to pair the alarm with. While this light is flashing, press the button on the transmitter once to wake it up, then hold for 5 seconds until the alarm light on the transmitter flashes orange. Then release the button on the transmitter.



Step 5: The pager will beep to signal a successful connection. Then it will re-enter programming mode in case you need to connect other wireless devices. If you do not need to connect more devices press and hold both the UP and DOWN arrows on the front of the pager to exit programming mode.




Step 6: Test the connection by stepping on the sensor mat. The alarm should sound on the pager and the connected channel number should flash in red while the alarm sounds. To silence the alarm press the STOP button on the pager. If this does not occur follow these steps again from the beginning or contact us for help.




Note: The wireless pager sounds the alerts given from the sensor mat. The pager is sold separately.


NOTIFICATION MODES FOR THE WIRELESS PAGER



Step 1: Set notification mode for the WPG.

The WPG has three notification modes.

Notification Mode	Notification	Sign
Audible notification	Beeps	
Vibration notification	Vibration	
Audible + Vibration notification	Audible and Vibration	

To set the notification modes, press and hold down both   Button on the front of the WPG for more than 5 seconds until the  image is blinking.



To select another notification mode, press and release the UP Button or DOWN Button on the front of the WPG. The next notification light  will start blinking. Continue to press and release the Up or Down buttons until the desired notification light is selected. Once you've landed on your desired notification light, press and release the  Button, you will hear a 1 second beep and/ or feel a 1 second vibration. This means the current notification mode is recorded to the WPG.

Step 2: Exit programming mode: Press and hold both the UP and Down buttons on the WPG for more than three seconds until you see the number light turn off. The WPG has now exited programming mode.


ADJUSTING THE TRANSMITTER TIME DELAY

Press and release the time delay button on the front face of the transmitter until desired time delay setting is selected. Options are 0 seconds, 2 seconds and 4 seconds.

Note: This function can reduce false alarms if the patient momentarily steps on the mat within the set delay time.


TEST THE FLOOR ALARM TRANSMITTER

It is important to test the floor alarm before using with a patient. To test, apply pressure on the floor mat and release. If working correctly, the alarm light will flash one time and the programmed receiving devices will alarm. Press the button on the front of each receiving device to silence alarm.



CHANGING THE TRANSMITTER

The usable life of this floor alarm transmitter is 24 months. Actual life depends on how frequently the transmitter is used. When the transmitter needs replacing, the red "Change Alarm" light will flash. To replace the transmitter, remove the transmitter from the pad and connect a new transmitter to the pad.



WARNINGS

This device will not prevent a fall or injury. It is a supplemental part of an overall fall management program.

Failure to follow all product instructions may result in the failure or improper function of the product.

The WPG is able to communicate up to 150 feet away from the patient in open field. The actual communication distance or wireless range will vary from building to building. Test the wireless range in your building prior to putting the system into service. Always keep the pager located within the wireless range.

The mat, transmitter and any other devices should be checked prior to each use for proper functioning.

TROUBLESHOOTING

When your system is not properly functioning, please follow these steps.

Step 1. Change the batteries if the “Change Alarm” light is lighted

Step 2: Change the transmitter if the “Change Alarm” light is lighted.

Step 3. Check that the mat is properly plugged into transmitter

Step 4: Check on how long the WPG and the mat have been in use. They may have run out of usual life.

Step 5: To determine which piece of equipment is at fault, follow the steps below:

- a. Take a new working transmitter and connect with the original mat.
- b. If the system works properly, the problem is with the original transmitter.
- c. If the system still does not work properly, the problem is with the original mat.

ACTIVATE YOUR GUARANTEE



You are protected by Vive Health's industry leading guarantees and customer service:

If you did not purchase through vivehealth.com, please register at vivehealth.com/register to validate your guarantee.

Product Code: LVA1074FLR

✉ service@vivehealth.com

☎ 1-800-487-3808
Monday - Friday 9am - 5pm EST