

# 尺寸：128\*168mm



#### [Product Details]

Product name: LED Photon Facial Mask  
Product model: TM-MR-115  
Voltage: 5V  
Red light wavelength: 630nm ± 30nm  
Blue light wavelength: 430nm ± 30nm  
Near infrared wavelength: 830nm ± 30nm  
Orange light wavelength: 600nm ± 30nm  
Size: 415 \* 205 \* 4mm  
Product weight: about 325g  
Main materials: electronic biological light source, medical grade silicone

#### [Product introduction]

Use high-tech to control the red, blue, orange and near-infrared light waves of specific wavelengths, stimulate the skin in pulse mode, and control the pulse intensity by the respiratory frequency of the skin. Realize various cosmetic functions of LED biological active light:

1. Red light: stimulate collagen regeneration and tighten skin
2. Blue light: inhibit acne and repair acne damage;
3. Orange light: whiten skin, lighten spots and improve fine lines;
4. Near infrared light: soothe the skin, expel toxins, inhibit and Anti wrinkles;
5. The light lasts for 10 minutes and is automatically turned off at regular intervals.

The beauty effect of promoting skin collagen regeneration, Make skin become more brighting and whitening. Promote and reduce fine lines and wrinkles. Shrink pores, improve skin brightness, make skin smooth and elastic. Dilute stains & red blood, accelerate skin metabolism and repair after sunburn. Promote facial blood circulation and lock skin moisture. to achieve delicate skinmake the face cells activate energy to make the skin more tender and firmingImprove acne and detumescence, and improve all skin problems. make your skin texture better and better

#### [Instructions]

1. First wash your face, then fix the mask on your face and fasten belt. please be careful to wear the mask gently and buckle belt on the ear to avoid falling
  2. The first charge should last for 12 hours and last for 10 minutes each time use.
  3. Please close your eyes when using the mask
  4. If want to turn off in use; Please press the button for 3 seconds;
  5. It is recommended to irradiate at least once a day or more according to the skin absorption
  6. It will have significant effect if it is used continu for 7-14 days
1. Clean and keep dry the face
  2. Use LED photon mask directly or use with traditional facial mask
  3. Lie down in a comfortable position, clip the photon mask belt on both sides of the ear, and adjust comfortable position, close your eyes, turn on switch, and use last for 10 minutes
  4. Remove the photon mask after10 minutes
  5. Gently massage the face for 5-15 minutes

#### [Product List]

Skin rejuvenating mask X1  
User Manual X1  
Controller X1  
Strap X1  
Charging cable X1  
Warranty card x1

#### [Dear users]

If this product does not work properly, please fill in this form carefully. Be sure to contact us and send it to the company, we will archive it so that we can provide more considerate services.

1. From the date of purchase, in case of failure in normal use, free of charge for six months, our company does not provide freight
2. The following conditions are excluded in the warranty scope.
  - 1) No warranty is available.
  - 2) Damage caused by user's operation error.
  - 3) Products are transformed without permission.
  - 4) Spare parts consumption.
  - 5) Change of product appearance.

Product	TM-MR-113	
Name		Telephone
Address		
Purchased		Purchase date
The basis for purchasing this product is: <input type="checkbox"/> Fashion <input type="checkbox"/> Advertisement <input type="checkbox"/> User recommendation <input type="checkbox"/> Online purchase <input type="checkbox"/> Recommended by friends		
User survey:		

Note: This product has been carefully inspected and quality controlled.  
The warranty only applies to this product

#### [Precautions]

1. When using the phototherapy mask, do not charge it immediately when the light is changing dark to avoid damaging the battery;
2. Please do not put the phototherapy mask and charger in a damp place;
3. Please connect all connections correctly before turning on the power;
4. Use a clean towel to wipe the facial mask, and do not directly contact the water of the mask body;
5. The circuit stripes in the mask cannot be squeezed or folded. Please use a suitable flat plate, do not put heavy objects on it;

#### [Warranty matters]

1. From the date of purchase, the entire product is guaranteed for six months free of charge;
2. Product damage not caused by quality problems shall be repaired by the customer's expense;
3. The nature of spare parts is not within the scope of warranty;
4. If the warranty card is lost or changed, the warranty will not be free

内1

内2