tria.

HAIR REMOVAL LASER 4X

Instructions for Use

EN

SKIN TONE CHART

CHECK YOUR SKIN TONE to find out if you can safely use the Tria Hair Removal Laser 4X. Match the skin in the area you wish to treat to the range of skin tones shown on this chart. If your skin tone matches any of the skin tones labelled "Not Safe" for treatment, then **DO NOT USE** the Tria Laser 4X.

INTENDED SKIN TONES:

NOT SAFE:

FOR USE ON:

Legs; Arms; Bikini Line; Underarms; Feet; Hands; Chest; Back; Stomach; Women's Facial Hair (cheek line downward)

READ THESE INSTRUCTIONS BEFORE USE

Before you use your new Tria Laser 4X, read this booklet for important information:

- Confirm that the Tria Laser 4X is safe for you.
- Learn how to use the Tria Laser 4X.
- View our treatment results timeline.
- Read our tips for best results.
- Review contraindications and warnings, which are marked with the symbol below:

 This **Warning Symbol** appears next to information about possible safety risks. Be sure to read and follow all the warnings and safety information.

We recommend that you save these instructions for future reference. If you have questions or concerns, please contact Tria Beauty Customer Service. You can also find the *Quick Start Guide* on our website (see back cover).

TABLE OF CONTENTS

INTRODUCTION	2
GETTING TO KNOW YOUR TRIA LASER 4X	3
WHAT THE TRIA LASER 4X IS USED FOR (INDICATIONS)	5
WHO SHOULD NOT USE THE TRIA LASER 4X (CONTRAINDICATIONS)	5
WARNINGS	5
BENEFITS AND RISKS	7
TRIA LASER 4X TREATMENT OVERVIEW	8
WHAT TO EXPECT	9
BEFORE YOUR FIRST TREATMENT	9
BEFORE EVERY TREATMENT	10
HOW TO TREAT	12
TAKING CARE OF YOUR TRIA LASER 4X	14
TROUBLESHOOTING	14

CONGRATULATIONS ON YOUR NEW TRIA LASER 4X

The Tria Hair Removal Laser 4X gives you permanent results at home, unlike temporary solutions such as shaving and waxing. The Tria Laser 4X works in the same way as professional laser hair removal. It treats the hair follicle, not just the hair, giving you permanent hair reduction after a series of treatments.

LET'S GET STARTED!


INTRODUCTION


HOW DOES LASER HAIR REMOVAL WORK?

The Tria Laser 4X deactivates the follicles in your skin that produce hair. It emits a pulse of laser light, which heats up the dark pigment inside the hair and deactivates the hair follicle. That's why it can deliver permanent results and works only on naturally brown or black hair.

To be most effective, the Tria Laser 4X needs to treat each hair follicle during its active growth period. Hair grows in cycles, so multiple treatments are necessary. For best results, treat every two weeks for up to 3 months or until you are satisfied with the amount of hair reduction.

GETTING TO KNOW YOUR TRIA LASER 4X

WHAT THE TRIA LASER 4X IS USED FOR (INDICATIONS)

The Tria Hair Removal Laser 4X is used to remove unwanted hair. It is an over-the-counter system intended for adjunctive use with shaving for hair removal sustained with periodic treatments. It is also intended for permanent reduction in hair regrowth defined as long-term, stable reduction in the number of hairs regrowing when measured out to 6, 9 and 12 months after the completion of the treatment regime.

It is intended for use only on light to medium skin tones (see the Skin Tone Chart on the inside front cover) and naturally brown or black hair. It will not work on white, grey, blonde or red hair.

WHO SHOULD NOT USE THE TRIA LASER 4X (CONTRAINDICATIONS)

DO NOT use the Tria Laser 4X if your skin tone matches one listed as “Not Safe” on the Skin Tone Chart on the inside front cover of this booklet. If you use the Tria Laser 4X on skin deemed “Not Safe” for treatment, the skin may absorb too much laser light and be injured (burned, discoloured or scarred).

PLEASE NOTE: The skin sensor is a safety feature that unlocks the Tria Laser 4X. If it is not safe for you to use the Tria Laser 4X, you won't be able to unlock it with the skin sensor.

DO NOT use the Tria Laser 4X if you are under the age of 18 or use it on anyone else under the age of 18. The use of the Tria Laser 4X on or by children has not been studied and could result in serious injury.

DO NOT use the Tria Laser 4X if you have a lighter hair colour, such as blonde, white, red or grey. In order for the laser to work, you must have naturally dark hair where you are treating. Lighter hair pigments will not absorb enough of the laser's energy to disable the hair follicle.

WARNINGS

⚠ DO NOT use the Tria Laser 4X on or around the eyes, eyebrows or eyelashes. Doing so can cause permanent eye injury.

⚠ DO NOT use the Tria Laser 4X if the Skin Tone Chart shows that your skin tone is not suitable. If you have African, Indian, Native American or Pacific Island heritage, you may have a skin tone that is not suitable—be sure to check the Skin Tone Chart. Darker skin tones can absorb too much laser light, which can injure the skin.

 **DO NOT use the Tria Laser 4X** anywhere on the ears, nipples, genitals or around the anus. These areas may have more sensitive skin, a darker skin tone, and/or a greater hair density, and using the Tria Laser 4X there may cause serious skin injury.

 **DO NOT use the Tria Laser 4X** on parts of your skin that are darker than your usual skin tone: tattoos, dark brown or black spots (such as large freckles, birthmarks or moles) or skin that is red or swollen from a previous treatment. Dark skin may absorb too much laser light, which can injure the skin.

 **DO NOT use the Tria Laser 4X** on a man's face, jaw or neck. The hair in a man's beard area is too dense, and using the Tria Laser 4X in these areas may result in skin injury.

Hair removal by laser can cause increased hair growth in some individuals. Based upon currently available data, the highest risk groups for this response are females of Mediterranean, Middle Eastern and South Asian heritage treated on the face and neck.

 **DO NOT use the Tria Laser 4X** on damaged skin, including open wounds, cold sores, sunburns or areas where you have had a skin peel or other skin resurfacing procedure within the last 6-8 weeks (check with your doctor). These conditions can make the skin more susceptible to injury.

 **DO NOT use the Tria Laser 4X** if you are pregnant or breast-feeding. The Tria Laser 4X has not been tested on pregnant women. Hormonal changes could increase skin sensitivity and the risk of skin injury.

 **DO NOT use the Tria Laser 4X** if it appears damaged, and **DO NOT** attempt to open or repair the Tria Laser 4X. These actions increase the risk of skin or eye injury. Please contact Customer Service for assistance.

 Always follow the basic safety precautions for electrical devices: **DO NOT** get the Tria Laser 4X or the charger damp or wet. **DO NOT** use the Tria Laser 4X while bathing or put it next to a bath or sink. **DO NOT** operate a charger with a damaged plug or cord. **ALWAYS** keep the cord away from heated surfaces. Failure to follow these precautions can cause electrocution or electric shock.

 **DO NOT use the Tria Laser 4X** for any purpose other than hair removal.

 The Tria Laser 4X is for personal use. The Tria Laser 4X is not for use in a medical office, clinic, salon or other commercial setting.

 **DO NOT use the Tria Laser 4X** on any area where you may someday want your hair to grow back. The results may be permanent and irreversible.

 Invisible Laser Radiation. Do Not View Directly with Optical Instruments. Class 1 Laser Product.

 Viewing the laser output with certain optical instruments (for example eye loupes, magnifiers, and microscopes) within a distance of 100 mm may pose an eye hazard.

 No modification of this equipment is allowed.

 This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.

BENEFITS AND RISKS

A medically supervised clinical study has shown that the Tria Laser 4X is safe and works when used as directed.

Seventy-nine (79) users performed a 6-week course of treatment (3 treatments) and all side effects were observed. The number of hairs at each treatment site was counted before each treatment and at 1, 2, 3, 6 and 12 months after the last treatment. No additional treatments were performed after the third treatment.

BENEFITS

When used as directed, the Tria Laser 4X gives you permanent results. The results vary from person-to-person and depend on which body area is treated and how closely treatment instructions are followed. The clinical study demonstrated the following benefits:

- 61% average reduction in hair count 3 weeks after the first treatment
- 70% average reduction in hair count 3 weeks after the second treatment
- 60% average reduction in hair count 1 month after the third treatment
- 33% average reduction in hair count 12 months following the third treatment

In the study, 69 of 79 users (87%) reported finer hair at 1 month after the last treatment than before any treatment. And 55 of 79 users (70%) reported **lighter** hair at 1 month after the last treatment than before any treatment.

These results were obtained with only 3 treatments. **For best results, treat every 2 weeks for 3 months** or until you are satisfied with the amount of hair reduction.

RISKS

 It is important to use the Tria Laser 4X properly. Make sure that you understand and follow all warnings (see pages 5-7). The skin sensor will prevent use of the device on contraindicated skin types (see page 5).

If you use the Tria Laser 4X as directed, the risks are low:

- In the clinical study, almost all users reported mild to moderate pain or discomfort during their first treatment, with a reduction in sensitivity reported in their later treatments. This is normal. You may feel warmth, burning, tingling or itching; some users described feeling the same level of discomfort as when a rubber band snaps against the skin. If the pain is too intense or persists after a treatment, stop using the Tria Laser 4X and check with your doctor before using it again.
- Your skin may become red right after treatment or within 24 hours. This redness generally clears up within a few hours. See your doctor if the redness does not go away within 2-3 days.
- In rare cases, some users may see an increase in the number of hairs in the treated area after the initial treatment. If this increase persists over several treatments, stop using the Tria Laser 4X and see your doctor.

The following symptoms affected no users in the clinical study (0 out of 88) but have been reported with other hair removal lasers:

- Skin swelling, especially around the hair follicle. If the swelling doesn't go away within 2-3 days, stop using the Tria Laser 4X and see your doctor.
- Temporary or permanent skin discolouration. If your skin colour changes, stop using the Tria Laser 4X and see your doctor.

- Blisters. This is very unlikely for people with a light or medium skin tone. If you do develop a blister, stop using the Tria Laser 4X and don't let the blister become infected. Consult your doctor if needed.

TRIA LASER TREATMENT OVERVIEW

BEFORE YOUR FIRST TREATMENT

- 6 weeks before your treatment, stop waxing or plucking hair and do not get a skin peel.
- At least 1 day before treatment, perform a patch test on your skin (see page 10).

BEFORE EVERY TREATMENT

- Charge your Tria Laser 4X (see page 9).
- Unplug it from the battery charger.
- Cleanse, shave and dry your skin (see page 10).
- Turn on the Tria Laser 4X (see page 11).
- Unlock it with the skin sensor (see page 11).
- Set the treatment level (see page 11).

BETWEEN TREATMENTS

- Shave, but don't wax or pluck hair.

TREAT UNTIL YOU'RE DONE

- Be consistent and follow treatment instructions.
- For best results, treat every two weeks for up to 3 months or until you see the desired results.

WHAT TO EXPECT

WHILE USING THE TRIA LASER 4X

- The sensation varies, depending on the person, treatment area and treatment level. Some users may feel the sensation of a rubber band snapping on the skin. Others will feel a warm, prickling sensation.
- As the amount of hair in the treated areas decreases, the sensation will also decrease. For this reason, you may be able to use a higher treatment level after your initial treatments.
- Some users experience mild redness that goes away within 24 hours.

DURING YOUR FIRST 3 MONTHS OF TREATMENT

- After your first treatment, there will be no visible change to the hair. Laser hair removal does not cause hair to immediately fall out.
- With additional treatments, there will be a gradual reduction in the number of hairs with each treatment.
- The hair that does regrow will also begin to lighten and become finer.


AFTER 3 MONTHS OF TREATMENT

- Hair follicles that have been deactivated by the Tria Laser 4X will not grow back.
- Do occasional touch-ups or additional treatments if necessary.

BEFORE YOUR FIRST TREATMENT

CHARGE THE TRIA LASER 4X

Charge the Tria Laser 4X for at least 2 hours before you use it or anytime the battery indicator is "empty" and blinking. If you are covering a large area or have the Tria Laser 4X set to a high energy level, you may need to charge the Tria Laser 4X in the middle of a treatment. This is normal.

TO CHARGE:

1. Plug the battery charger into the Tria Laser 4X charger socket and into a power outlet.
2. The battery indicator will show it is "filling up" while the Tria Laser 4X is charging.
3. Once fully charged, the battery light will stay on and you may unplug the charger. The Tria Laser 4X can remain plugged in for extended periods of time without damage.

NOTE: To prevent damage to your Tria Laser 4X, use only the included Tria battery charger. If you lose the charger, please contact Customer Service.

NOTE: You may charge the laser in Quiet mode (charges without the sound of cooling fan), but it will take longer. To change charging modes, press **and hold** the power button while charging **until you hear a tone**.

TEST THE LASER ON A PATCH OF SKIN

Before your first full treatment on a new area of your body, we recommend you test the Tria Laser 4X on a patch of skin there and try different energy levels.

1. Prepare your skin and unlock the Tria Laser 4X (see page 10). The Tria Laser 4X will be set to Level 1 (low).
2. Place the treatment window flat against your skin, and test one laser pulse.
3. If Level 1 was tolerable for you, then move to another test area and increase to Level 2, testing a patch at that level with one laser pulse. To adjust the Tria Laser 4X setting up or down, press the power button (see page 11).
4. If Level 2 was tolerable for you, change test locations and increase by one level until you find

the highest setting that is still tolerable. Treating at the highest level possible will maximise results.

5. Wait 24 hours and then examine the test areas. If the skin appears normal, start your first treatment at the highest level you are comfortable using.

Retest your skin if your skin tone changes dramatically (due to tanning).

 If you notice blistering, burns, changes in your skin tone or any other side effects, or if you are not sure about the outcome of the test, **DO NOT USE** your Tria Laser 4X. Contact Customer Service or your doctor.

BEFORE EVERY TREATMENT


PREPARE YOUR SKIN

1. Cleanse the treatment area, removing all cosmetics, lotions and creams that may block the laser light.
2. Shave the treatment area. This removes the surface hair, which improves comfort and allows the light to better reach the hair in the follicle.
3. Dry the treatment area.

UNLOCK THE LASER

The Tria Laser 4X will lock when it's not in use. To unlock the Tria Laser 4X, use the skin sensor, which helps confirm that the Tria laser is safe for your skin tone. To unlock:

1. Unplug the Tria Laser 4X from the battery charger and press the power button to turn it on. The display will show the lock and skin sensor icons.
2. Hold the skin sensor (red light at the base of the Tria Laser 4X) against the skin in the area you wish to treat. The skin sensor will automatically sense and check your skin. If the Tria Laser 4X is safe for your skin tone, the Tria Laser 4X will emit a series of beeps and unlock.


3. If the lock icon is still visible and you do not hear a series of beeps, then it may not be safe for you to use the Tria Laser 4X. Or, you may be using the skin sensor incorrectly and need to try again. To try again, lift the skin sensor from your skin and follow Steps 1 and 2 above. If you still cannot unlock the Tria Laser 4X, please contact Customer Service for assistance.

SET THE TREATMENT LEVEL

The Tria Laser 4X has different energy levels. In a clinical study, higher levels were found to be 2-3 times more effective. **So for best results, we recommend you use the highest setting that you can tolerate.**

When the Tria Laser 4X is on and unlocked, press the power button to select your desired level. Each time you press the power button, the level will increase. Once the Tria Laser 4X is at the maximum level, pressing the power button again will return it to the lowest level.


HOW TO TREAT

1. TURN ON THE LASER, UNLOCK IT AND SET THE TREATMENT LEVEL (SEE PAGE 11).

2. TOUCH AND HOLD THE TREATMENT WINDOW TO THE SKIN UNTIL YOU HEAR TWO BEEPS.

- If you do not hear beeps, make sure that the treatment window is flat against your skin. The laser will only pulse when the whole window is touching your skin.
- If you hear a buzzing tone instead of the second beep, you lifted the device from the skin before the pulse was over.
- Do not move the treatment window between the first and second beep because this will disrupt the laser pulse and reduce its effectiveness.


3. LIFT THE LASER AND PLACE IT BACK ON THE SKIN SO THAT IT OVERLAPS THE PREVIOUS SPOT.

- Each time you reposition the treatment window, move it about half the distance of the window (about 1/4 inch or 5 mm). This creates an overlapping pattern in the treatment area, which ensures you cover all the hair follicles. This works out to about 25 laser pulses per square inch or 40 pulses per 10 square cm.

⚠ DO NOT treat repeatedly in the same spot.

Always move the treatment window to a new area after each laser pulse. Failure to do so may result in discomfort or skin heating, which could cause skin injury.

4. KEEP GOING UNTIL YOU HAVE COMPLETELY COVERED THE AREA YOU WANT TO TREAT.

Use the Tria Laser 4X Pulse Counter to help get the best results. See the table below as a guide to how many pulses you may need to cover each area per treatment. These pulse counts and times are approximations – your actual treatment area and number of pulses may vary. Always follow

the treatment instructions and warnings provided in this booklet.

NOTE: When treating larger areas like legs and back, you may need more than 600 pulses to cover the area using an overlapping pattern. If so, please recharge the Tria Laser 4X and continue treating until you've covered the desired treatment area.

PULSE COUNTER GUIDE

BODY PART	APPROXIMATE # OF PULSES PER TREATMENT	APPROXIMATE MINUTES
Underarms	100 per side	4-5 per side
Bikini	200 per side	8-10 per side
Upper Lip	25	2-3
Upper Leg	600 per side	25-30 per side
Lower Leg	600 per side	25-30 per side

Pulse counts and times are approximations.

5. TURN OFF THE LASER BY HOLDING DOWN THE POWER BUTTON.

Repeat treatment every 2 weeks for up to 3 months or until you reach the desired results.

NOTE: The cooling fan may go on and off during your treatment. This is normal and nothing to worry about.

TAKING CARE OF YOUR TRIA LASER 4X

The Tria Laser 4X doesn't require special care, but we do recommend that you keep it clean, especially the treatment window and skin sensor.

CLEANING INSTRUCTIONS

- Make sure that the Tria Laser 4X is turned off and unplug the battery charger from the Tria Laser 4X and the power outlet.
- To clean the Tria Laser 4X and battery charger, moisten a soft cloth with water and mild soap, and wipe the treatment window and the exterior of the device. Be careful not to damage the treatment window. Allow to dry completely before charging or turning on your Tria Laser 4X.
- To clean the skin sensor located at the base of the Tria Laser 4X, wipe it gently with a dry, lint-free cloth.

TROUBLESHOOTING

Please contact Customer Service for additional help.

YOUR QUESTION	IT'S POSSIBLE THAT	HOW TO TROUBLESHOOT
My Tria Laser 4X won't turn on.	Your Tria Laser 4X is not charged.	Charge your Tria Laser 4X. (see page 9)

 **DO NOT** get the Tria Laser 4X or battery charger damp or wet. If it gets wet, **DO NOT** use the Tria Laser 4X or the battery charger. Contact Customer Service. **Performing a treatment with a wet Tria Laser 4X (or component) can cause electrocution or electric shock.**

LONG-TERM STORAGE

If you do not plan on using your device for a long period of time, we recommend that you recharge the device every 3 months to maintain the life of the battery.

RECYCLING/DISPOSAL

The Tria Laser 4X contains lithium batteries that cannot be replaced. To recycle or dispose, please follow the regulations where you live or contact Customer Service for information.

YOUR QUESTION	IT'S POSSIBLE THAT	HOW TO TROUBLESHOOT
While I'm charging my Tria Laser 4X, the battery indicator on the display doesn't "fill up."	Your Tria Laser 4X is not connected properly.	Confirm that your Tria Laser 4X is firmly plugged into the charger and the wall outlet.
	Your Tria Laser 4X is already fully charged and ready to use.	Unplug the battery charger from the wall and from the Tria Laser 4X.
Even after I charge my Tria Laser 4X, the battery indicator is still "filling up."	Your Tria Laser 4X isn't fully charged.	Charge your Tria Laser 4X. (see page 9)
When I put my Tria Laser 4X against my skin, nothing happens.	Your Tria Laser 4X is locked.	The Tria Laser 4X will lock when it's not in use. Unlock your Tria Laser 4X with the skin sensor. (see page 11)
When I put my Tria Laser 4X against my skin, nothing happens.	Your Tria Laser 4X is plugged into the battery charger.	Your Tria Laser 4X will not work when it is plugged in to the battery charger. Unplug your Tria Laser 4X.
	The treatment window isn't completely in contact with your skin.	Confirm that the treatment window is completely flat against your skin.

YOUR QUESTION	IT'S POSSIBLE THAT	HOW TO TROUBLESHOOT
When I put my Laser 4X against my skin, nothing happens.	The area you're trying to treat is too small, curved or bony to treat.	Try again, pulling the skin flat. If you can't get the whole treatment window to touch the treatment area, then the Tria Laser 4X won't work.
I can't use my Tria Laser 4X and an image of a lock shows on the display.	Your Tria Laser 4X is locked.	The Tria Laser 4X will lock when it's not in use. Unlock your Tria Laser 4X with the skin sensor. (see page 11)
The display is flashing quickly.	Your Tria Laser 4X is not working properly.	Please contact Customer Service for assistance.
Dirt or debris is stuck to the treatment window.	Your Tria Laser 4X needs to be cleaned.	See cleaning instructions on page 14. If you can't get the area clean, contact Customer Service for assistance.
I tried to unlock the Tria Laser 4X, but it's still locked.	It may not be safe for you to use the Tria Laser 4X.	Match your skin tone to the Skin Tone Chart on the inside cover of this booklet, and contact Customer Service with any questions.
	You didn't use the skin sensor correctly.	Try again, holding the skin sensor flat and firm against your skin. (see page 11)

YOUR QUESTION	IT'S POSSIBLE THAT	HOW TO TROUBLESHOOT
When I reposition the Tria Laser 4X, it seems to take a long time before the next pulse of laser light.	The room is warm. The Tria Laser 4X may pulse at a slower rate in a warm environment.	Use the Tria Laser 4X in a cooler environment.
My Tria Laser 4X was working properly, but then it stopped working. The display is off and does not turn on even when I press the power button.	The battery needs to be charged.	Charge your Tria Laser 4X. (see page 9)
When I'm using my Tria Laser 4X, I can hear a fan working.	This is normal.	Continue your treatment.

TECHNICAL SPECIFICATIONS

LASER TYPE: Diode laser

LASER CLASS: 1
According with IEC 60825-1:2014

OUTPUT WAVELENGTH: 810 nm

OUTPUT FLUENCE: 7-20 J/cm²

BATTERY CHARGER INPUT:
100-240 V - , 50/60Hz, 29 VA

BATTERY CHARGER OUTPUT:
5.0 V DC, 2.0 A

EXPECTED SERVICE LIFE: 3 Years minimum

PHYSICAL

PHYSICAL OPERATING TEMPERATURE:

5°C to 30°C

OPERATING RELATIVE HUMIDITY:

Non-condensing 15-85%

OPERATING ATMOSPHERIC PRESSURE:

700 hPa to 1060 hPa

TRANSPORT AND

STORAGE TEMPERATURES:

-20°C to 50°C

TRANSPORT RELATIVE HUMIDITY:

Non-condensing 0-93%

SYMBOL GLOSSARY

Warning Symbol

Laser Warning Symbol

WEEE Symbol

Type BF Symbol

Serial Number Symbol

CE Mark

Date of Manufacture Symbol

Manufacturer Symbol

Temperature Limitation Symbol

Accompanying Documentation
Must be Consulted Symbol

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Protected by patent(s). See tribeauty.com/legal. All rights reserved. Tria, Tria Beauty, and the Tria logo are trademarks of Tria Beauty. / Protegido por patente(s). Ver tribeauty.com/legal. Todos los derechos reservados. Tria, Tria Beauty y el logotipo de Tria son marcas registradas de Tria Beauty. / Protégé par un ou plusieurs brevets. Voir tribeauty.com/legal. Tous droits réservés. Tria, Tria Beauty et le logo Tria sont des marques de commerce de Tria Beauty. / Protetto da uno o più brevetti. Vedi tribeauty.com/legal. Tutti i diritti riservati. Tria, Tria Beauty, e il logo Tria sono marchi di fabbrica di Tria Beauty. / Geschützt durch Patent(e). Siehe tribeauty.com/legal. Alle Rechte vorbehalten. Tria, Tria Beauty und das Tria-Logo sind Marken von Tria Beauty. / Chronione patentem/-ami. Patrz: tribeauty.com/legal. Wszystkie prawa zastrzeżone. Tria, Tria Beauty i logo Tria są znakami towarowymi Tria Beauty. / Beschermd door octrooi(en). Zie tribeauty.com/legal. Alle rechten voorbehouden. Tria, Tria Beauty, en het Tria-logo zijn handelsmerken van Tria Beauty.
EC REP MDSS GmbH, Schiffgraben 41
Tria Beauty, a company of Channel Investments, LLC