

www.CandyLipz.com

www.CandyLipz.com

CONTENTS

Thank You Note		1
About CandyLipz® Technology		2
CandyLipz® AMAZING Styles		2
Caution Before Use		4
Who Should Not Use CandyLipz®		5
Lip Plumper Styles and Sizes		6
Xtreme Lip-Shaper® Set BENEFITS		8
What Is In Each Box For Model A, B, C, or D?		9
About CandyLipz® Premium Pump Materials		9
Product Accessories		9
How To Use The Blocker J. J. in		10
How To Use Mouth Cover		11
Safety Precatuions		12
Warnings		12
Care & Storage		13
What to do After You Get Your Lip Pump		14
Instructions for Use		14
Frequently Asked Questions		20
Beauty Tips		30
Before and After Gallery		32
Return Policy		34
Contact		35
Disclaimers		36
www.CandyLi	$\overline{nz.c}$	om

vww.CandyLıpz.com



Clinically proven to last for up to 2 hours!

THANK YOU FOR YOUR PURCHASE!

You are now a proud owner of the world's most OBSESSED over lip plumping beauty tool. Get ready to experience the thrill of your life! Go from ordinary to magazine perfect in minutes. You decide how huge you want your lips to be! The only lip plumping system that knows no limit - "The Rolls Royce of Lip Plumpers." ~ Splash Magazine, Hollywood.

ABOUT CandyLipz® Technology

CandyLipz® Xtreme Lip-Shaper® System is made possible by its revolutionary breakthrough technology for enhancing the appearance of lips using dual elements combining hand-operated suction with an advanced built-in, patent-pending Lip-Shaping technology. CandyLipz won ten prestigious awards internationally and is the first lip enhancement, home-use cosmetic tool clinically proven to increase lip volume and safe for intended use.

CandyLipz[®] AMAZING Styles

There isn't anything like CandyLipz[®] in the world. You can now shape and contour the appearance of your lips the way you wanted. CandyLipz gives you the option to enhance both the upper and lower lip together or separately. It is the only system on the market to create the luscious Single-Lobed or Double-Lobed lip styles.

GET BREATHTAKING LIPS INSTANTLY



CandyLipz is portable. It can be used anytime and anywhere. It is designed to be used hands-free so you can still work with your hands. Your lip enhancement session takes just two minutes and the look lasts for up to two hours- clinically proven! You can look magazine perfect without the side effects of plastic surgery, injections, drugs, or chemicals.

CAUTION BEFORE USE:

Temporary hickey marks are likely to appear for 7 days or less on the lips, around and/or inside the mouth on first-time users. The discolorations may vary from light pink to bluish purple. This is a common side effect. With continuous use of the lip plumper daily, the hickey marks will not return. If you stop using the lip plumper for a few weeks and come back to plump your lips again, the hickey marks may come back again. It is best to continue using the lip plumper at least once daily. Some firsttime users may get a blood blister. In this case, wait for it to heal completely before you plump again.



WHO SHOULD NOT USE CANDYLIPZ

Individuals who have sensitive skin; blood disorders, hypertension, diabetes, cancer, that bruise easily or are otherwise in poor health.

DO NOT APPLY CANDYLIPZ SUCTION PLUMPER ON YOUR NECK OR ANYWHERE ELSE ON THE BODY. IT IS ONLY INTENDED FOR THE LIPS/MOUTH.

Lip Plumper Styles and Sizes

STYLE	SIZE
	S-M
Model A	Accommodates horizontal lip widths under 2 inches for thin to average lip size
	S-M
Model B www.Canc	Accommodates horizontal lip widths under 2 inches for thin VLUP to average lip size
Model C	L Accommodates horizontal lip widths 2 inches and above for full lips, wide lips, and arched teeth
Model D	L Accommodates horizontal lip widths 2 inches and above for full lips, wide lips, and arched teeth

SHAPE	SINGLE-LOBED	DOUBLE-LOBED
(at a	No	(P)
		No DZ.COM
	No	
		No

www.CandyLipz.com

Xtreme Lip-Shaper[®] Set BENEFITS

Xtreme Lip-Shaper[®] Set gives amazing appearance of:

- 1. A fuller and poutier upper and/or lower lip
- 2. Luscious and pouty Single-Lobed or Double-Lobed lips
- 3. A fuller and poutier shape between the nose and the upper lip
- 4. A more contoured Cupid's bow
- 5. A shorter distance between the nose and the upper lip

- 6. Lifted corners of the mouth
- 7. Fewer lip wrinkles
- 8. Reduced lip drooping
- 9. Reduced lip pits and grooves
- 10. Restored aging lips
- 11. A natural rosy color to your lips

WHAT IS IN EACH BOX FOR MODEL A, B, C, OR D?

- Lip Pump
- Blocker
- Instructional Guidelines
- Mouth Cushion

ABOUT CANDYLIPZ® PREMIUM PUMP MATERIALS CandyLipz.com

- ✓ All pump parts that touch your lips are 100% food grade materials.
- ✓ All materials are eco-friendly, non-toxic, and ultra-hypoallergenic.
- ✓ Recycle all parts when you retire your pumps.

9

PRODUCT ACCESSORIES

HOW TO USE THE BLOCKER

Blocker	\bullet	Enhance the upper lip only
Blockér WU	v. 🕑 i	Enhance the lower lip only
Without a Blocker	Ð	Enhance both the up and lower lips simultaneously with a strong emphasis on the upper lip enlargement

HOW TO USE THE MOUTH COVER

The mouth cushion reduces red markings around the mouth, promotes suction strength, and enhances a tighter grip around the mouth.



SAFETY PRECAUTIONS

- Do not use on swollen, bruised, damaged, open wounds, or ulcerated lips and/or mouth.
- 2. Do not use if you have conditions which are prone to bleeding.
- 3. Do not use if you have lip implants or lip fillers.
- Do not share your lip enlargement pump with anyone for health and hygiene reasons!

www.CandyLipz.com WARNINGS

The mouth is the dirtiest part of the human body. Keep your lips, lip pump, and lip pump accessories such as mouth cover, blocker, and lid clean. If an infection occurs after your lip enlargement session, do ensure proper first aid and/or seek medical attention as appropriate.

CARE and STORAGE

- 1. Store your lip pump in a clean and cool, dry place.
- 2. Keep out of reach of children.
- 3. Wash your lip pump with mild dishwasher soap in cool water. Use a small toothbrush to clean the pump and mouthpiece. To remove water from inside the pump, roll up a piece of napkin and funnel the napkin into the pump. Allow the napkin to soak up the water then pull out the napkin. Repeat until all the water is soaked up and then allow the pump to/air dry before storage. A clean pump produces a better suction than a greasy, dirty pump. Wash your lip pump daily before use.
- Replace your mouth cover as frequently as possible. Wash your mouth cover in warm water with mild dishwasher soap daily.

CAUTION: Do not wash with alcohol as it can irritate your mouth.

WHAT TO DO AFTER YOU GET YOUR PUMP?

After you purchase your lip pump, remove the blocker from the top chamber and wash all lip pump parts with mild dishwasher soap in cool water before you use your lip pump.

INSTRUCTIONS FOR USE

STEP 1: Your lips must be clean and free of oil or grease! Remove all makeup on your lips: lipstick, lip-gloss, lip balm, dirt, gunk, etc. Massage your lips and mouth for a few minutes before you start.

STEP 2: Before you use the lip pump, make sure your lips are hydrated, smooth, soft, and pliable. Dry, cracked, or inflamed lips prevent the lip pump from working properly and having the desired effect. If your lips have lesions, sores, bruises, or cuts, they must be treated and healed before using the pump. **STEP 3:** Put the lip pump to your lips and test to see if the lip pump mouth opening is aligned to the center of your lips. If you are using the lip pump with a creaser mouth opening, you need to make sure that the creaser is centered at the midline of your lower lip



www.CandyLipz.com

STEP 4: Next, pucker your lips as far out as possible and hold them together. Learn not to drip saliva into your pump. Do not breathe into the pump or let air enter it. Breathe through your nose to avoid the problem.

STEP 5: Compress the lip pump as shown here (flip method), pucker and wiggle your lips as far into the lip pump as possible and place the middle platform between your upper and lower lips. Release the lip pump slowly, Relax your puckered lips while they are being pulled into the lip pump. If the suction is perfect, the lip pump will stick to your lips even when you are not holding it. You should be able to use the lip pump hands-free.

Tips: Here are a few more tips to help you create a good suction: 1) Dip the mouth piece into water as it helps create more suction (for practice only) 2) Use the mouth cushion or 3) Swallow some air to create a stronger suction.



STEP 6: If your suction is not strong enough and the lip pump falls off your lips, repeat step 4 and 5 until you achieve THE PERFECT SUCTION. The major reason why the lip pump falls off is because your lips and the skin surrounding the lips are not pushed deep enough into the mouthpiece and you did not create

a good suction. Try Again. Don't get discouraged. Practice makes better!

STEP 7: Once you have created the perfect suction, leave the lip pump in place for 1 to 2 minutes. If you lose suction during the session, try again. **IMPORTANT:** If you are a first-time user, read YOUR 10 DAYS CONDITIONING PERIOD.

STEP 8: If you need more time to fully enlarge your lips, take a 15-30 second break before you continue for an additional time. During your break, massage the applied area. Most people need only about two minutes to get the appearance that they want. The more time you add, the larger your lips will get.



www.CandyLipz.com

STEP 9: IMPORTANT: To remove the lip pump, compress while it is still attached to your lips by suction. Slowly and gently wiggle the pump until it comes off easily. Do not yank it off your mouth.

STEP 10: After your lips are enlarged, you can apply lip ointment, moisturizer, lipstick, or lip balm.

STEP 11: Your lips will gradually deflate in between one and one-half to two hours. You can repeat your 2 minutes lip enlargement application throughout the day to keep your lips plump and luscious.

IMPORTANT TIPS:

- If you prefer to work on just the lower lip alone, add the Blocker before you follow the rest of the steps indicated above.

- To create the correct suction, make sure the inner rim of the mouthpiece grips the area surrounding the lips firmly.

FREQUENTLY ASKED QUESTIONS

1. WHAT IS THE LIP CONDITIONING PERIOD?

It is a period where you must give your lips at least 7 to 10 days to adjust and adapt to the pulling suction of the lip plumper. During this period, you will get temporary hickey marks which is a common side effect on first-time users.

IMPORTANT: Some hickey marks may not appear right away. They may appear within 1 to 12 hours. Use less time and go at your own pace.

Day 1: Do 4 minutes total with 1 minute duration and 1 minute break in between.

Day 2 to day 3: Rest

Day 4 to 6: apply the product 1 session per day, each session will last for a total of 4 minutes with 2 minute duration and a 2 minute break in between.

Day 7 to 9: apply the product 1 session per day, each session will last for a total of 6 minutes with 2 minute duration and a 10 minute break in between. By this time, your hickey marks should be all gone. All individuals are different. Some individuals resolve their hickey marks in 4 days and some take up to 10 days. It is all normal. If the hickey marks linger on for 2 to 3 weeks, you need to stop using the lip plumper completely.

Day 10 and forward: You can plump your lips as often as you like. If your lips are sore, please rest a few days before you start again. Bruising should not occur again. However, if minor marks appear, they may be the last marks of your conditioning period.

TIP: During the break and rest time, gently massage the affected areas as frequently as you can to improve circulation.

2.DOES IT HURT?

You will feel a slight pinch and a slight suction pulling on your lips. The pain is minimal but your lips will 21

become accustomed to the suction after about one week. Thereafter, you will feel little or no pain whatsoever. If you do not like the strong pull of the suction on your lips, compress the pump at your comfort level until your lips are conditioned to the suction of the pump. Once you are accustomed to the suction, you can compress the lip pump at its maximum capacity. Thus, the stronger the suction, the larger your lips will be in a shorter period of time.

3. WILL THE LIP ENHANCEMENT PROCESS CAUSE LIP BRUISING?

95% of first-time users experience mild to moderate lip discolorations, hickey marks, or bruising. In the case that your lips bruise, you will need to give your lips a rest until the bruises go away. The bruises normally disappear within seven days. A great method to get rid of your discolorations quick is to massage the affected areas frequently throughout the day. Press moderately on the affected area with a knuckle using a back and forth motion. If your skin is a bit sore, give it a rest. Drink plenty of water with fresh squeezed lemon juice, raw vegetable juice, and eat more iron rich-foods during your conditioning period.

4. IS CANDYLIPZ SAFE FOR LIP PLUMPING?

Yes. Candylipz is clinically proven safe for lip plumping and is dermatologist approved

OUR CLINICAL TRIAL SHOWED THAT LIP VOLUME INCREASED BY 36 PERCENT IN 60 DAYS OF USE TWICE DAILY FOR 2 MINUTES EACH TIME.

5. HOW LONG DOES THE LIP ENHANCEMENT LAST?

According to the clinical trial, the appearance of enlarged lips may last up to two hours. You may use your lip plumper as frequently as you like throughout the day to maintain your pout. Note: Individual results may vary due to diet, health, and lifestyle.

6. HOW IMPORTANT IS NIGHTTIME APPLICATION?

Believe it or not, the best way to enhance the appearance of your lips and to help your lips look LUSCIOUSLY FULLER is to do a daily nighttime application! Before you fall asleep, make your lips look as large as you can by using the Single-Lobed Lip-Shaper® pump. Your lips will look enlarged until the morning because you are not talking or moving. Yes. Talking and moving can deflate your lips much quicker than when you are at rest. This is why lip injections might last just a few

weeks for some of you who have tried them. If you are a talker, you may have to plump your lips more frequently than individuals who talk less.

7. HOW MANY TIMES CAN YOU USE THE LIP PUMP PER DAY?

You can use the lip enhancement pump as many times per day as you like after 7 to 10 days of conditioning. Give your lips at least a 15-30 seconds rest period after each application. Remember to check your lips after one to two minutes to determine their size. Because the pump can be used hands-free, it is very easy for you to forget to remove the lip pump from your lips. **WARNING:** The appearance of your lips can become EXTREMELY LARCE!

8. HOW DO YOU KNOW IF YOUR LIP PUMP IS LEAKING AIR?

There are two methods to test your plumper. You can add a mouth cushion and swallow some air while you plump. If the plumper gradually loses suction, you may have a defective plumper.

9. IS DRINKING ENOUGH WATER IMPORTANT TO HELP THE LIPS APPEAR ENLARGED?

Yes. The appearance of your lips will stay enlarged longer if you are hydrated. If the appearance of your lips is not enhanced after two minutes, you are dehydrated. Stay hydrated. Take your weight in pounds and divide by two in order to determine how many ounces of water you should drink per day. Water, not coffee, caffeinated drinks or soda beverages! 20% of the water you consume will come from foods.

10. IS IT NORMAL THAT MY LIPS ARE BLUE AFTER I TAKE OFF THE PUMP?

Yes. The lack of oxygen will turn your lips blue. This symptom will disappear in a few seconds. If your lips stay blue/purple for a longer period of time (such as days) after you use the lip plumper, add iron rich foods and lemon juice to your diet until your bluelips symptom goes away. Stay hydrated. You can still plump during this period. Do half the time you normally would do until your blue lips restore to their natural color. Massage the lips and mouth tissues frequently to improve circulation.

11. DOES IT TAKE A LONGER TIME FOR ELDERLY INDIVIDUALS IN THEIR LATE 50s TO OBTAIN THEIR RESULTS?

Yes. It does take a longer time for elderly individuals in their late 50s to obtain their results because their lip skin is not as hydrated, pliable, and flexible as in younger individuals. However, with continuous use of the lip enhancement pump, the lip skin may appear to

rejuvenate and become more elastic. When this occurs, it may take less time to enlarge the appearance of lips in elderly individuals.

12. IF MY BLOCKER FALLS INTO THE UPPER CHAMBER OF THE PUMP, HOW DO I GET IT OUT?

You can use tweezers to pull the Blocker out.

13. ARE THERE DRUGS THAT CAN SLOW MY LIP ENHANCEMENT PROCESS OR CAUSE THE LIPS TO DEFLATE FASTER?

Yes. Anti-inflammatory drugs or pain killers (even caffeine and tobacco smoke) will deflate your lips very quick! Avoid if possible.

To learn important tips and tricks, go to www.blog.candylipz.com

14. ARE TRAINING VIDEOS AVAILABLE?

To get your full training course on how to use the Xtreme Lip-Shaper[®] Set, use your smart phone to scan in the QR code or you can visit www. CandyLipz.com/videos.



BEAUTY TIPS

- Some individuals may experience a little pinch from 1 using the Xtreme Double-Lobed Lip-Shaper® Pump. To make your experience more pleasant, use the Xtreme Single-Lobed Lip-Shaper® Pump to enhance the appearance of your lips first. Once you are satisfied with your lip size, use Double-Lobed Lip-Shaper® Pump to add the Double-Lobed lip appearance. Experiment and see what you like. With continuous use of the Double-Lobed lip plumper, you will get used to it and may not feel the pinching. To maintain the vertical crease on the lower lip, you have to practice as frequently as you can. The more you use this style, the longer the appearance of the double-lobe look may stay.
- 2. If you like your lips to appear extremely large, use the Xtreme Single-Lobed Lip-Shaper[®] Pump first to enhance the appearance of your lips. Once you are satisfied with your lip size, use Double-Lobed Lip-Shaper[®] Pump to add the Double-Lobed lip appearance.

- Drinking cold water may cause your lips to shrink faster. Drink warm water to help your enlarged lips appearance last longer.
- 4. If you want your upper lip to be larger only, use the Blocker to block suction to the top chamber and revert the lip plumper before you use as shown in the image here.



5. Always plump your lips before going to sleep.

BEFORE AND AFTER GALLERY All images in gallery are unedited.



www.CandyLipz.com



RETURN POLICY

All sales are final due to health and hygienic reasons. If your plumper is defective within 30 days of your purchase, please send support@candylipz.zendesk.com an email and we will send you a replacement lip plumper. Thank you so much for using CandyLipz.

www.CandyLipz.com

CONTACT

CandyLipz[®] customer service is available from 9 A.M. to 5 P.M. PST Monday to Friday to answer any questions and concerns you may have. If the information available to you via our instructional guide, website, and/or training videos do not answer all your questions, please send us an email to support@candylipz.zendesk.com and one of our CandyLipz[®] representatives will get back to you within 24 hours via email or a live phone appointment call back. Please list your questions and indicate whether you want to schedule an appointment or you just want a returned email with the answers to your questions. Due to high call and email volumes, we will not respond to your requests if the answers are available through the sources listed above. We truly appreciate your time. Thank you for your support.

DISCLAIMERS

The safety of this product has been clinically proven safe for intended use. It is Dermatologist Approved. Please follow the instructions carefully in order to achieve the results you desire without any adverse effects. Consult professional medical advice when appropriate. Xtreme Lip-Shaper[®] is not a medical device. Statements made by CandyLipz[®] LLC. have not been evaluated by the Food and Drug Administration (FDA). CandyLipz[®] LLC. products are not intended to diagnose, treat, cure or prevent any illness or disease.

Country of Manufacture: China Manufactured for and Distributed by CandyLipz® LLC. San Francisco, CA 94127 Website: www.CandyLipz.com IG: @mycandylipz YT: www.youtube.com/user/MyCandyLipz/videos FB: www.facebook.com/mycandylipz Email: support@candylipz.zendesk.com Tel: (844)CANDYLIPZ







