LUXE.

LED FACE REJUVENATION

LUXE LED MASK device designed for home use. Its purpose is to emit a consistent, narrow-band of light for the treatment of facial cosmetic concerns. The mask is worn on the face and secured with an adjustable Velcro

LUXE LED MASK has been designed to deliver light directly to the skin, maximizing the absorption by the target cells and minimizing waste. It utilizes 7 clinically proven wavelengths of light, each shown to stimulate collagen and elastin production.

become pregnant, or breastfeeding.

• Do not use the LUXE LED MASK if you have Lupus erythematosus, photosensitive eczema, or albinism. Using the LUXE LED MASK

SPECIFICATIONS LED Light: 7 LED colours • Near Infrared Light: 810-850 nm LED Lamps: 121 lamps Working Time Settings: 10-30 minutes

 AC/DC Adapter Output: 5V, DC 1000mA Light Wavelengths and Effects: • Near Infrared Light (810-850 nm): Anti-aging and firming, reduces to remove these products may cause a skin reaction or diminish

treatment results. • Red Light (640-670 nm): Reduces fine lines and wrinkles, anti-aging
• Do not use the LUXE LED MASK if the mask, controller, charging Green Light (490-540 nm): Soothing and purifying effects plug, or cord is damaged. • Blue Light (440-470 nm): Targets acne and skin inflammation, • Stop using the LUXE LED MASK immediately if it feels hot or if anti-inflammatory • Yellow Light (565-570 nm): Decomposes pigment, improves rough |
• Do not use the LUXE LED MASK near water or in environments where water might come into contact with the device.

 Purple Light (380-440 nm): Repairs acne scars, acne Cyan Light (460-495 nm): Gradually enhances cellular energy Laser Light: Accelerates cell activity, decomposes dark spots

FIRST-TIME SETUP

WARNINGS

headaches.

contact with water.

become pregnant, or breastfeeding.

with these conditions may cause a severe skin reaction.

Do not use LUXE LED MASK if you experience light-induced

Do not use LUXE LED MASK if you have any genetic eye conditions.

Always remove creams, lotions, or makeup before using the

LUXE LED MASK. Ingredients in these products may cause light

Do not use the LUXE LED MASK if you suspect it has come into

• Do not attempt to modify the LUXE LED MASK in any way.

sensitivity or reduce the effectiveness of the treatment. Failure

conditions may result in a severe skin reaction.

your skin becomes uncomfortably warm.

The LUXE LED MASK is a wearable, rechargeable LED phototherapy

• Do not use LUXE LED MASK for any conditions other than those

Unpacking and Charging: listed under the indications for use. The LUXE LED MASK has not • Before using the LUXE LED MASK, it is important to fully charge environment. been tested for other conditions than those listed, and the the battery. associated risks are unknown. • The indicator is red light on the control to show that the device is damp and humid areas. • Do not use LUXE LED MASK if you are pregnant, planning to

photosensitive eczema, or albinism. Using the LUXE LED MASK | OPERATION

 Do not use the LUXE LED MASK if you have any photosensitive products. disorder or sensitivity to light. Using the device under these 2. Put on the mask (recommended use is at least once a week). to the manufacturer or an authorized representative at the end 3. Start the treatment using the controller.

Charging takes approximately 4-5 hours.



ON/OFF Press power button for 1-2 seconds Near Infrared Light Button **E** LED Intensity 3 levels of LED intensity Setting Treatment Time + Setting Treatment Time -Treatment times 10 / 15 / 20 / 25 / 30 minutes CARE AND MAINTENANCE

 Store the LUXE LED MASK at room temperature in a dry Avoid storing the mask in dusty conditions, direct sunlight, or

charging. The battery is fully charged when the light turns blue.

• Before each use clean the mask following the cleaning instructions provided, ensuring the device is switched off during this process. Do not use any cleaning products that are not recommended, as

they may damage the mask's materials. 1. Cleanse and exfoliate your skin to remove any traces of skincare . Do not dispose of any part of the LUXE LED MASK with regular household waste. To protect the environment, return the device of its life, or follow local guidelines for safe disposal.

> WARRANTY The LUXE LED MASK is warranted for 1 year under normal use, and

the battery is also warranted for 1 year. • This warranty applies only to the original purchaser. If the product becomes faulty within the warranty period, contact your provider. • The authorized provider reserves the right to repair a defective product. The authorized provider also reserves the right to replace a

defective product with the most compatible product currently available. This warranty does notcover damage caused by:

 Accidents, misuse, abuse, or alterations. Servicing by unauthorized persons. Connection to incorrect power sources.

Use with unauthorized accessories.