


# bio roller<sup>®</sup> G4




The brand of BIO roller integrate the advantages of the existing roller on the market. It makes significant innovation and improvement. BIO roller is aimed to improve the skin problem in order to bring you beauty and make your life full of joy.


# If you have these troubles



Bald Patch



Wide Crevice



Skin Care

## Advantages



### Multi-directional Portable Operation

Easy to use longitudinally, vertically and diagonally, with even force and flexible rotation




Integrated cover, ready to use



Turn over → sterilizer tank



1200 pins



1200 stainless needles create thousands of microchannels that activate the body's natural healing process, increase blood circulation and stimulate sleeping hair follicles.




# For Hair Growth & Beard Growth

G4 can penetrate the epidermis and create thousands of microchannels to initiate skin tissue repair, stimulate cell regeneration, and promote scalp renewal. At the same time, it stimulates hair follicles, revitalizing them and providing a healthy microenvironment. Use with hair growth fluid to double absorption and continue to provide active energy and nutrients for maximum protection against hair loss.




## Hair Growth



## Scalp Renewal

Delay scalp aging, accelerate epidermal cell renewal, and restore scalp vitality



## Beard Growth



# For Skin Care

G4 is great for large area skin care. It penetrates the epidermis and creates thousands of microchannels to initiate skin tissue repair, stimulate cell regeneration, rejuvenate the skin, shrink pores and restore skin elasticity.



Stretch mark on belly




Large areas of acne



Stretch mark on thighs

# Product Parameters



## ● Specification

Brand	Bio Roller	Model	G4
Shell Material	ABS	Number of Needles	1200 stitches
Needle Material	Stainless steel	G.W.	7.2*7.2*7.2 cm/ 1 pc/ 100 g
Needle Length	0.25mm、0.5mm、0.75mm、1.0mm、1.5mm、2.0mm	1 carton	55*37*25 cm/105 pcs/11.5 kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box ② G4 (1 pc)

### How to use :

1. Disinfect the roller with 75% alcohol.
2. Clean the target skin regions thoroughly.
3. Apply skincare product or serum upon the skin region and start rolling gently over the skin.
4. Clean and disinfect the roller after use.