



# User's Guide

Thank you for purchasing the Body Back Vibe Massager! This easy to use massager is great for either home or professional use. The comfortable, oscillating vibration of this massager allows you to ease muscle tension, fatigue, and soreness by yourself or with a partner. The adjustable speed of the Body Back Vibe gives you the control. Let your Body Back Vibe do the work for you and feel pain and discomfort melt away!

# Body Back Vibe Massager

---

The Body Back Vibe Massager is ideal for either professional or home use. Its soothing and comfortable vibrating effect quickly relieves tight muscles and overall tension. This powerful tool has many therapeutic advantages. It may increase circulation, improve digestion when used around the abdominal area, reduce cellulite, balance hormones and much more!

- Helps Calm Deep Muscular Pain
- Relieve Muscle Spasm
- Deluxe Oscillating Movement Increases Blood Circulation
- May Improve Digestion and Reduces Stress
- Perfect for Professional or Home Use

## Operating Instructions

---

1. Plug electrical cord into outlet.
2. Turn on massager using switch.
3. Adjust speed to desired setting.
4. Place massager in contact with body - solo or with a partner.
5. Do not exceed 20 minutes of use to avoid overheating massager.
6. To turn off, move switch into the neutral, center position.
7. Unplug massager

# Great For The Whole Body

---



**Back**



**Shoulders**



**Glutes**



**Arms**



**Hands**



**Feet**



**Chest**



**Abdomen**

## Accessories

---

A strap is provided to safely and securely position the Body Back Vibe Massager in place if you are self-massaging. Fasten the strap under your back to massage your abdomen or anywhere else you need stability.



# Important Safety Instructions

---

Please read all safety and usage instructions before operating massager.

- Use this massager only according to the operating instructions listed in this manual. Unintended use may result in injury to the user or damage to the massager.
- DO NOT USE NEAR WATER OR WET SURFACES
- Close supervision is advised when operating the massager, especially around children, elderly and disabled persons.
- Never leave the massager unattended while in use.
- Unplug when not in use
- Do not operate if there is any damage to the massager, cord, or plug-in.
- Do not operate if there is an abnormal sound from the massager.
- Do not insert or attach any foreign objects in or on the massager.
- Do not operate under a blanket, bedspread or in any place where airways could become blocked.
- This massager is intended for indoor use only - do not use outdoors.
- Operate away from heat sources.
- Do not use aerosol products near the massager or operate close to where oxygen is being administered
- Let the massager cool after use
- Do not use for more than 20 minutes per session.
- Clean with a dry cloth or towel after use.
- Check with physician before use especially people who have hypertension, heart or kidney problems, or who are pregnant.

## Warranty

---

We offer a one year warranty on the Body Back Vibe Massager. If you should have a problem with the massager for any reason, please contact us at [shipping@bodyback.com](mailto:shipping@bodyback.com). Should your massager be in need of replacement or repair, please be prepared to show proof of purchase.

All misuse, unauthorized repair, or damage from abuse will not be covered by the warranty.

## Product Specifications

---

Voltage: 110V

Measurement: 235x120x185 mm

Rotate Speed: High 3500 round/min +/- 10%, Low 2800 round/min +/- 10%

Weight: 3.5 kg



Body Back Company  
280 E. Hersey St. #19  
Ashland, OR 97520  
[www.bodyback.com](http://www.bodyback.com)

### WARNING:

Use caution when wearing loose clothing, jewelry, or loose/long hair when operating this massager as they may get caught in the massager while in use.