

cellpod®

bubbolux

Silicone Vibration Brush &
Pore Refining Deep Cleansing Foam EX



cellpod is a beauty device brand by Ecowell, designed to make aesthetic care easily accessible at home.

cellpod



It's different when cellpod makes it

cellpod is the beauty device brand of Ecowell Co., Ltd. a research-oriented company specializing in skin technology. We have developed technologies to address customers' skin concerns through systematic research and innovation led by various experts with doctoral degrees.

With our expertise and technological prowess, we not only offer solutions under the cellpod brand but also provide customized solutions for industry partners to implement their own specialized beauty devices.

Through collaboration with industry partners, we aim to grow together and provide optimized skincare solutions to a broader range of customers.





FOR THOSE WHO

- ✓ Individuals seeking to thoroughly remove makeup residue.
- ✓ Those searching for an easy-to-use daily cleansing tool.
- ✓ Anyone wishing to soothe their skin after extensive outdoor activities.
- ✓ Individuals with sensitive skin in need of gentle cleansing.
- ✓ Those struggling with the effective absorption of cosmetics or skincare products.
- ✓ Anyone looking for an easy and convenient home care solution.

With soft and delicate vibrations for a cleaner and more vibrant finish

- ✓ Makeup residue removal
- ✓ Pore management
- ✓ Skin texture improvement




SUPERIOR WATERPROOFING AND CLEANSING EFFICACY



Gentle and delicate vibrations

Delivers approximately
12,000 vibrations per minute
(about 12,000 RPM) to gently remove
impurities deep within the pores.



**12,000
RPM**

cellpcd

SGS TEST REPORT

Report No: SGS-R23-2612-EN00
Page : 1 / 7

Report No. : SGS-R23-2612-EN00
Order No. : DTRC2311002375EN
Applicant : ecowell
Address : 1107, 11th floor, Hanhwa Biz Metro Bldg(1Cha), 551-17, Yangcheon-ro, Gangseo-gu, Seoul, Korea
Product : SILICONE VIBRATING CLEANSER
Model : EBI-E11114
Date of Test : October 26, 2023.
Location of Test : Permanent testing Lab. On Site Testing
(Address : 20, Dongtansandan 10-gil and 12, Dongtansandan 10-gil, Hwaseong-si, Gyeonggi-do, Korea)
Standard / Method : IEC60529:1989+A1:1999+A2:2013
Degrees of protection provided by enclosures (IPX8)
Test Result : Refer to page 7.
Use of Report : Quality evaluation

This is certified that the above mentioned products have been tested for the sample provided by client.

Affirmation	Tested by Name : Son, Joon-hyung	Technical Manager Name : Kim, Hak-ju

■ No part of this document may be duplicated or reproduced by any means without the express written permission.
■ The test report is not related to our KOLAS accreditation.

November 30, 2023.

SGS Korea Co., Ltd. Dongtan Laboratory



This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at www.sgs.com/eng, and, conditions that attention is drawn to the inclusion of results, identification and production issues defined therein. Any holder of this document is advised that information contained herein reflects the Company's findings at the time of the inspection only and within the limits of Client's instructions. If any. The Company's sole responsibility is to its Client and the document does not constitute proof of a transaction from any party at the time and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company, any unauthorized alteration, forgery or modification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated the results shown in this test report refer only to the samples tested and such samples are retained for 90 days only.

RGP 009-P02 (02)
SGS Korea Co., Ltd. | 20, Dongtansandan 10-gil, Dongtan-myeon, Hwaseong-si, Gyeonggi-do, Korea | T+82 (0)31 340 9011 | F+82 (0)31 340 9009 | www.sgs.com | www.sgs.com/kr
Member of SGS Group (Societe Generale de Surveillances)

SMART&COMPACT



Smart operation

Automatically activates when a finger is inserted into the strap, eliminating the need for an on/off switch.

Exceptional portability

A compact size that fits comfortably in one hand (61mm70.5mm28mm). Comes with a cellpod pouch for convenient carrying.



Soft silicone brush

A **hygienic** and gentle **non-toxic** fine bristle silicone brush that does not irritate the skin.



Optimized waterproof internal structure

Finger support band with automatic detection sensor



Lithium coin battery case




Secondary silicone waterproof cover



COMPARISON OF CLEANSING EFFICACY WITH COMPETITOR PRODUCTS

The comparison between Bubbleux Silicone Vibration Cleanser and FOREO LUNA MINI showed similar performance, with both products demonstrating twice the cleansing speed compared to manual washing.




REASONABLE PRICE RANGE



Mass-produced Chinese products, typically priced between approximately \$6 to \$25, are available on various open markets and direct purchase channels from China. However, these products are generally perceived to have poor quality.



Including a dedicated foam cleanser, a reasonable price range can be formed from the late \$23 to the early \$32.



Based on web searches for well-known brand products, the price range is formed from approximately \$36 to \$266.

OPERATES WITH A SINGLE CR2032 BATTERY, ELIMINATING THE NEED FOR CUMBERSOME CHARGING

Bubbleux Silicone Vibration Cleanser is powered by a single CR2032 battery.



Based on a one-minute daily use, it can last for approximately 60 days, and for two-minute daily use, around 30 days.


*Actual duration may vary depending on individual usage frequency.

competitor	OA	JUVERA	FOREO
Charging Time	150 mins	300 mins	120 mins
Usage Time	40 mins	70 mins	90 mins

How to Replace the Battery

① Upper part ② Lower part ③ Battery cover

1. Pull apart the bottom of the product to separate it.
2. Remove the battery cover and replace the battery.
3. Ensure the positive side is facing up, then reassemble in reverse order.



PORE REFINING DEEP CLEANSING FOAM EX



A high-moisture cleanser with bouncy foam and firm elasticity

A nourishing, resilient texture with elastic foam and cushioning effect to minimize skin irritation.

Gentle and smooth cleansing

Gently exfoliates with cellulose, removes residual impurities with soft and dense foam.

Natural-derived ingredients for moisturizing and soothing effects



Extracts of Muhuza



Extracts of Soapberry,



Extracts of Quillaja



Extracts of Yucca

Cleansing pollutants penetrating into the skin



Extracts of Macadamia



Fermented extract of Water hyacinth



Extract of chrysanthemum

HOW TO USE



1

Apply an appropriate amount of cleansing foam onto the face along with a moderate amount of water.

Gently massage the skin from the inside to the outside in circular motions using your fingers wrapped around a finger strap.



2



3

After use, remove your fingers from the finger strap and rinse with water.

SET COMPOSITION



1 silicone vibrating brush



cellpod portable pouch



2 units of Pore Refining Deep Cleansing Foam EX

*Please direct inquiries regarding silicone vibrating brushes or foam cleansing products to the sales department of EcoWell Co., Ltd. through their contact person, email, or our company's website.

SILICONE VIBRATING BRUSH & PORE REFINING DEEP CLEANSING FOAM EX REBORN UNDER YOUR BRAND

**Manufactured as OEM/ODM and distributed
and sold under your own brand**

Customization options are available to match your brand, including printing your brand logo, changing product and package colors, and adjusting the formulation of the foam cleanser to align with your brand.

For further details, please reach out to the designated person at the Ecowell Sales Department via phone, email, or by submitting an inquiry through our website.

EMAIL: sales@ecowell.co.kr

Homepage: ecowell.co.kr



PRODUCT SPECIFICATIONS



Item	Details
Battery	CR2032 / 1EA
Continuous usage time (actual usage time)	50 mins (75 mins)
Motor	Bar-type vibration motor
Function	Micro-vibration
ON/OFF	Magnetic switch
Waterproof rating	IPX8

FAQ

Q. Where is the power button located?

A. Bubblux does not have a separate power on/off switch. When you insert your finger between the strap and the device, the vibration will start automatically.

Q. What is the battery life of this product?

A. This product uses one CR2032 battery. Based on usage of 1 minute per day, you can expect to use this product for approximately 60 days with one CR2032 battery. With 2 minutes of daily usage, it would last for around 30 days. However, actual usage time may vary depending on individual usage frequency.

Q. After using the product, how should it be stored?

A. After use, wipe off any moisture and store the product with the finger strap facing downwards. This helps prevent battery drainage and keeps the silicone brush cleaner.