cellpod®

# ELB

EMS LED BOOSTER

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


출처: Yan, Chadakan & Ng, Janice. (2022). Efficacy of handheld iontophoresis device in enhancing transdermal vitamin C delivery: A split-face clinical trial. Journal of Cosmetic Dermatology. 21. 10.1111/jocd.14702.


출처: Effect of At-home Beauty Device Equipped with Microcurrent and Electrical Muscle Stimulation Functions on the Elasticity and Lifting of Facial Skin So-Young Kwak, Da-Som Shin, Yu-Kyung Hur, Su-Gyeong Hong, and Ye-Ji Choi\* Researcher, APR Corporation R&D Center

출처: Effects of a skin-massaging device on the ex vivo expression of human dermis proteins and in-vivo facial wrinkles Elisa Caberlotto, Laetitia Ruiz, Zane Miller, Mickael Poletti, Lauri Tadlock


출처: Application of 630-nm and 850-nm Light-emitting Diodes and Microcurrent to Accelerate Collagen and Elastin Deposition in Porcine Skin

출처: Consider of Micro-Current's effect to variation of Facial Wrinkle trend, Randomized Clinical Trial Study. Fatemeh Saniee, Hamid Reza Ghafarian Shirazi, Khosroo Khademi Kalantari, Parviz Yazdanpanah, Asghar Reza Soltani, Nader Dabiri, Nazanin Ghafarian Shirazi, Farzad Karimpour

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