



VANITY PLANET

Leda

Blue LED Sonic Facial Brush

USER GUIDE

Cleanse sensitive, blemish-prone skin with Leda, the sonic powered face brush with silicone bristles and therapeutic LED Blue light therapy. Cushion your complexion with soft, deep-cleansing massage and the reparative nature of FDA-approved blue lights. Remove pore-clogging impurities with daily treatments of blue light wavelengths to kill blemish-causing bacteria without harmful effects to the skin.

Includes

- Blue LED Sonic Facial Brush
- 2 Interchangeable Silicone Treatment Heads
- USB Cable
- User Guide



Features

SILICONE FOR SENSITIVE SKIN

Gentle silicone cleansing brush designed for sensitive skin

SAFE AND EFFECTIVE

Blue LED lights kill blemish-causing bacteria without harmful side effects

FDA-APPROVED LIGHT THERAPY

Blue LED light is proven to help treat mild-to-moderate acne

3 SONIC BRUSH MODES

Customizable massage modes deep clean for a nourished complexion

Instructions

LIGHT SENSITIVITY TEST

- Before initial treatment, test yourself for light sensitivity
- Press the power button to power up Leda, and place Leda on or within 1/4" of your skin for 3 minutes
- If the treated area has turned red and the redness lasts for more than 2 hours, your skin is light sensitive
- Discontinue use of Leda and contact our customer service

Instructions

SONIC CLEANSING MODE

- Make sure Leda's battery is fully charged before usage
- Remove makeup with a cleansing towelette or wipe before using the attachments to avoid staining the bristles
- Attach the silicone brush head
- Dampen your skin with warm water
- Apply cleanser directly to your skin or the brush head. Avoid using exfoliators with large, jagged particles
- Press the power button to power up Leda, then press the power button to select a mode
- Gently glide Leda over your skin's surface in a rotating motion

LIGHT TREATMENT MODE

- Make sure Leda's battery is fully charged before usage
- Clean and dry treatment area prior to use
- If using serums, attach the silicone massage head
- Press the power button to power up Leda, place directly on the treatment area, and Leda will automatically shut off after 3 minutes
- Move Leda to the next treatment and repeat
- Press and hold the power button to turn Leda off when finished
- Use daily for best results, but don't exceed 3 treatments a day for each treatment area

Charging

- Twist and remove light head to charge
- Using the provided USB cable, connect the USB port to a 5V power supply or computer
- Charge your Leda for 2-4 hours

Upkeep

CLEANING

- Turn Leda off and unplug from any power source
- Remove brush head, rinse and dry
- Wipe device surface with a damp cloth or alcohol wipes

MAINTENANCE

- Leda is a medical device and it is important to keep it clean
- Leda has no user serviceable parts
- All settings are set at the factory

DISCLAIMER

Leda is not intended to diagnose, cure or prevent any disease. We makes no claims, representations or warranties regarding the ability of Leda to cure any physical, skin or mental condition. A qualified health professional should be consulted with regards to any condition requiring medical attention. This product contains a lithium-ion battery that must be disposed of properly, according to local laws and regulations, and not as household waste. When the product reaches its end of life, check with your local state and municipal waste disposal regulations, to learn about disposal and recycling options.

Safety

- Stop using Leda if you think it is defective. Using a broken device could result in injury
- Store Leda out of direct sunlight to prevent damage to the plastic casing
- If Leda overheats or becomes uncomfortably hot to your skin, stop usage
- Do not exceed recommended treatment times or usage
- Do not use Leda on infants or young children without physician approval
- Do not shine the device lights directly into your eyes
- Consult with a doctor prior to use if you suspect that you have a silicone allergy
- Contact your doctor prior to using if you are pregnant or suspect you may be pregnant or have suspicious or cancerous lesions

Contraindications

- Do not use this device if your skin is light-sensitive or you are currently using medication that may cause skin sensitivity to light. Doing so may result in photosensitive reaction.
- If unsure whether your medication causes sensitivity to light, confirm with your physician prior to using this device.
- Discontinue use immediately if you experience any sensitivity to using this device.

Warnings

- Do not use if the device has been exposed to water or other liquids, if it is defective or is damaged in any way. Doing so could cause an electric shock or result in injury. Please contact Customer Service at info@vanityplanet.com if the device has been damaged in any way.
- Do not attempt to take the device apart. Any modification or attempted service will void the warranty.
- Do not store the device in the sun or any location with high temperatures as this may cause damage to the silicone and plastic casing.
- Discontinue use and contact Customer Service if the device overheats or becomes unconformably hot to your skin.
- Do not exceed recommended treatment times or usage.

Precautions

- Do not use this device on infants or young children without physician approval.
- Do not shine the device lights directly into the eyes. Doing so may cause temporary vision discomfort.
- This product contains medical-grade silicone. If you have or suspect to have allergies to silicone, consult your doctor prior to use.
- Contact your doctor if you have any medical questions regarding use of this device.
- Consult your physician prior to using this device if you:
 - Are pregnant or suspect you may be pregnant.
 - Have suspicious cancerous lesions.

One-Year Warranty

Visit vanityplanet.com/warranty to register your product and activate your warranty.

This Vanity Planet product is warranted to be free of manufacturer's defects in materials or workmanship for one year from the date of purchase. Damage or wear resulting from an accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair are not covered under this warranty.

Should this product require warranty service (or replacement at our discretion), please contact client services. Proof of purchase is required. Products returned without prior authorization from Vanity Planet will not be accepted and the sender will not receive a refund or replacement.

vanityplanet.com
info@vanityplanet.com

There are no express warranties except as listed above. This warranty gives you specific legal rights and you may have other rights which vary from state to state.

Please do not return this product to the retailer. Contact us directly for warranty service.



VID

Distributed by Dastmalchi, LLC
Newport Beach, CA 92660
Made In China