

# Veil

LED Acne  
Spot Treatment

What if you could stop acne before it started? Veil features both blue and red LED light to eliminate bacteria that causes breakouts while reducing inflammation and promoting skin health.

## Contraindication

- Contact your doctor prior to using the device if you are or suspect you are taking medications known to cause light sensitivity.
- Do not use this device if your skin is light sensitive or if you are currently using medication that may cause skin sensitivity to light. Doing so may result in a photosensitive reaction.
- Discontinue use immediately if you experience any sensitivity while using this device.

## Warnings

- Do not operate the device where it could come into contact with water. Doing so could cause an electric shock or result in injury.
- Stop using the device immediately if you think it is defective. Using a damaged device could result in injury.
- Do not attempt to take the device apart. Any modification or attempted service will void the warranty.
- Do not store the device in the sun or a location with high temperatures as this may cause damage to the plastic casing.

## Precautions

**Do not use this device on infants or young children without physician approval.**

- Do not shine the light directly into your eyes. Doing so may cause temporary vision discomfort.
- Contact your doctor if you have any medical questions regarding the use of this device.
- Contact your doctor prior to using Veil if you:
  1. Are pregnant or suspect you may be pregnant.
  2. Have suspicious or cancerous lesions.

## Light Sensitivity Test

1. Power on your device.
2. Place the device on or within 1/4" of the skin on your forearm for three minutes.
3. After three minutes, check the area. If it has turned red and lasts for more than two hours, your skin is light sensitive.
4. Discontinue use.

## Includes

- Veil LED Acne Spot Treatment
- User Guide
- Strap

## Instructions

1. We recommend washing and drying the treatment area prior to use.
2. Insert one AAA battery (sold separately) and gently twist to close the cover.
3. Turn on device by gently pulling the two halves apart (without twisting).
4. Place the device directly on the skin over the desired treatment area for three minutes. Do not move or wave the device across the treatment area. The device will automatically shut off after the 3 minute treatment is finished.
5. Move the device to the next treatment area and repeat step four.
6. Turn off the device after treatment by gently pushing the two halves back together.

**Use the device daily for best results.**

# Veil

## LED Acne Spot Treatment

User guide

### Disclaimer

This product is not intended to diagnose, cure or prevent any disease.

Vanity Planet makes no claims, representations or warranties regarding the ability of this product to cure any physical, skin or mental condition.

A qualified health professional should be consulted with regards to any condition requiring medical attention.

### Upkeep

#### Note:

1. Because this is a medical device, it's important to keep it clean.
2. There are no user-serviceable parts in the device.
3. All settings are set at the factory.

#### Cleaning:

1. Turn the device off.
2. Wipe clean with a damp cloth or alcohol wipes.
3. Allow it to completely dry prior to use.

### Troubleshooting

#### Veil won't turn on:

- Is the battery fully charged?
- Is the battery aligned correctly?
- Is the device turned on?
- Is the device cap (at the battery end) fully closed?

### Warranty

Visit [vanityplanet.com/warranty](http://vanityplanet.com/warranty) to register your product and activate your warranty.

This Vanity Planet product is warranted to be free of manufacturer's defects in materials or workmanship for one year from the date of purchase. Damage or wear resulting from an accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair are not covered under this warranty.

Should this product require warranty service (or replacement at our discretion), please contact client services. Proof of purchase is required. Products returned without prior authorization from Vanity Planet will not be accepted and the sender will not receive a refund or replacement.

[vanityplanet.com](http://vanityplanet.com)

[info@vanityplanet.com](mailto:info@vanityplanet.com)

1.888.870.2132

There are no express warranties except as listed above. This warranty gives you specific legal rights and you may have other rights which vary from state to state.

Please do not return this product to the retailer. Contact us directly for warranty service.