

JP51/51M/200/210/250 Solo Quick Start Guide

JETPIK
Power Floss System



Step by Step

STEP 1

Insert one Floss Capsule into the nozzle



STEP 2

Place the nozzle into the handle completely



STEP 3

Connect the water supply hose to the handle and fill cup with water



STEP 4

Press On/Off button to let the water blow the floss out into nozzle (prime the floss)



STEP 5

Put the nozzle toward your teeth and slightly close your lips. Lean over the sink and press On/Off button to start cleaning



Keep the tube and travel lid dry if you will not use it for a while will hold your system free of germs and mold build up!

JP51/51M/200/210/250 Solo Quick Start Guide

JETPIK Power Floss System



JETPIK web site: www.jetpik.com