JETPÌK



STEP 1

Insert one Floss Capsule into the nozzle



STEP 4

Press On/Off button to let the water blow the floss out into nozzle (prime the floss)

STEP 2

Place the nozzle into the handle completely



STEP 5

Put the nozzle toward your teeth and slightly close your lips . Lean over the sink and press On/Off button to start cleaning



STEP 3

Connect the water supply hose to the handle and fill cup with water

Use JETPIK UV Sanitizer regularly to kill bacteria and keep nozzle and toothbrush clean



JETPIK web site: www.jetpik.com

Quick Start Guide

JETPÌK Power Floss System

JETPIK web site: www.jetpik.com



