


fancii

Good-bye  
Boring  
**Beauty  
Routine!**

# LILY

EMS/RF Wand with  
Light Therapy

- 1 / Massage Head
- 2 / LED Color Light
- 3 / Gear Light
- 4 / EMS/RF Switch
- 5 / Color Lighting Switch



## FEATURES

**2-IN-1 FACIAL TONING** - Elevate your skin treatment with Electric Muscle Stimulation (EMS) that increases facial muscle vitality, eliminates puffiness, improves facial firmness, and restores the elasticity of your skin and Radio Frequency (RF) which helps firm, tighten, and lift the skin while stimulating and promoting collagen production.

**LIGHT THERAPY** - Target and treat a variety of skin conditions using five different light emitting diode (LED) wavelength options.

**SAFE, PAINLESS, AND EFFECTIVE** - Heal and treat your skin in the comfort of your own home and enjoy zero recovery downtime. Suitable for all skin types including sensitive skin.

**ECO-FRIENDLY POWER** - Lily is designed with a rechargeable battery that lasts up to 3 hours of use with a full charge.

## WARNINGS

- ⚠ Read **Warnings** carefully before using the device.
- ⚠ Children must not play with this device
- ⚠ Cleaning and user maintenance may not be carried out by children without supervision
- ⚠ Use this product only for its intended use as described in this manual.
- ⚠ Discontinue use immediately and contact Customer Service if you notice any malfunction or unusual reaction to this product.
- ⚠ Do not use this product in a defective state as this could result in injury.

## DISCLAIMER

⚠ Fancii makes no claims, representations, or warranties regarding Lily's ability to diagnose, prevent or cure any physical or mental condition. Consult a qualified health professional regarding any condition that requires medical attention.

⚠ Conduct a light sensitivity test before usage of the device. See **Let's Begin** for more information.

## SAFETY PRECAUTIONS

1 / The following persons should not use this device unless approved by a medical professional:

- Persons who are currently undergoing treatment, have acute medical conditions, tubercular conditions, or have recently undergone surgery.
- Persons with heat-sensitive skin or skin diseases such as rashes, atopic dermatitis, etc.
- Persons with conditions such as light sensitivity, photophobia, skin hemorrhages, vitiligo, or other severe skin diseases.
- Women who are pregnant, breastfeeding, or are currently having their menstrual period. Hormonal fluctuations can make skin more sensitive and reactive to light therapy.
- Persons with skin injuries such as open wounds, irritation, swelling, severe sunburns, etc.
- Persons experiencing overfatigue.
- Persons with capillary dilatation caused by long term use of steroid hormonal drugs or a disorder of liver function.
- Persons who had an operation on their eye(s)/face; who have a disease of the eye(s)/face; those who are currently undergoing treatment for a disease of the eye(s)/face.
- Persons taking medication that induces drowsiness.

2 / Avoid direct eye contact with the device due to the LED light therapy. Please protect your eyes by keeping them closed throughout the process or using protective eye gear such as goggles (not included).

3 / Use the product with caution. Do not hold the device while powered on one point of the face for an extended period of time to avoid skin damage or bruising. Ensure that the device is in constant motion across the face throughout use.

4 / Keep the product away from children.

5 / Always unplug the device once completely charged. Do not leave the device unsupervised.

6 / Keep the device away from wet surfaces. Do not use near a bathtub, sink, or any household furniture containing water.

7 / If the item has fallen into water or other liquids, do not reach for the machine. Unplug immediately.

8 / Keep the product and the power cord away from hot surfaces.

9 / Clean the device regularly to avoid the accumulation of grease and other residues.

## **CHOOSING THE RIGHT MODE:**

### **ALL-IN-ONE MODE (EMS + RF + LIGHT THERAPY)**

Enhance your facial toning and sculpting treatment with the powerful combination of EMS and RF, further improved by incorporating light therapy to address particular skin issues like acne, aging, and dullness, among others. Maximize the benefits of every treatment session.

### **TONE MODE (EMS + RF)**

Experience advanced and effective facial toning with the combination of EMS and RF. Elevate your skincare routine with this next-level treatment, which helps increase facial muscle vitality, eliminate puffiness, improve firmness, and restore skin elasticity. Additionally, RF helps to firm, tighten, and lift the skin while stimulating collagen production, providing even more benefits for your skin.

### **TARGET MODE (LIGHT THERAPY)**

Five different LED light styles to target specific skin concerns such as acne, redness, aging skin, etc. and help boost skincare absorption and overall treatment effectiveness.

## LED LIGHT THERAPY

- 🌈 Mixed Light - Enjoy a comprehensive rejuvenation treatment by using a mixture of three therapy lights for healthier, smoother skin.
- 🔵 Blue Light - Fight acne-causing bacteria and reduce oil production that leads to clogged pores with blue light therapy.
- 🟡 Yellow Light - Helps combat redness and irritation while boosting the skin's blood circulation. Yellow light may also be used to treat sunburn and radiation damage from ultraviolet (UV) rays.
- 🔴 Red Light - Improves your skin's appearance by stimulating collagen production and rejuvenating skin cells to treat wrinkles, redness, acne, scars, and other signs of aging to produce smooth, supple skin.
- 🌟 Flashing Mixed Light - For use after skincare application to enhance the absorption of products and boost their benefits.

## LET'S BEGIN

### LIGHT SENSITIVITY TEST

While light therapy is generally safe for most users because it does not contain ultraviolet (UV) light, we recommend a light sensitivity test or check with your physician before treatment.

1 / Power Lily up by pressing the (⦿) button to activate the light therapy feature.

2 / Test the device on a small patch or area of your skin (e.g. on your forearm) for no more than three minutes.

3 / Turn the device off and check the treated area to see if your skin has turned red. If the redness does not go away after 2 hours, your skin may be sensitive to the light.

4 / If your skin is sensitive to the light, discontinue the use of Lily and contact our Customer Service (note: Fancii stands by a 30-day return/ refund policy for orders placed via fancii.com).

### USING LILY

Mode	Function	Time	Frequency
<b>All-In-One</b>	Electric Muscle Stimulation (EMS) Radio Frequency (RF) LED Light Therapy	Approx. 5 mins per session	2-3 times a week
<b>Tone</b>	Electric Muscle Stimulation (EMS) Radio Frequency (RF)	Approx. 4 mins per treatment area	2-3 times a week
<b>Target</b>	Mixed Light Blue Light Yellow Light Red Light Flashing Mixed Light	Approx. 3 mins per treatment area	3-5 times a week



1 / We recommend charging Lily when you first receive it to ensure the battery is fully powered. Follow the instructions in the **Stay Powered** section to read how.

2 / Ensure your skin is clean, dry, and free of makeup.

3 / Choose and apply only skincare serum or essence evenly on your face.

## ALL-IN-ONE MODE

1 / Press the (♣) button to turn the device on at level 1 intensity. The device has 4 levels of intensity. Press the button again until you've reached your desired level.

2 / Press the (♣) button to turn on the LED light therapy option. Choose starting from mixed light, blue, yellow, red, and mixed flashing light. Press the button again until you've reached your desired light setting.

3 / Hold the device head flat against your face and gently massage in a circular motion starting from the forehead down to the neck.

4 / Turn the device off by pressing the buttons for 2-3 seconds. Alternatively, the device will automatically turn off after 10 minutes of continuous usage.

## **TONE MODE**

1 / Press the (♣) button to turn the device on at level 1 intensity. The device has 4 levels of intensity. Press the button again until you've reached your desired level.

2 / Hold the device head flat against your face and gently massage in a circular motion starting from the forehead down to the neck.

3 / Turn the device off by pressing the buttons for 2-3 seconds. Alternatively, the device will automatically turn off after 10 minutes of continuous usage.




## **TARGET MODE**

1 / Press the (♣) button to turn on the LED light therapy option. Choose starting from mixed light, blue, yellow, red, and mixed flashing light. Press the button again until you've reached your desired light setting.

2 / Hold the device head flat against your face and gently massage in a circular motion starting from the forehead down to the neck.

3 / Turn the device off by pressing the buttons for 2-3 seconds. Alternatively, the device will automatically turn off after 10 minutes of continuous usage.

## MASSAGE AREAS



## **STAY POWERED**

1 / Use the USB-C cord to charge the battery.

2 / Attach cord to USB-C port on the bottom of Lily and connect the other end to a power source. The battery takes approximately 4 hours to fully charge.

3 / When charging, the light on the right side of Lily will turn red. When fully charged, it will turn green.

## **CLEANING, STORAGE AND MAINTENANCE**

After use, unplug the device and store it in a dry, safe place. Use a clean, dry cloth to wipe down the device. Wipe the device head with a clean, damp cloth before storing away.

Ensure that the device is never exposed to direct sunlight or submerged underwater.

## **BATTERY DISPOSAL**

Note: Opening and disassembling the device will void its warranty. The process is not reversible. This action must only be undertaken when the device is ready to be disposed of. All batteries should be recycled or disposed of as per state or local guidelines, and should not be thrown away in household waste. To remove the battery, you need to disassemble the unit and carefully extract the lithium battery. Wear gloves during this process for your safety. Do not dispose of it in a fire as the battery will likely explode and cause injury.

## Here To Help

We pride ourselves on providing premium beauty tools which is why all of our products are backed by a 24-month warranty. If you have any questions or need further assistance, you can send us an email, connect with us on Live Chat or give us a call! Our friendly support team is here to help.



**1 · 800 · 285 · 0183**



**support@fancii.com**



LIVE ONLINE CHAT AT:

**fancii.com**

fancii