

SGS


Clinical study to determine the efficacy of 5-in-1 beauty device for wrinkle reduction, firming, elasticity and improving dermal density.

1 month study on 33 healthy volunteers using one cosmetic device

SUMMARY DOCUMENT

MAGNITONE

FACEROCKET CLINICAL STUDY



Summary

SGS

- This study was designed to assess the impact on the appearance of facial skin with regular use of a 5-in-1 beauty device (MAGNITONE FaceRocket) over 4 weeks.
- Instrumental and imaging-based measurements were made, together with volunteers' self-assessment of the efficacy attributes of the device.
- 33 volunteers aged 33 to 50 were recruited. The volunteers had visible signs of ageing on the face, including lines and wrinkles and skin laxity.
- Assessments were made after a single use, after 14 days and after 28 days of regular use.

MAGNITONE

FACEROCKET CLINICAL STUDY



Summary of results



With regular use, after 28 days the device is clinically proven to improve:

- Skin firmness by 29%
- Dermal density by 35%
- Wrinkle reduction by 10%
- Skin elasticity by 7%

Average result based on 33 participants

- 97% said their skin looked more lifted + toned
- 97% said their skin looked more radiant
- 100% said their skin looked more youthful.



MAGNITONE

FACEROCKET CLINICAL STUDY


Results

Skin Elasticity - Cutometer®


A Cutometer® ‘sucks’ the skin, ‘pinching’ it to measure its elasticity. The higher the measured value, the better the skin elasticity.

There was a decrease in skin elasticity, indicating a firming and tightening of the skin. With regular use over 4 weeks, this firming and tightening became statistically significant ($P < 0.05$).

Skin elasticity increased by an average of 7.24% after 28 days.



Results



Skin Firmness - Cutometer®

A Cutometer® ‘sucks’ the skin, ‘pinching’ it to measure its elasticity. The smaller the measured value, the better the skin firmness.

There was a decrease in skin elasticity, indicating a firming and tightening of the skin. With regular use over 4 weeks, this firming and tightening became statistically significant ($P < 0.05$).

Skin firmness increased by an average of 29.2% after 28 days.

Results


Derma Density - Ultrasound

High-frequency ultrasound allows an objective assessment of the cutaneous structure: the thickness of the epidermis and dermis as well as the dermal density, which are useful parameters for the quantification of the cutaneous regeneration process*. The higher the measured value, the greater the skin density.


With regular use over 4 weeks, the increase in dermal density became statistically significant ($P < 0.05$).

Derma density increased by an average of 35.1% after 28 days.

*Crisan D, Crisan M, Moldovan M, Lupsor M, Badea R. Ultrasonographic assessment of the cutaneous changes induced by topical flavonoid therapy. Clin Cosmet Investig Dermatol. 2012;5:7-13. doi: 10.2147/CCID.S25840. Epub 2012 Jan 10. PMID: 22291475; PMCID: PMC3267407.



Results







Wrinkle Volume - Canfieldsci Optical 3D Skin Measuring

Canfieldsci is an image-based skin analysis system. Both left and right-side images were analysed to measure wrinkle volume. The smaller the measured value, the fewer the wrinkles.

There was a decrease in wrinkle volume at the nasolabial fold. With regular use over 4 weeks, this became statistically significant ($P < 0.05$).

Wrinkle volume decreased by an average of 10.7% after 28 days.

Results





Evaluation Parameter	Time Point	Evaluation Sample		Parameter Remarks
		Rate of change	Difference Analysis (Comparing with D0)	
Skin Elasticity Value R2	After application for 14 days	7.59%	Significant difference	The increase of measured value means the increase of skin elasticity level
	After application for 28 days	7.24%	Significant difference	
Skin Firmness Value F4	After application for 14 days	-3.52%	No significant difference	The decrease of measured value means the increase of firmness level
	After application for 28 days	-29.20%	Significant difference	
Dermal Density	After application for 14 days	18.16%	Significant difference	The increase of measured value means the increase of skin density
	After application for 28 days	35.17%	Significant difference	
Wrinkles Volume	After application for 14 days	-10.07%	Significant difference	The decrease of measured value means the decrease of wrinkles volume
	After application for 28 days	-10.07%	Significant difference	

Examples of subjects

Dermal Density - Ultrasound

Wrinkle Volume - Canfieldsci Optical 3D Skin Measuring



Participant survey



Evaluation Item	Evaluation Standard
1. Are you satisfied with the cleanliness of facial skin after using this product?	<p>Evaluation criteria of Q1-Q11: 1 – Very Satisfied, 2-Satisfied, 3– Relatively Satisfied, 4 - Less Satisfied, 5-Dissatisfied, 6-Very Dissatisfied.</p> <p>Statistical figures <4 is the valid number, Satisfaction=Valid number/Total number×100%.</p> <p>Evaluation criteria of Q12-Q13: Yes or No.</p> <p>Evaluation criteria of Q14-Q17 see 4.11 subjects' self- assessment for details</p> <p>Evaluation time point:D28</p>
2. Are you satisfied with the effect of promoting the absorption of essence and improving the moisture after using this product?	
3. After using this product, are you satisfied with the effect of lifting and tightening the facial skin and eliminating edema?	
4. Are you satisfied with the brightening effect of facial skin after using this product?	
5. Are you satisfied with the smooth and tender effect of facial skin after using this product?	
6. Are you satisfied with the effect of reducing dry lines and fine lines of eye skin after using this product?	
7. Are you satisfied with the effect of fading facial wrinkles after using this product?	
8. Are you satisfied with the effect of moisturizing and strengthening skin barrier after using this product?	
9. Are you satisfied with the anti-aging effect after using this product?	
10. After using this product, are you satisfied with the effect of cooling and soothing facial skin and repairing after sunburn?	
11. Are you satisfied with the appearance design and the comfort of this product?	
12. Do you have any skin allergy after using this product?	
13. When you use this product, is the electric current of EMS micro-current acceptable?	
14. Are you satisfied with the effect of fading dark circles under your eyes after using this product?	
15. Are you satisfied with the effect of relieving acne inflammation of facial skin after using this product?	
16. What is the most commonly used mode of this product?	
17. Would you like to recommend this product to your friends?	

Subjects were given an evaluation form at 14 and 28 days to assess their experience using the device. A variety of evaluation criteria were assessed to determine the results.

After 28 days, key takeaways include:

- 97% said their skin looked more lifted + toned
- 97% said their skin looked more radiant
- 100% said their skin looked more youthful.

