

# We Equip the World to Heal.









# **BuffEnuff® Power Massager**

BE Fit, BE Flexible, BE Ready

# **User Manual**

# **BuffEnuff® Power Massager User Manual**

#### **General Information**

The **BuffEnuff®** Power Massager is a portable cordless handheld percussion massager. It features convenient hand held grips with a soft application surface used for pain relief, recovery from exercise, and general well being. The **BuffEnuff®** sends gentle percussive healing massage to tissue and muscles for relaxation and pain relief. The power for the device is controlled by the easy grip trigger button on the main handle and the application pad is applied using desired technique. The intensity, duration, and technique are applied according to your personal needs or the requirements of the client receiving application care.

#### Included in this BuffEnuff® Kit

#### \*BuffEnuff® Power Massager

- \*1 14.4 V Lithium Ion Battery
- \*1 Wall Charging Unit
- \*2 Plush Crowns
- \*BuffEnuff® Carry Bag

#### **Features**

#### BuffEnuff® Power Massager (For a schematic of the BuffEnuff® Power Massager visit www.percussionfitness.com)

- 1) 3300 BPM Motor
- 2) Cordless, Rechargeable and Portable
- 3) Dimensions: Height 12 " / Depth 7.5" / Width 6"
- 4) Weight: 3.20 lbs.
- 5) Grip: The **BuffEnuff**® has many handgrip options, more than other massagers.
- 6) Foam Disk Pad: 6" wide, 1" thick durable and soft foam.

#### 1.1 Overview

Do you suffer from Muscle Soreness, Pain, or Stress?

The ultimate healing device, the **BuffEnuff®** "**BE" Power Massager** delivers powerful percussive therapy to tight muscles, trigger points and scar tissue and can quickly relax tense muscles, provide instant relief from tension, and relieve everyday aches and chronic pain.

Pain is your body's attempt to alert you there could be a problem. Pain is important because it helps you stay alert to changes within your body. Ignoring pain may cause damage or injury. Sometimes pain is prolonged and requires medical attention to determine the cause and possible ways to help alleviate it.

The perception of pain does not begin until the messages traveling from the body to the brain are received and decoded. The pain messages travel from an injured area (both chronic and acute) along small nerve pathways and into the spinal cord to the brain. In the brain, the messages are detected and a signal is sent to the original area where the first messages arose. Your awareness of pain is felt immediately in some cases and delayed in other instances. Pain is neither good nor bad but mainly a signal for you to interpret and provide appropriate response.

#### **Explanation of Percussion Massage**

Percussion Massage is a non-invasive drug free method for controlling pain and improving function. The **BuffEnuff® Power Massager** uses gentle, yet powerful, waves of mechanical motion to gently assist muscles to improve blood flow, lymphatic activity, and flexibility.

In some cases, the physiological actions created by application with the **BuffEnuff® Power Massager** can lead to pain relief. Percussion Massage does not cure any physiological problem. Instead, it enhances your own body's recuperative capacity to heal itself. Individual responses will vary depending on condition, age, genetics, activity level, etc. However, many people experience benefits from regular use of the **BuffEnuff® Power Massager's** percussive care application.

### 1.2 Indications for use

The **BuffEnuff® Power Massager** may be used for the following conditions (these statements have not been evaluated by the food and drug administration.):

- 1) Symptomatic relief of burdensome chronic pain.
- 2) Recovery from exercise.
- Prevention of exercise and movement related injury.
- 4) Relaxation of hypertonic musculature

- of morecome blood now onodiation.
- 6) Treating myofascial dystonia.
- 7) Post traumatic pain.
- 8) Maintaining or increasing range of motion.
- 9) Post-surgical pain.
- 10) Scar tissue treatment.
- 11) General tension and anxiety management.
- 12) Prevention of disuse atrophy.
- 13) Integrated wellness program.

#### Benefits of Use

- Great for Self-Care. Clinical Use and Rehabilitation
- Peak Performance
- Treating Scar Tissue
- Treating Trigger Points and Sore Muscles
- Ask us about our Affiliate Program!
- Warm-up, Pre-Workout and Boost Post-Workout Recovery

#### IMPORTANT SAFETY INFORMATION

Read the instruction manual before operating the **BuffEnuff® Power Massager**. Failure to comply with all "Contraindications", "Warnings", "Cautions", and "Adverse Reactions" can cause harm to user or device.

#### 1.4 GENERAL SAFETY RULES FOR ALL BATTERY OPERATED TOOLS

**WARNING!** When using battery operated tools, basic safety precautions should always be followed to reduce the risk of fire, leaking batteries and personal injury, including the following: Read all these instructions before operating this product, and save these instructions.

- Disposal of battery: ensure the battery is disposed of safely; it is advised to contact a lithium ion battery recycling agent.
- **Do not abuse the supply cord.** Never carry the tool by the cord or yank it to disconnect from the socket. Keep the cord away from heat, oil, and sharp edges.
- Keep your work area well lit. Cluttered shelves and dark areas invite accidents.
- Do not operate power tools in explosive atmospheres, such as in the presence of flammable liquids, gases, or dust. Power tools create sparks, which may ignite dust and fumes.
- Keep unnecessary bystanders, children, and pets while operating a power tool. Distractions may cause you to lose control.
- A battery operated tool with integral batteries or a separate battery pack must be recharged only with the specified charger for the battery. A charger, which may be suitable for one type of battery, may create a risk of fire if used with another battery.
- Use battery operated tool only with specified battery pack. Use of any other batteries may create risk of fire.
- WARNING! NEVER allow any parts of the battery pack or power supply to come into contact with water and other liquids, oils or lotions. Exposure to water or other liquids, oils or lotions can result in electric shock, burns, or fire. Exposure to water and other liquids, oils, or lotions can result in the corrosion of the battery pack, battery leakage, electric shock, burns, or fire.

#### 1.5 GENERAL SAFETY RULES (cont.)

#### **Personal Safety**

- Stay alert, watch what you are doing, and use common sense when operating a power tool. Do not operate the tool while tired or under the influence of drugs, alcohol, or medications. A moment of inattention while operating power tools may result in serious bodily injury or property damage.
- Dress properly. Do not wear loose clothing or jewelry. Keep hair, clothing, jewelry, fabric and gloves away from moving parts. Loose hair, clothing, jewelry, fabric, and gloves can become caught in the moving parts.
- Avoid accidental starting. Make sure switch is in the locked or OFF position when inserting battery pack.
- Operate power tool with clean hands only. Avoid gripping tool with oily hands. Oily hands may result in loss of grip and loss of control of the power tool.
- Do not overreach. Keep proper footing and balance at all times. Proper footing and balance enables better control of the tool in unexpected situations.

### 1.6 GENERAL SAFETY RULES (cont.)

#### **Tool Use and Care**

- **Do not force tool.** Using excessive pressure load may result in injury or tool malfunction.
- Do not use tool if switch does not turn it on or off. Any tool that cannot be controlled with the switch is dangerous and must be repaired or replaced.
- Turn the power switch to the OFF position and disconnect the battery pack before making any adjustments, changing accessories, or storing the tool. Such preventative safety measures reduce the risk of starting the tool accidentally.
- Store idle tools out of the reach of children or untrained persons. Tools are dangerous in the hands of untrained users.
- When battery pack is not in use, keep it away from other metal objects like; paperclips, coins, screws, metal tools, or other small metal objects that can make a connection from one terminal to another. Shorting the battery terminals together may cause, sparks, burns, or a fire.
- Maintain tools with care. Properly maintained tools are less likely to bind, are easier to control, and have a longer lifespan.
- Check for misalignment, or binding of moving parts, breakage of parts, and any other condition that will affect the tool's operation. If damaged, have the tool serviced before using. Poorly maintained tools cause many accidents.

#### 1.7 GENERAL SAFETY RULES (cont.)

- Use only accessories that are recommended for use by the manufacturer for your model. Accessories that are appropriate for one tool may become hazardous when used with another tool.
- Do not attempt to service this power tool. If this tool is broken or challenged in any way, discontinue use and consult the warranty information. Servicing this tool with foreign parts may result in fire or injury.
- DO not attempt to polish, wax or sand objects with this power tool. This tool is intended for massage only. Although the construct and appearance of this tool is similar to other power tools used for manual jobs, this tool should only be used for its intended purpose: massage

parts. They contain chemicals that can damage, weaken, or destroy plastic.

#### 1.8 BATTERY AND CHARGER SAFETY RULES

- 1) Charge only 14.4V battery pack with your charger supplied in this kit, other types may burst causing damage or personal injury.
- 2) Do not charge battery when temperature is below 10°@ (50°F) or above 40°@ (104°F).
- 3) Do not attempt to use a step-up transformer, an engine generator or DC power receptacle.
- 4) Do not short the battery pack; a battery short can cause a large current flow, overheating, possible burn, and even a fire.
- 5) Do not allow anything to cover or clog the charger vents.
- 6) Do not touch the terminals with any conductive material.
- 7) Do not store or carry battery pack in a container, bag or pocket with other metal objects, such as nails, coins, paper clips, fingernail clippers, keys, etc.
- 8) Do not expose battery pack to water or rain.
- 9) Do not store the tool in a location where the temperature may exceed 50°@ (122°F).
- 10) Do not incinerate the battery pack even if it is damaged or no longer holds a charge. The battery pack can explode in a fire.

#### 1.9 Rechargeable Battery Power Tool Safety Rules

- 1. Battery powered tools are generally in an operative condition. Be aware of possible hazards. Before using charger and tool, read all instructions and markings on tool, charger and battery pack. Before connecting the charger to the battery pack (plugging into receptacle or outlet), be sure voltage supplied is same as power requirements specified on charger. A voltage source greater than specified can **result in serious injury to the user**. When in doubt, **do not plug in the charger**.
- 2. **Do not** charge the battery pack in wet or damp environment or expose to rain or snow.
- 3. **Do not** operate charger or tool if damaged in any way. **Do not** operate charger with damaged cord. Return product to factory. Incorrect assembly may result in electrical shock or fire. To prevent damage to cord when disconnecting from power source, grip plug body when removing it from a receptacle. Make sure cord is not subjected to unnecessary abuse.
- 4. Unplug charger from outlet when attempting maintenance or cleaning.
- 5. **Do not** substitute charger of any other device or use an attachment not recommended by manufacturer. Doing so may result in risk of electrical shock, fire, or risk to persons.
- 6. If batteries leak, avoid contact with the liquid. Should contact occur, wash immediately with soap and water and seek medical attention.
- 7. **Do not** alter or modify tool. Doing so may result in dangerous conditions.
- 8. Keep charger free of foreign matter.
- **9. Do not** run tool while carrying it. Operate in a stable, stationary position.
- 10. The use of any other accessories not specified in this manual may create a hazard.
- 11. Disassembly or tampering with housing of device is not advised and will void warranty.
- 12. Do not use force when inserting battery pack. To insert fully, make sure both latches on battery pack have snapped into place.

#### 2.0 Contraindications

- 1) This device should not be used to treat symptomatic local pain unless causes and conditions of the pain are established and a pain syndrome has been diagnosed by a qualified healthcare professional.
- 2) This device should not be used when cancerous tumors are present in the treatment area.
- 3) **Percussion massage should not** be applied over skin eruptions (e.g. varicose veins, thrombophlebitis, phlebitis, etc.), swollen areas, inflamed tissue, infections, or injuries such as sprains, broken bones, bruises and burst blood vessels.
- 4) **Do not** use this device over the head, face and neck regions. With special training and techniques some minimal application of the device may be appropriate for these regions. Consult a trained technician or licensed health care provider for guidance.
- 5) **Do not** use this device if patient has a cardiac pacemaker or any implanted defibrillator.
- 6) Patients with sepsis or any type of bacterial toxemia should not use this device.
- 7) This device should not be used over poorly enervated areas.
- 8) **Do not** use this device if patient has blood clotting, serious arterial circulatory problems in the lower limbs, or is currently taking blood thinning medication.
- 9) This device should be not be used by people with epilepsy.
- 10) Patients with abdominal or inguinal hernia should not use this device.
- 11) **Do not** use this device if you have heart disease without consulting your physician.
- 12) Do not use this device if you have kidney disease without consulting your physician.

### Warnings, Cautions, and Adverse Reactions

# 2.1 WARNINGS

- 1) This device is optimally used under the continued supervision of a licensed physician or certified practitioner.
- 2) The long-term effects of percussive massager stimulation are unknown.
- 3) Percussion Power Massager devices do not have any curative value.
- 4) Percussion massage is a symptomatic treatment, and, as such, suppresses the sensation of pain, which would otherwise serve as a protective mechanism.
- 5) Safety has not been established for the use of percussive massage during pregnancy. Do not use during pregnancy unless directed by your physician.
- 6) Percussion massage is not effective for pain from central origin, such as headache.
- 7) Electronic monitoring equipment (such as ECG equipment, and ECG alarms) may not operate correctly when percussion massager stimulation is in use.
- 8) Never use the percussion massager in wet conditions such as bath, shower, pool, hot tub, lake, pond, ocean, etc.
- 9) Caution should be used in applying percussion massager stimulation to patients suffering from heart disease. Further clinical data is needed to show there are no adverse results.
- 10) Never use percussion massagers near the heart. Caution should be used when applying the percussion massager pad anywhere on the front of the thorax (noted by the breastbone and ribs). Take extreme caution when placing on or near the two large pectoral muscles. Here it can increase the risk of ventrical fibrillation and lead to cardiac arrest.
- 11) Percussion massager stimulation should never be placed over the ears, nose, eyes, throat, mouth, near the genitals, or internally.
- 12) Never use percussion massager stimulation on areas of skin that lack normal sensation.
- 13) Apply percussion massager pad only to clean, unbroken, healthy skin.

- 17/ 100p the perodocion modelyer out of redoit of children
- 15) Percussion massagers are not made for use on animals or pets.
- 16) Consult your doctor if you have any doubt about using the percussion massager whatsoever.

#### 2.2 CAUTIONS

- 1) Percussion massagers are for single patient use only.
- 2) Bacterial infections, viruses and substances such as dirt, oil, blood and secretions can be transferred from one person to another if percussion devices (particularly the cloth crowns that wrap around the disk) are shared between people.
- 3) Keep yourself informed of the contraindications.
- 4) This percussion massager is not intended for unattended, personal use by patients who are noncompliant, emotionally disturbed, mentally handicapped, suffering from dementia, or has mental faculty challenges.
- 5) Thoroughly read and understand the warnings, cautions and operating instructions before using the percussion massager device. Know the limitations and hazards associated with using the percussion massager device. Always follow the operating instructions provided by your healthcare provider.
- Instructions of use are listed; any improper uses may be dangerous.
- 7) Do not use a percussion massager device on undiagnosed pain syndromes until consulting a physician.
- 8) Patients with an implanted electronic device, such as a cardiac pacemaker, implanted defibrillator, or any other electronic device, should not operate the percussion massager without first consulting a doctor.
- 9) Patients with internal components such as plates, screws and mesh should use caution when operating the percussion massager over these areas.
- 10) Patients with blood clotting, blood thinning, blood pressure issues, heart disease, kidney disease, liver failure, epilepsy, cancer, or any other health condition should not use this device without first consulting a physician.
- 11) Some patients may experience skin irritation or hypersensitivity, possibly due to the frictional application of the percussion device and its dermabrasion action, or as a result of allergic reaction to the cloth bonnet material. If rash develops and pain persists, discontinue use and consult a doctor.
- 12) Percussion massager device should only be placed on healthy skin. Isolated cases of skin irritation may occur at the site of continual long-term percussion massager focus. If this occurs, discontinue application and consult a physician.
- 13) Skin irritation, bruising and soreness may occur if excessive pressure is applied. If this occurs, discontinue application and consult prescribing healthcare provider.
- 14) This device and associated symptomatic treatment applications have yet to be certified by the food and drug administration and users are warned to use at your own risk.
- 15) Percussion massager disk placement should be based on the educational information provided herein, with the guidance of a health care professional or prescribing practitioner.
- 16) If application of percussion massager becomes uncomfortable due to pressure or length of focus, reduce pressure or length of time applied to an area, and contact a doctor is the problem persists.
- 17) Effectiveness is variable and depends on many cofactors such as activity, health, age, proximity to injury, etc. A spectrum of outcomes from applying the percussion massager may occur according to these cofactors. Considering cofactors, always use best judgment before operating device
- 18) This device should not be used while driving, operating machinery, close to water, or any activity in which involuntary muscle contractions, distraction, and increased relaxed states may put the user at undue risk of injury.
- 19) Use caution if using chemical-based or essential oil based lotions, gels or oils in conjunction with the percussion massager. Burns or skin irritation may occur, and if so, consult a doctor.
- 20) The percussion massager has not been approved by regulatory agencies for use on the body. Though similar products that exist in the marketplace are used on the body and have demonstrated benefit and approval by the general public. Ongoing research is recommended to demonstrate a positive correlation between use of devices like this one and favorable outcomes. Use of this device is on a voluntary basis and user assumes all risk for application.
- 21) Long hair, loose clothing and fabrics, and jewelry can become tangled and caught in the spinning head of the percussion massager. This may result in hairloss, ruined clothing and jewelry, or device failure. Follow operating instructions to ensure this does not happen. If tangling happens, stay calm and immediately power off the device before carefully untangling.
- 22) Make sure the cloth crown is centered properly on the disk before operating the percussion massager. Turn the massager off before removal of the cloth crown.
- 23) The **BuffEnuff® Power Massager** should only be used with the rechargeable lithium-ion batteries and cloth crowns recommended for use by the manufacturer.

# 2.3 ADVERSE REACTIONS

- 1) Skin irritation and excessive painful dermabrasion are adverse reactions, which have been associated with prolonged application, excessively applied pressure, and rare allergic reactions to the cloth crown material. If skin irritation or muscular pain occurs, discontinue use and consult your physician.
- 2) If percussion massager stimulation levels are too intense, reduce applied pressure to a comfortable level, and consult a physician if problems persist.

### **INSTRUCTIONS FOR USE**

# 2.4 Applications for BuffEnuff® Power Massager

The **BuffEnuff® Healing Technique** is created to help you achieve the best results from your **BuffEnuff® Power Massager**. The techniques are simple to learn and teach others. No special training is required to use the **BuffEnuff® Power Massager** and it is easy to use on yourself too!

Additional advanced training is available and certification courses are provided for fitness and health professionals. Courses are fun and highly interactive and will provide health and fitness professionals with an in depth understanding and application of this powerful therapeutic healing massager.

DECORDING THE LAD CONTACT ON THE BODT

- **FLAT**® When the flat part of the pad surface is used in full contact with the body.
- TILT® When contacting the body with a slight tilt in the Buff Enough and the pad is at about a 30-45 degree angle to the body.
- EDGE® When contacting the body with the side of the pad and may by used for deeper penetrating healing.

#### THE APPLICATION OF BUFFENUFF® TECHNIQUE

<u>HOVER</u>® - Using a flat contact hold the **BuffEnuff**® softly over the body with light or no pressure applied and move gently. HOVER® helps patient/client get used to the percussion action of the **BuffEnuff® Power Massager**.

<u>WRAPPING</u>® - Using a flat contact rest the buffer gently on the body and apply light downward pressure while moving the **BuffEnuff®** around the body surfaces from one side to other. Imagine wrapping a large bow around the body as you trace a line back and forth carefully massaging all body surfaces along the back, buttocks, and legs. Avoid the head during this technique. WRAPPING® is an excellent technique to use following a HOVERING®.

<u>CLEARING</u>® - Using a tilt contact, apply a moderate amount of downward pressure toward the tilted edge of the **BuffEnuff**®. Then move the **BuffEnuff**® along the muscle surface up and down like you were ironing a garment. Using one pattern of movement (ie. top - down) is the most effective way to clear the muscle of waste products that can cause muscle pain. CLEARING® is helpful after exercise to reduce muscle soreness.

**STRIPING**® - Using an edge contact, imagine you're painting a stripe along the body working perpendicular to the muscle fibers. Make long and short motions across the body and apply pressure to tolerance of patient/client. Always get feedback about the pressure you're using on the **BuffEnuff**®. Note: Edge contact delivers more concentrated percussion action into muscle tissue.

<u>DRUMMING</u>® - Using a flat contact, hold the **BuffEnuff**® gently over the treatment area in one location. Now gently lift the **BuffEnuff**® 1-3 inches off the treated area and then gently tap the pad up and down like you're using the **BuffEnuff**® as a drumstick and the body is a drum. Gentle up and down movements with minimal force are all that is required.

**GROUNDING**® - Using a flat contact, hold the **BuffEnuff**® firmly over the treatment area and apply firm pressure. While holding the pressure firmly, lighten the pressure briefly without allowing the pad to lose contact with the body. Then apply deep firm pressure again. Repeat as needed.

**SQUARING**® - Using the **BuffEnuff's**® edge apply moderate pressure for 3-5 seconds. Then remove and apply in the same location except move the edge contact perpendicularly to so the healing waves are distributed at a different angle to the tissue.

<u>HINGING</u>® - Using a flat contact or a slight variation of the tilt or edge, hold the buffer over the joint that will receive the hinging technique. Maintain a moderate depth contact with the **BuffEnuff**® and instruct patient to move his/her joint through its normal range of motion. It is very important to maintain full contact with joint being HINGED® to maximize the therapeutic results.

<u>JUMPING</u>® - This technique is similar to DRUMMING®, except you will use the BuffEnuff® edge rather than the flat contact. Technique can be used anywhere but is most effective on tendon attachment points of muscles like the achilles tendon.

**SPARKING**® - Using a flat contact, hold buffer lightly and then vigorously move the buffer back and forth and in a circular motion like you're rubbing out a sore muscle in the body or applying lotion to the skin. The patient will feel additional generated warmth from the application of SPARKING® applied with the rotation of the pad head on the **BuffEnuff**®. This technique is great following SQUARING® for trigger point care.

# **BUFFENUFF® TECHNIQUE COPYRIGHT 2016**

# 2.5 Battery Operation

- Your new battery pack is not fully charged when received; therefore, charging is necessary prior to use.
- There are two conditions when battery packs will not accept a full charge: when it is a new battery pack, and when the battery pack is not used for an extended period of time. These are normal conditions and do not indicate a problem. However, after approximately 4 charge and discharge cycles, the battery will be able to accept a full charge.
- When not in use, the battery pack can be left on charge (indicated by the red light) without any adverse effect; the battery cannot be "overcharged." However, when the battery pack is not going to be used for an extended period of time, the charger should be unplugged from the power source, and the battery pack disconnected from the charger.
- If you need to charge multiple battery packs with the same charger, leave 15 minutes of rest time between charges. This will result in the charger lasting longer.

#### 2.6 Charging the Battery

battory. I requent use may require auditional onarged batteries for the user's convenience

- 1) The battery is located at the bottom of the handle.
- 2) To release the battery, squeeze the two grip tabs, located on either side of the battery, and gently pull battery away from the massager.
- 3) Plug the power supply cord into a wall outlet.
- 4) Insert the power supply cord adapter into the battery connection port located on the top of the battery near the connection stem.
- 5) When the battery is charging, a red light is displayed on the power supply.
- 6) Charging a dead battery requires approximately one hour of recharge time. Recharging times may vary.
- 7) When the lithium-ion battery is fully charged, the display will change to green.
- 8) When the battery is fully charged, disconnect from the power supply and insert battery into the power massager.
- 9) A fully charged battery will supply power to the massager for approximately 30-40 minutes of non-stop massaging. Times may vary depending on application, loads applied to the massager, and frequency of use.

#### 2.7 Crown

The cloth crown is the elastic fabric cap, which is centered around the foam disk on the massager. The crown is meant to separate the spinning foam disk of the massager from contacting the skin of the user/receiver. The crown is plush and soft.

- 1) Lotion Application: The crown is safe to be in contact with the skin, and can even be used with oils, lotions and gels.
- 2) Sanitation: Because the crown is the only actual part, which is in contact with the skin of the user/receiver, it is important to wash the crown periodically, as it will pick up dirt and oils. Crowns can be machine washed and dried. To prolong the lifespan of the crown, wash your crowns in the sink with hot water and soap. Hang to dry. Crowns may also be sanitized periodically using fluorescent lighting certified as cleaning agents to kill bacteria and other unsanitary material. This does not replace thorough cleaning in a washing machine or hand washing.
- 3) Extra Crowns: The crown is the only part of the BuffEnuff® that will need to be routinely replaced. Each crown will sustain many massage sessions. Frequent users may want to purchase a package of crowns for convenience.

#### **CAUTION**

In rare cases, users have experienced skin irritation, possibly due to allergic reactions to the fabric of the crown. If this happens, try using the massager over a layer of clothing, a sheet, or other fabric. Discontinue use completely if the problem persists.

#### 2.8 CLEANING AND CARE

#### Massager

- 1) Periodically wipe down with a moist antimicrobial wipe.
- 2) Extra cleaning is necessary when used by sweaty people in the gym, and when using massage lotions and gels.

#### Crown

- 1) Because the crown is the only actual part, which is in contact with the skin of the user/receiver, it is important to wash the crown periodically, as it will pick up dirt and oils.
- 2) Wash your crowns in the sink with hot water and soap. Hang to dry.
- B) Use certified fluorescent lights that have their primary function as cleaning agents to periodically sanitize in between washings.

#### 2.9 STORAGE

- 1) It is best to keep the massager, battery and power supply, crown, extra batteries, extra bonnets in your **BuffEnuff® Carry Bag** for convenience.
- 2) The **BuffEnuff® Power Massager** is small, portable, lightweight, and made to be taken to the gym, locker room, in the car, to work, or easily stored in a desk drawer.
- 3) Store in a dry space away from water and excessive heat.

#### 3.0 WARRANTY

Please contact your dealer in case of claim under the warranty. If you have to send the percussion massager to dealer, please enclose a copy of your receipt and state what defect(s) exist. The following warranty terms apply:

The warranty period for the **BuffEnuff® Power Massager** is one year from date of purchase. In case of warranty claim, the date of purchase must be proven by dated sales receipt or invoice.

- ) Replacements under the warranty period do not extend the warranty period for a returned massager.
- 2) Warranty is provided on a prorated basis based on length time massager has been in use from purchase date.
- 3) The following is excluded under the warranty:
  - all damage from improper use
  - all damage from mishandling, inadvertently dropping the unit, etc.
  - all damages from tampering with the mechanical functions of the unit
  - all damages resulting from transport in gym bags, suitcases, etc.
  - accessories which will experience normal wear and tear
  - liability for direct or indirect consequential losses caused by the percussion massager is excluded even if the damage to the unit is accepted as a warranty claim

www.percussionfitness.com

CALL (833) BUF-ENUF (833) 283-3683

For questions about our Affiliate Program or to order supplies.

**Quantity Discounts Available**