

#### Please Note:

- The total time required for treatment with the HairMax LaserComb Advanced 7 is 15 minutes.
- It is important to move the LaserComb very slowly. Remember you are bathing the scalp with laser light. If you move it too fast, the treatment will not achieve the optimal results to promote hair growth.
- While there are no negative effects from overuse, it is recommended that you use the LaserComb according to directions.
- The LaserComb has advanced electronics to make sure that the device is working properly. If for some reason you are experiencing difficulties with your device, before sending it back to us for repair, please refer to the Troubleshooting section of your User Manual. For further assistance please call +1-561-417-0200.

**HAIRMAX**<sup>®</sup>  
The science of hair growth.

## QUICK START GUIDE

### Advanced 7 LaserComb



#### REGISTER

To validate your warranty please register online at [www.hairmax.com/warranty](http://www.hairmax.com/warranty) or call Customer Service toll-free at 1-866-527-3726

#### CONTACT

Lexington International, LLC  
777 Yamato Road, Suite 105  
Boca Raton, FL 33431  
+1-561-417-0200  
[HAIRMAX.com](http://HAIRMAX.com)



The HairMax LaserComb® is a very effective device when used correctly. It is manufactured using precision components and should last 10-15 years under normal usage. **Please take time to read the entire User Manual before operating the HairMax LaserComb.**

## Advanced 7 LaserComb



Use the LaserComb for 15 minutes, 3 times a week. The treatments should be spaced every other day, for example Monday, Wednesday, Friday. Move the LaserComb slowly 1/2 inch every 4 seconds.

It is recommended but not essential that the LaserComb be used on a clean scalp after washing. Your hair can be wet, towel dried or dry. The LaserComb should be used prior to applying any topical treatments, gels, spray or concealing products.

	<p><b>A</b> Your Advanced 7 LaserComb comes with four AA batteries already installed. Prior to use, remove the cardboard strip located in the battery compartment (Figure A). For more detailed instructions on how to remove the cardboard strip and to replace the batteries please refer to your User Manual.</p>
	<p><b>B</b> The LaserComb can also be powered directly by the optional power cord (not included). Plug the power cord into the end of the LaserComb. Plug the other end into an electrical wall outlet (Figure B).</p>
	<p><b>C</b> To turn on the LaserComb slide the power/laser button to the middle or top position. Middle position activates the laser and vibrate feature. Top position is laser only (Figure C).</p>

	<p><b>D</b> Before beginning treatment, be sure to test the LaserComb for sensitivity (see page 5 of your User Manual). To start your treatment place the LaserComb flat on your scalp so that both rows of teeth are touching the scalp (Figure D).</p>
	<p><b>E</b> Start at your hairline. Move the LaserComb from front to back of scalp in a slow combing motion, stopping every 1/2 inch (1cm) and pausing for 4 seconds (Figure E). Do not lift the LaserComb off your scalp while combing from front to back.</p>
	<p><b>F</b> Next, lift the LaserComb off your scalp and start at a new position at the hair line. Continue in this manner from front to back and then side to side (Figure F) until you cover your entire scalp with multiple passes. When the LaserComb is on it will vibrate every four (4) seconds. This means that it is time to move the LaserComb to the next area of the scalp. Repeat these passes within the required treatment of 15 minutes.</p>
	<p><b>G</b> If you have long hair, use your fingers to lift your hair out of the way while using the LaserComb to be sure that the laser energy reaches your scalp (Figure G).</p>
	<p><b>H</b> After treatment, turn off the LaserComb by sliding the power/laser button to the lowest position. If using the optional power cord, remove the power cord from the electrical wall outlet (Figure H). Store the device in the cradle/stand for future use.</p>