


## IS HAIRMAX® RIGHT FOR YOU?

HairMax is most effective in women and men in mild to moderate stages of hair loss and thinning. If your hair loss looks like the images below, HairMax is recommended for you.

### FEMALE



### MALE



## DOCTOR RECOMMENDED

Physicians worldwide recommend HairMax Laser Light Therapy to their patients. HairMax is supported by a Medical Advisory Board of established leaders in the fields of dermatology and hair restoration.

## FREQUENTLY ASKED QUESTIONS

### How often do I need to use my Prima 7 LaserComb?

It is recommended to use your Prima 7 three times a week for 15 minutes, with at least a day between each treatment.

### Will results improve if I use my Prima 7 LaserComb more than 3 times a week?

Three times a week has been found to be the optimal treatment time. Using your Prima 7 more than 3 times a week has not shown any additional benefits.

### How long will I need to continue hair loss treatments with my Prima 7 LaserComb?

To keep seeing benefits use your Prima 7 as directed on an ongoing basis. If you do not use the device as recommended, there is a strong possibility that you will not see benefits or your hair loss will regress to where you were prior to using the Prima 7.

### Can I use my Prima 7 LaserComb on wet or dry hair?

You can use your device on dry or damp hair. It is important that your scalp be clean and free of any styling products. A build-up of topicals such as oils, gels, sprays, etc. can block the laser light from effectively reaching the hair follicles.

## HAIRMAX® THINNING HAIR CARE

From roots to ends HairMax brings new life to your thinning hair with the den•sity collection, a full line of specially formulated hair care products, including shampoo, conditioner, scalp treatments, supplements and hair fibers.



Lexington International, LLC  
1040 Holland Dr, Boca Raton, FL 33487  
Tel. +1.561.417.0200 | U.S. Toll Free 1.800.9REGROW  
FAX. +1.561.892.0747

[www.HairMax.com](http://www.HairMax.com)

LIT-CONBROLEX-PRM7-02 01/17 © Copyright 2017 Lexington Intl., LLC. All rights reserved.



## THE ULTIMATE SOLUTION FOR HAIR LOSS


HairMax® Prima 7 LaserComb



**HAIRMAX®**  
The science of hair growth.™

## HOW LASER LIGHT THERAPY PROMOTES HAIR GROWTH

HairMax® delivers therapeutic, nourishing light energy to your hair follicles to activate them for stronger growth. This boost of light energy revitalizes and enlivens weakened hair follicles to stimulate growth factors and reverse the thinning process. This will help extend and restore your hair's natural growth cycle, resulting in hair growth with increased density, fullness and vibrancy.



LASER PHOTOTHERAPY FOR HAIR GROWTH

## 3 PART SYNERGISTIC SYSTEM

The intensity of the laser light, the density of the laser light and the hair parting teeth work in unison to provide a highly efficacious hair loss treatment.

*"I am writing to sing the praises of the HairMax LaserComb. I first learned of your product a few years ago as I began experiencing female thinning hair. The HairMax LaserComb has been a godsend, as it has kept my hair shiny, helped with new hair growth and appears to have even thickened it. All I can say is thanks for making such a great product available and affordable for personal use." - Ruth Stone*



## SEEING IS BELIEVING

RESTORE, REVITALIZE AND REGROW YOUR HAIR!



**Female:** Significant new hair growth in the temple areas.



**Male:** Significant new hair growth in the frontal areas of the scalp.

## CLINICALLY PROVEN

HairMax has been proven safe and effective in 7 clinical studies with 460 participants. Men and women in these studies experienced significant hair growth with an average increase in hair count of 129 additional hairs per square inch. HairMax clinical study results have been published in 6 peer review medical journals.

# The Ultimate Solution For Hair Loss & Thinning Hair

CLINICALLY TESTED

FDA CLEARED\*

DOCTOR RECOMMENDED

15 minutes

3 times per week

## HAIRMAX® PRIMA 7 LASERCOMB

This sophisticated laser device safely and effectively delivers nourishing laser light to energize and revitalize your hair follicles to promote hair growth. The patented hair parting teeth incorporated in the device parts your hair to ensure the light energy reaches your scalp.

Patented Hair Parting Teeth

7 Medical Grade Lasers

- Stimulates Hair Growth
- Treats Hair Loss
- Reverses the Thinning Process
- Increases Density and Fullness
- Revitalizes Damaged Hair
- No Harmful Side Effects
- Drug-Free Treatment

NOURISH. STIMULATE. GROW.

Over **90%** of participants in clinical studies experienced significant hair growth

\*HairMax laser devices are indicated to treat Androgenetic Alopecia, and promote hair growth in males who have Norwood Hamilton Classifications of IIa to V and in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I to IV. Laser Light - Avoid direct eye exposure.

