

# GAIAM WELLBEING™

## HAND RELIEF KIT

### Setup, Workout Tips and Care & Safety Guide

Please read before setting up or using your Hand Relief Kit.

#### **SETUP**

1. Remove Hand Relief Balls from box and check for shipping damage (scratches, punctures or gouges).
2. Your Hand Relief Balls are now ready for use.

#### **WORKOUT TIPS**

1. Follow the exercise instructions presented on the guide.
2. Ensure your hands are clean and dry before use.
3. Use a smooth and steady motion when performing each exercise.
4. Perform exercise routines to the best of your ability without strain.
5. Avoid holding your breath while exercising.

#### **CARE**

1. Dry Hand Relief Balls with a towel after use.
2. Avoid exposing the Hand Relief Balls to heat or excessive sunlight.
3. Store in a cool, dry place.

GAIAM  WELLBEING™

# HAND RELIEF KIT

## Setup, Workout Tips and Care & Safety Guide

Please read before setting up or using your Hand Relief Kit.

### IMPORTANT SAFETY CAUTIONS!

1. Not all exercise is suitable for everyone. To reduce the risk of injury, consult your doctor before beginning this or any exercise program or before using this or any other exercise equipment. The instruction presented herein is in no way intended as a substitute for medical counseling.
2. If you suffer from heart disease, high blood pressure or any other disease or condition, consult your doctor before beginning this workout.
3. Perform exercises in a slow and controlled manner.
4. Stop and rest if you feel dizzy or short of breath.
5. Use product only as intended and demonstrated on the exercise guide.
6. This product is sold without any warranties or guarantees of any kind and manufacturer and distributor disclaim any liability, loss or damage caused by its use.

Learn more at <http://life.gaiam.com>

©2016 Gaiam and the "Flower of Life" logo design are registered trademarks and Wellbeing is a trademark of Gaiam Americas, Inc. or its subsidiaries.  
Manufactured and distributed by Gaiam Americas, Inc., Louisville, CO 80027-2452.

MADE IN CHINA.

SKU 05-52108CARE