HOT & COLD THERAPY FOOT WRAP

Setup, Workout Tips and Care & Safety Guide

Please read before setting up or using your Hot & Cold Therapy Foot Wrap.

SETUP

1. Remove the Hot & Cold Therapy Foot Wrap from the box and check for possible shipping damage (scratches, punctures or gouges).
2. Your Hot & Cold Therapy Foot Wrap is now ready for use.

INSTRUCTIONS FOR USE

COLD THERAPY

1. Begin by taking out the gel core from the center back of the wrap.
2. Place the gel core in the freezer for 4–6 hours to freeze.
3. Remove the gel core from the freezer and slide back into place.
Product will stay cold for 15–20 minutes. Refreeze as necessary.

HEAT THERAPY

1. Begin by taking out the gel core from the center back of the wrap.
2. Place the gel core into 120-degree (very hot) tap water for 7–10 minutes; test before applying to skin.
3. Slide the gel core pack into the therapy foot wrap. The heat is retained for approximately 10–15 minutes. Reheat as necessary.

Do not use other heating methods—conventional or microwave ovens or stove tops, for example—as overheating may cause burns or damage the Hot & Cold Therapy Foot Wrap. We do not recommend placing the Hot & Cold Therapy Foot Wrap in boiling water, as this may damage the seals on the product and/or cause burns.

CARE

1. Dry the Hot & Cold Therapy Foot Wrap with a towel after use. Check for nicks or tears.
2. Spot-clean the Hot & Cold Therapy Foot Wrap with a damp cloth. Air-dry.
3. Avoid exposing the Hot & Cold Therapy Foot Wrap to sharp or abrasive surfaces, direct heat or excessive sunlight.

HOT & COLD THERAPY FOOT WRAP

Setup, Workout Tips and Care & Safety Guide

Please read before setting up or using your Hot & Cold Therapy Foot Wrap.

SETUP

IMPORTANT SAFETY CONSIDERATIONS

1. Self-treatment programs are not suitable for everyone. To reduce the risk of injury or complications, consult your doctor before use of the product.
2. This product is not intended to cure any medical conditions and is in no way intended as a substitute for medical counseling.

CAUTION: Use only as directed. Consult with your doctor or physical therapist before beginning any self-treatment program. Suitable for ages 13 and over. Keep out of reach of younger children. Do not use if you have diabetes, or if you have open sores or cracked or blistered skin in the treatment area. If symptoms worsen, decrease pressure; if adverse symptoms persist, discontinue use and contact your doctor.

VISIT GAIAM.CA/hotandcoldtherapyfootwrap FOR A DOWNLOADABLE EXERCISE GUIDE.

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