

TEI Spa OxyDerm

At-home Ozone Facial Tool for clear, clean blemish-free skin

KIT INCLUDES:

- OxyDerm Tool
 - OxyCalming Serum (30mL)
 - Travel Pouch
-

ABOUT OXYDERM HIGH FREQUENCY OZONE BEAUTY DEVICE:

The OxyDerm Professional Beauty Device has been developed by the professionals at TEI Spa to aid in the purification and sterilization of skin, fight against inflammation, and accelerate the skin's metabolism. TEI Spa professional technology has made it possible for those with blemish prone skin to experience maximum anti-acne results in the comfort of their own home, in just minutes a day. The OxyDerm not only helps clear up acne and blemished skin, but visibly brightens skin tone and complexion for health, radiant skin.



FEATURES & BENEFITS:

- Uses ozone to sterilize the skin and kill acne-causing bacteria
- Increases skin metabolism to healthy cell function and promote skin healing
- Reduce acne, scars and redness with the professional level at-home tool proven to give you clear, radiant skin

- Strengthens cell functions to improve the clarity and quality of your skin, giving it a more beautiful and youthful appearance



HIGH FREQUENCY TREATMENTS:

High frequency treatments are used to stimulate, sanitize and heal the skin. Used for years, this type of treatment is designed to provide high frequency currents that generate warmth in the skin and stimulate nerve endings, in turn increasing blood circulation, nourishing the skin and improving product penetration within the skin. By increasing skin metabolism, high frequency ozone facials encourage healthy cell function and promote skin healing. Direct high frequency offers the results of treating dry and dehydrated skin by increasing the activity in the sebaceous and sudoriferous glands and causing a skin warming effect. Direct high frequency ozone provides benefits for those with blemish prone and acneic type skin by killing bacteria and eliminating oils. In addition, high frequency facial treatments result in an improvement of the appearance of the skin as well as a reduction in fine lines and wrinkles.

HOW IT WORKS:

The OxyDerm At Home Ozone Facial Device works by generating tiny pulses of electricity, much like static electricity. When this static electricity travels through the air it ionizes the air it's passing through creating ozone (or three molecules of Oxygen) which possesses strong oxidation or sanitizing abilities. By expelling carbon dioxide into skin, your skin is able to absorb more oxygen, which allows better circulation. By strengthening cell functions you improve the clarity and quality of your skin, giving it a more beautiful, clear and youthful appearance.

DIRECTIONS:

Removal all metal Jewelry before use. Not for use by individuals with braces.

Be sure to use the OxyDerm only on clean and dry skin. It is recommended to use the OxyDerm nightly after removing all makeup, cleansing and drying the face.

For All Skin Types:

- For Oily Skin: Produce a drying effect that kills bacteria and diminishes excess sebum. Use daily, or as needed, on T-zone and U-zone in small circular motions for 1-2 minutes.
- For Combination Skin: Use 3-4 times per week, or as needed. Focus use on oily areas and use minimally on sensitive areas.
- For Dry Skin: Produce a warming effect that increases activity in the glands and optimizes skin nourishment. Use 2-3 times per week, or as needed, on T-zone in small circular motions for up to 1 minute.

For Full Body Use:

- Treat the Back & Body by using tool in large circular motions for 1-2 minutes. Ideal for treating acne, marks and redness on the body.
- Treat the Chest & Neck by using tool in large circular motions for 1-2 minutes. Ideal for treating redness, irritation and sensitivity on the body.
- Treat Specific Blemishes by using the tool in a small circular motion at a slower pace on affected area or blemish. Use for up to 30 seconds per target area.

OXYCALMING SERUM:

TEI Spa's unique OxyCalming Serum relaxes and moisturizes to soothe the skin after OxyDerm treatments. Natural, organic, botanical extracts, such as aloe, green tea and cucumber, aid in

the relief of irritation as they reduce redness and inflammation to instantly deliver lasting results and comfort. The OxyCalming serum is ideal for use on sensitive skin or for those who experience rosacea, burning or stinging.

The OxyCalming Serum was designed for use on the face and body in order to relieve any irritation, specifically after OxyDerm treatments.

Recommended Use: Apply 1-2 pumps of OxyCalming Serum on the hands and gently rub into treatment area. The OxyCalming Serum is ideal for use as a nightly skin serum on dry, irritated or burning skin. For best results do not wash off or wipe skin after use and allow serum to seep into the skin until dry.



INGREDIENTS: deionized water, apple fruit water, aloe barbadensis leaf juice, green tea extract, xanthan gum, oat extract, hyaluronic acid, cucumber fruit extract, allantoin, grapefruit seed extract, grape seed extract, thyme extract, thyme extract, apple extract, matricaria flower extract, azulene extract, pentylene glycol.

STORAGE & STERILIZATION:

If using the device on more than one person, unplug unit and gently sterilize the head with a disinfectant wipe or non-abrasive solvent after each use, being careful not to soak in water or alcohol.

If using the device on one person only, unplug unit and gently clean the head with a clean towel after each use.

Always secure protective cap after use and store in convenient travel pouch in a clean, dry place.

SPECIFICATIONS:

Power/Dial Control

I- Low Setting= 0.0003 PPM

II- Medium Setting =0.0004 PPM

III- High Setting= 0.0013 PPM

VOLTAGE INFORMATION:

Power Supply: 110V-120V/ 60Hz AC Adapter

Power: 1 Watt

Amperage: 0.001A

TOOL SIZE:

Tool Size: 8" H x 1.5" W

Cord Length: 72 inches

OZONE LIGHT:

Oscillation Rate: 100,000 Hz/Cycles per second

BEAUTIFUL SKIN BEGINS WITH A CLEANER FACE:

Traditional methods of acne prevention don't get you the results you'd really like. Fortunately, the lead inventor of TEI Spa discovered a breakthrough in how we sterilize our skin and rid of and prevent acne. Using High Frequency Technology, OxyDerm glides through the face while producing Ozone to kill bacteria that causes acne, redness and inflammation. Gently remove dirt, oil and shine while getting the radiant skin without harsh chemicals and drying cleansers.

- Removes more dirt, oil and bacteria without drying the skin or harming sensitive skin
- Works with the skin, not against it, to gently banish blemish causing bacteria stuck in the pores
- Prepares the skin for a vibrant and healthy glow
- Gentle enough to use 2x a day