

How to Use?

- Pull the insulating tape out of the battery compartment when using for the first time
- Press the start button to turn on the device and wait for 1 minute till the room temperature appears on the display
- Switch to °C/°F as your desired unit for temperature measurement
- Set the buzzer by pressing the ON/MEM button for 3 seconds

For Forehead measurement mode

- Pull off the probe lens cap to remove it
- Once the forehead icon appears, hold the unit with the probe lens 1 to 3 cm away from the center of the forehead and press "START" button. The measurement is completed in 1 second with a long beep
- Move the unit away from the forehead and check the measurement result
- Press and hold the ON/MEM button until "OFF" appears on the display to turn off the unit. The unit also turns off automatically if no operation is performed for 1 minute (the room temperature will still appear on the display)

For Memory function

- Press the ON/MEM button to turn on the unit
- Press the ON/MEM button again. The memory number appears on the display
- Release the ON/MEM button for the most recent result
- Press the ON/MEM button repeatedly to view the older results
- Press and hold the ON/MEM button until "OFF" appears on the display to turn off the unit (the room temperature will still appear on the display)

For Surface measurement mode

- Press the ON/MEM button to turn on the unit
- Press and hold the ON/MEM button and, press the START button simultaneously (ensure the surface icon appears when taking temperature)
- Move the unit/probe lens close to the object (1 cm apart is recommended) and press the START button. It measures continuously until the START button is released
- **P.S.-** Surface measurement mode is not intended for body temperature measurement or medical use
- Press and hold the ON/MEM button until "OFF" appears on the display to turn off the unit. The unit also turns off automatically if no operation is performed for 1 minute (the room temperature will still appear on the display)

REMEMBER

- Do not use the unit in presence of a strong static or electromagnetic field
- Do not apply the probe lens to the object directly
- Follow the instructions of your doctor as a high or prolonged fever needs medical attention
- For more information please refer to the product package list

Uploaded by www.heyzindagi.in