

How to Apply?

- Take an underpad out of the packaging and unfold it in half, lengthwise
- Help the patient move slightly to make space for laying the underpad on the bed
- Unfold the underpad entirely and spread it out evenly ensuring that there are no creases or wrinkles
- Remove the silicone release tapes and press the pad onto the fabric to ensure a light grip
- Help the patient to rest comfortably on the underpad

REMEMBER

- Please do not let your hands or other objects come in contact with the inner material of the underpad to avoid any accidental transfer of infection - for yourself/patient. Always use protective gloves during application and disposal

- Plan for disposal in advance - do not throw the underpads with regular garbage. They must be held in a separate container/bags for final dispatch with the garbage truck