## How to use Modo de empleo

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1.

- Sit/ lie down in a comfortable position.

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2.


- Place the Heating Pad on the affected body part requiring hot fomentation.

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4.

- Switch on the pad from hand controller.

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5.

- Choose the high, medium or low temperature mode by pressing the selector button.

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6.

- After use, Switch off, remove the pad and place it in the box.

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## Heating Pad Ortho

English

Heating Pad is a user friendly device to provide controlled hot fomentation. Hot fomentation of the injured/inflamed area is particularly helpful for sub acute and chronic injuries, it induces vaso dilation and tissue healing. The warmth decreases muscle spasm, relaxes tense muscles, relieves pain and can increase the range of motion. Along with pain relieving drugs, hot fomentation has a synergistic effect.

### CONSTRUCTION:

Made out of Teflon coated heating element, covered with a soft fleece fabric. The layered design ensures uniform & appropriate heating. Multiple safety features like timed & thermal cutoff are inbuilt for user safety. User friendly controls allow easy operation.

### APPLICATIONS:

- Muscular and Joint Pain,
- Backache
- Sprains and spasms
- Rheumatoid or other arthritis
- Geriatric care
- Sore or over worked muscles
- General comfort in extra cold conditions
- Abscess and Boils

### CONTENTS:

Polyester, Ni-chrome wire, PTFE, Cotton, PVC, Polyamide, Copper

### ▲ CAUTIONS:

- Disconnect main supply after use.
- Avoid sharp fold of the Pad and the main lead.
- Don't pull the main lead.
- Avoid contact with sharp edges.
- Do not sleep on the Heating Pad.
- Check for any mechanical damage before use.
- Don't use for infant or mentally challenged
- Don't use on persons with sensitive skin or ischemic limbs.
- Monitor temperature at regular intervals.

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