tynor

Ortho Heating Gel Bag (Electric)

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English

INTRODUCTION:

Ortho Heating Gel Bag (electric) is a user friendly, portable device to provide effective heat therapy for long duration. The therapy of the injured and painful area is particularly helpful for sub-acute and chronic injuries, it induces vaso dilation and tissue healing. The warmth decreases muscle spasm, relaxes tense muscles, relieves pain and can increase the range of motion. Along with pain relieving drugs, heat therapy has a synergistic effect.

Quick electric charging ensures product is ready to use in minutes. Cordless use ensures safety and unrestricted access.

CONSTRUCTION:

The gel pouch is made from three layered PVC, coated with soft feather like fabric. Multiple safety features like thermal cutoff, thermal fuse and LED indicator ensure better life & safety.

APPLICATIONS:

- Muscular and Joint Pain
- Backache
- Sprains and spasms
- Rheumatoid or other arthritis
- Geriatric care
- Sore or over worked muscles
- General comfort in extra cold conditions
- Abscess and Boils

FEATURES:

- Long Heat retention time
- Portable
- Auto-cutoff LED indicator
- Multi-layered insulation for safety and proper distribution of heat.
- Soft and flexible

▲ CAUTION:

- Don't use product while charging
- Disconnect the electric lead before use
- Avoid sharp fold of the Pad and the main lead
- Don't pull the main lead
- Avoid contact with sharp edges
- Do not sleep on the Heating Pad
- Check for any mechanical damage before use
- Don't use for infant or mentally challenged
- Don't use on persons with sensitive skin or ischemic limbs.

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English

POLICY FOR DEFECTIVE GOODS:

All products manufactured by Tynor Orthotics (P) Limited carry six months warranty against material or manufacturing defect. Defective units must be returned to Tynor Orthotics (P) Limited through the Super Dealer concerned for replacement, along with duly filled " Product Return Form" . Liability of Tynor Orthotics (P) Limited shall be limited to the cost of product only .

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Tynor Orthotics Pvt. Ltd.