Table of Contents

1	Device description	2
2	Intended use	3
3	Safety information	3
4	Package contents	5
5	General information about calluses	5
6	How to use	6
6.1	Power supply.....	6
6.2	Removal and attachment of roller	6
7	Operation	7
7.1	Application	7
7.2	Skin care after application	8
8	Cleaning and maintenance of the device	8
9	Storage	9

I. Silk'n Pedi Pro device description

1. Roller with diamond callus remover
2. Roller holder
3. On/Off switch
4. Silk'n Pedi Pro device
5. Battery compartment lid
6. Nail polishing head - rough (turquoise)
7. Nail polishing head - gentle (white)

Thank you for your trust!

Congratulations on your new Silk'n Pedi Pro callus remover and nail polisher. For safe handling of the device and in order to get to know the entire scope of features:

- Thoroughly read these operating instructions before using Silk'n Pedi Pro for the first time.
- Above all, observe the safety instructions!
- Keep these operating instructions for reference.
- We hope you enjoy your new Silk'n Pedi Pro!

2. Intended use

Silk'n PEDI Pro is intended for the removal of callus formation and for nail polishing. It is designed solely for use on humans.

It has been developed for private use and is unsuitable for commercial purpose.

Foreseeable misuse

- The device may not be used on reddened, injured or burnt skin. Similarly, healthy skin without calluses may not be treated with this device. Warts, moles and other protrusions (except calluses) must not be treated with this device. In all of these cases, this could lead to skin and/or deep tissue injury.
- Diabetics and persons with high blood pressure may only use this device after consultation with their general practitioner.

3. Safety information

Instructions for safe operation

- Before using Silk'n PEDI Pro for the first time, please read this User Manual in its entirety. Pay particular attention to sections on device use procedures, device operation, and after-use procedures. We recommend you re-familiarize yourself with this User Manual before each use of Silk'n PEDI Pro. Silk'n PEDI Pro is a powerful electrical device. As such, it should be used with special attention to safety. Please read all warnings and safety precautions before use, and strictly follow them when using Silk'n PEDI Pro.
- The device is for use by adults.

Keep Silk'n Pedi Pro and its packaging materials away from children!

WARNING! Risk of injury or skin irritation

- Do not use Silk'n Pedi Pro if there is visible damage to any of its parts or the roller.
- Do not attempt to make modifications, open or repair your Silk'n Pedi Pro™ device. Only authorized Silk'n™ repair centers are permitted to perform repairs.
- Do not file over previously damaged skin areas (e.g. through eczema or wounds), or over raised skin areas (e.g. moles or warts).

DANGER relating to batteries

- Do not subject batteries to extreme conditions, place them on radiators or expose them to direct sunlight.
- If battery acid has leaked, avoid contact with the skin, eyes and mucous membranes. In case of contact with the acid, rinse the affected area at once with plenty of clean water and seek medical assistance immediately.
- If a battery has been swallowed, medical assistance must be sought immediately.

WARNING! Avoid damage to the device

- Use only the original accessories.
- Avoid dropping the device.
- Use this device on human callused skin and human nails only. Other surfaces could damage the rollers.

4. Package contents

- 1 Silk'n Pedi Pro device
- 1 callus remover
- 2 nail polishing heads for rough and gentle finishes
- 1 set of operating instructions

5. General information about calluses

Callused skin is normal, especially on the feet. However, it should not become too thick or hard. Both of these conditions can cause pain.

Causes

An increase in callused skin can be caused by the wrong shoes, obesity, an incorrect gait, bone deformations or extremely dry skin. If you suffer from an increase in callused skin, try to determine the cause and eliminate it. Also consult your physician.

Prevention

- Regular foot baths refresh and relax your feet.
- Moisturize your feet regularly (especially after every foot bath).
- Make sure that your shoes fit well.

6. How to use

- Remove all packing material.
- Make sure all components are there.

You may use the callus remover in water, but use the nail heads only on dry nails and in a dry environment.

6.1 Power supply

1. Slide the battery compartment lid off the device.
2. If necessary, replace the old batteries with new ones.
3. Slide the battery compartment lid back onto the device.

6.2 Removal and attachment of roller

To remove, just hold the roller, gently press and pull out. To attach, push in until you hear a click.

Silk'n Pedi Pro has an RPM of 1700 +/-100.

7. Operation

The device comes complete with one callus remover and two nail polishing heads with different grain types:

- The coarse diamond grain roller is for use on thick calluses
- The rough nail polishing head (turquoise) is to even-out the nail surface
- The gentle nail polishing head (white) is for giving a smooth shine to the nail

Replace the nail polishing head/s if they appear worn.

WARNING!

Exert only light pressure. If too much pressure is exerted or if the treatment duration is too long, this can lead to skin injuries.

Stop the treatment immediately in case of pain.

7.1 Application

- You may use the callus remover in a wet environment (under the shower). Use the nails heads only in a dry environment.

Calluses

- Arrange for a surface (a towel, etc.) on which the removed callused skin can fall.
 1. Press On/Off switch. The roller starts to rotate.
 2. Apply the rotating roller slowly and carefully to the callus with as little pressure as possible.

3. Use the end of the roller on tough skin or if you want extra power.
4. To switch off the device, just press the On/Off switch.

Nails

1. Start with the rough polishing head.
2. Apply the rotating head slowly and carefully to the nail, with as little pressure as possible in order to even-out any ridges on the surface of the nail or rugged edges.
3. Follow up with the gentle polishing head to smooth out the nail and create a visible shine.

7.2 Skin care after application

- Apply body lotion, oil or hand cream to the skin.
- Do not apply any substances to the treated areas which could irritate the skin.

8. Cleaning and maintenance of the device

WARNING!

Do not switch the device on if the roller unit has been removed.

Do not use any astringent or abrasive cleaning agents.

The callus remover may be cleaned under water. Clean the nail heads with a soft brush in a dry environment only.

After Silk'n Pedi Pro Treatment

1. To switch off the device, press the On/Off SWITCH.

From time to time

- Wipe down the device with a damp cloth when necessary.
- Check whether the surface of the roller is undamaged.

9. Storage

DANGER! Keep the device out of reach of children.

NOTE: To prevent dust and other harmful substances from getting into the device, store the device with the attached roller.

If you are not going to use the device for a long period of time, remove the batteries.

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