Table of contents

1. Package contents/device parts	2
2. General information	2
3. Safety	3
4. Checking the device and package contents	5
5. Charging	6
6. Operation	6
7. Cleaning and maintenance	9
8. Storage.....	10
9. Troubleshooting	10
10. Customer service.....	11
11. Disposal	11
12. Specifications	12
13. Signs & symbols	12
14. Warranty	13

1. Package contents/device parts

EN

- | | |
|-------------------------|---------------------|
| 1 Device | 10 Body diamond tip |
| 2 LCD display | 11 Massage tip |
| 3 Mode button | 12 USB adapter |
| 4 Start/pause button | 13 USB cable |
| 5 Power button | 14 Filter, 30x |
| 6 Control buttons | 15 O-ring, 5x |
| 7 Stand/charging cradle | 16 Cleaning brush |
| 8 Fine diamond tip | 17 Pouch |
| 9 Regular diamond tip | |

2. General information

2.1 Reading and storing the user manual

This user manual accompanies this ReVit Prestige (hereafter referred to as the “device”), and contains important information on setup and handling.

Before using the device, read the user manual carefully. This particularly applies to the safety instructions. Failure to do so may result in personal injury or damage to the device.

When abroad, you must also observe country-specific guidelines and laws. Store the user manual for further use. Make sure to include this user manual when passing the device on to third parties.

2.2 Explanation of symbols

The following symbols and signal words are used in this user manual, on the device or on the packaging.

Warning: This symbol/word designates a hazard with moderate risk, which may result in death or severe injury if not avoided.

Notice: This symbol/word warns of possible damage to property.

3. Safety

3.1 Proper use

The device is exclusively designed to exfoliate your skin, stimulate the growth of healthy new cells and stimulate blood circulation in the skin for a smoother and fresher complexion. It is only intended for private household use and not suitable for commercial purposes.

Only use the device as described in this user manual. Any other use is considered improper and may result in damage to property or persons.

The manufacturer or vendor cannot be held liable for damages or injury incurred through improper or incorrect use.

3.2 Safety instructions

Risk of electric shock

-  Only connect the device if the line voltage of the socket corresponds to the data on the USB adapter.
-  Only connect the USB adapter to an easily accessible socket so that you can quickly disconnect it from the electricity supply in the event of a problem.
-  Only use the device with the enclosed parts.
-  Do not use the device if it is damaged or if the charging cradle, USB adapter or the USB cable is defective. Contact the manufacturer or customer service for replacement or repair to avoid risks.
-  Do not open the housing; instead, have a qualified professional perform repairs. Contact our customer service for this. Liability and warranty claims are waived in the event of repairs performed by the user, improper connection or incorrect operation.
-  Only parts that comply with the original device data may be used for repairs. This device contains electrical and mechanical parts which are essential for providing protection against sources of danger.
-  Do not use attachments not specified by the manufacturer.
-  Do not use the device while showering or bathing. Do not place the device such that there is a risk of it falling into a tub or sink.

-  Do not immerse the device or any of the parts in water or other liquids.
-  Never reach to retrieve an electrical device if it has fallen in water. In such a case, immediately disconnect the USB adapter.
-  Never insert any objects into the housing.
-  Always switch the device off and pull the USB adapter out of the socket when you are not using it, when you are replacing the parts, when you intend to clean it or in the event of a malfunction.
-  Never touch the USB adapter or USB cable with wet or damp hands.
-  Keep the device and all parts away from open flames and hot surfaces.
-  Lay the USB cable so that it does not pose a tripping hazard.
-  Do not kink the USB cable and do not lay it over sharp edges.
-  Only use the device indoors.

Danger for children and persons with impaired capacities, or lack of experience and knowledge.

-  The device is not intended for use by persons under the age of 16 or persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the device by a person responsible for their safety. Children must not play with the device. Cleaning and user maintenance must not be performed by unsupervised children.
-  Keep the device and all parts away from children. Small parts may cause choking if swallowed or inhaled.
-  Make sure that children do not play with the plastic wrapping. They may get caught in it when playing and suffocate.

Risk of damage

-  Do not use the device on the eye lid, inside the eye orbit, around areas that contain inflammatory acne, or areas that have rashes or wounds.
-  Do not focus the device on one spot on your skin or pass the device over the same area more than 3 times.

-  Never leave the device unattended while it is in operation.
-  Stop using the device immediately if you develop any abnormal reactions such as itching, uncomfortable inflammation, or scratches.

Consult your doctor before using the device if..

-  you suffer from allergies, skin irritation or skin diseases,
-  ~~you are pregnant,~~
-  you are not sure whether the device is suitable for you.

Risk of damage

-  Place the device on an easily accessible, level, dry and sufficiently stable surface. Do not place the device on the edge of the surface.
-  Handle with care. Dropping the device may damage the housing or motor.
-  Keep a filter in place at all times to prevent foreign substances from entering the motor.
-  Stop using the device if any parts exhibit cracks or breaks or are deformed. Only replace damaged parts with corresponding original spare parts.
-  Do not place or store the device in damp areas or rooms.

4. Checking the device and package contents

Risk of damage

-  If you are not cautious when opening the packaging with a sharp knife or other pointed object, you may quickly damage the device. Be very careful when opening it.
 1. Take all parts out of the packaging.
 2. Check to make sure that the delivery is complete (see **Fig. A**).
 3. Check whether the device or individual parts are damaged. If this is the case, do not use the device. Contact our customer service.

5. Charging

EN

Risk of injury

-  Do not charge when the device is in use or if it is wet. Make sure it is dry before charging.
-  Only charge the device with the enclosed charging cradle, USB adapter and USB cable.
-  Charge the device at temperatures between 0 and 40 °C.
-  Always unplug the device after charging has been completed.

It takes about 3 hours to fully charge the device. When the charging process is completed, the LCD display shows the charging status as full.

Option 1 (via mains supply):

1. Place the device **1** onto the charging cradle **7** by fitting the hole in its bottom to the pin of the charging cradle.
2. Connect the USB cable **13** to the charging cradle and connect its USB connector to the USB adapter **12**.
3. Plug the USB adapter into a mains socket.

Option 2 (via USB port):

1. Connect the USB cable **13** to the charging cradle **7** with the device properly installed.
2. Connect its USB connector to a compatible USB power source.

6. Operation

-  Clean all parts of the device before first use as described in the chapter "Cleaning and maintenance".
-  Make sure the filter, O-ring and the desired tip have been installed properly before use.
-  Before treating an area for the first time, perform a skin sensitivity test on a spot. If you have no adverse reaction, start your treatment.

- Always start on the lowest suction level and stroke no more than 2 times over each area for your first treatment, or until you have determined your ideal treatment setting.
- Do not use additional treatment products such as oils or ointments during treatments.
- The device will switch itself off after 9 minutes of continuous operation. ~~It is recommended not to perform a facial treatment for more than 4 minutes.~~
- When the power level is low, the device will alert with 2 beeps every 30 seconds and the battery symbol on the LCD display will blink. Refer to the chapter "Charging" to charge the device.

6.1 Treatment tips

- Use the device maximum 2 times a week. Treat your face for maximum 4 minutes per treatment. Treat the parts of your body by passing each area no more than 3 times.
- Select a suitable tip for your treatments:
 - **Fine diamond tip 8**: the finer treatment tip for the delicate and hard-to-reach areas of the face, such as around the eyes, nose and mouth.
 - **Regular diamond tip 9**: the standard treatment tip for all skin types for removing dead skin cells and stimulating the cell regrowth and the blood circulation. Suitable for face and neck.
 - **Body diamond tip 10**: the larger treatment tip for body parts such as arms and legs.
 - **Massage tip 11**: for massaging the skin and stimulating the blood circulation. Suitable for all skin types.
- Thoroughly wash your face or body part and pat dry with a soft towel to leave your skin bare before using the device.
- Gently move the tip across the skin using short and even strokes. For facial area, follow the below instructions:
 - **Cheeks**: Stroke from inside outwards.
 - **Nose**: Stroke downwards covering entire nose.
 - **Chin**: Stroke downwards along jawline.
 - **Forehead**: Stroke from centre outwards on each side.
 - **Eye area**: Stroke outwards between eye and ear. Do not stroke directly below the eyes or on the eyelids, and stay outside the eye orbit.

- **Neck:** Stroke downwards.
- After treatment, wash the treated area with cold water to remove remaining debris and gently pat dry. Apply a toner, serum or cream to minimize your pores and soothe your skin for the optimal efficacy of your treatments.
 - Light bruising and redness is completely normal after treatments and should subside within 24 hours. Effective dermabrasion treatments gently remove the outer layer of skin which may cause slight discomfort and redness.
-  Stop using the device immediately if you develop any abnormal reactions such as itching, uncomfortable inflammation, or scratches.

6.2 Attaching the tip

-  Make sure the device is switched off and unplugged before attaching or changing the tip.
 -  Always make sure a filter and O-ring are in place. Refer to the chapter “Replacing the filter” or “Replacing the O-ring” for replacement if necessary.
1. Attach the desired tip (8 / 9 / 10 / 11) to the device.
 2. Make sure the tip is firmly in place before use.

6.3 Operating in Auto mode

Auto mode is only for facial treatments.

1. Press the power button  to switch on the device.
2. Press the mode button  to select “Auto mode”.
3. Press the control buttons  to select the desired intensity level. Level 1 is the lowest level; while level 3 is the highest level.
4. Press the start/pause button  to begin the treatment.
 - Treat the facial area indicated on the LCD display. The LCD display also shows the countdown of the treatment time for each facial area. You will hear 2 beeps when the time elapses.
 - Move to the next facial area indicated on the display after each 2 beeps until the treatment cycle is completed.
 - During the treatment, you can press the Start/pause button and then press the control buttons anytime to switch to another facial area.

6.4 Operating in Manual mode

Manual mode is for both facial and body treatments. This mode allows you to follow your own treatment regime. It will not inform you about treatment times, areas or give you instructions on how to use it. The timer however, will keep track of the treatment time that has passed.

1. Press the power button **5** to switch on the device.
2. Press the mode button **3** to select "Manual".
3. Press the control buttons **6** to select the desired intensity level. Level 1 is the lowest level; while level 3 is the highest level.
4. Press the start/pause button **4** to begin the treatment.

7. Cleaning and maintenance

Risk of damage!

-  Always turn off and unplug the device before cleaning or replacing the spare parts.
-  Never immerse the device or any parts in water or other liquids. Make sure that no water or other liquids penetrate the housing.
-  Do not use any aggressive cleaners, brushes with metal or nylon bristles, or sharp or metallic cleaning utensils such as knives, hard scrapers and the like.

7.1 Cleaning the housing and the charging cradle

1. Wipe the device **1** and the charging cradle **7** (avoid contacting the pin of the charging cradle or the hole on the bottom of the device) with a damp cloth.
2. Wipe it completely dry with a dry cloth.

7.2 Cleaning the tip

1. Remove the tip (**8** / **9** / **10** / **11**) from the device and gently clean its inner side with the brush **16** after each use to remove dead skin cells and debris.
2. Wipe the outside of the tips thoroughly with a damp cloth.
3. Let it dry completely.

7.3 Replacing the filter

The filter should be replaced regularly for optimal results. Always make sure that the filter is in place to prevent particles such as dead skin cells or debris from entering the motor.

1. Gently remove the tip (8 / 9 / 10 / 11) and discard the filter 14.
2. Replace with a new filter and put the tip back into place.

7.4 Replacing the O-ring

The O-ring should only be replaced if it is damaged or worn. Always make sure that the O-ring is in place, otherwise the device will not work properly.

1. Gently remove the tip (8 / 9 / 10 / 11) and discard the O-ring 15.
2. Replace with a new O-ring and put the tip back into place.

8. Storage

- All parts must be completely dry before being stored in the pouch 17. You can also use the charging cradle 7 as a stand to hold the device.
- Do not wind the USB cable 13 round the device.

9. Troubleshooting

Some problems may be caused by minor faults that you can fix yourself. To do so, follow the instructions in the following table. If it is still not possible to resolve the problem, contact customer service. Do not repair the device yourself.

Fault	Possible cause and solution
The device won't turn on.	Make sure the USB adapter is securely plugged in the charging cradle, and plugged into a working mains socket for at least one hour. The battery symbol on the LCD display should indicate that device is charging. If the battery symbol does not come on, contact our customer service.

<p>The device has lost suction.</p>	<ul style="list-style-type: none"> • Check the tip closely for any hairline cracks, punctures, or spaces. • Remove the tip and ensure that the O-ring is in place. If you notice cracks or that the O-ring is missing, replace it with a spare ring. • Remove the tip and ensure that the filter is in place. If you notice that the filter has worn down or if you see a lot of dead skin cell and debris, replace it with a new filter. • If these solutions do not work, contact our customer service.
<p>The tip (8 / 9 / 10) is not coarse enough.</p>	<p>The tip will accumulate skin cell build-up over time, making it feel smoother. Clean the tip with a damp cloth to keep it clear of skin debris. The tips should be replaced every year, or earlier if they are worn. You can purchase replacement tips and filters through the retailer where you bought the device, or you can order directly from our website: www.silkn.com.</p>

10. Customer service

For more information about this or other Silk'n products, visit our website: www.silkn.com. This manual is also available as a PDF download from www.silkn.com.

Report unexpected operations or events, or any other problem with your device to Customer Service. Our contact details are:

Toll free number: 1-877-367-4556

Email: contact@silkn.com

11. Disposal

Dispose of the packaging separated into single type materials. Dispose of cardboard and carton as waste paper and film via the recyclable material collection service.

Do not dispose with household waste! The WEEE directive requires that you dispose of this electrical device (including all accessories and parts) at a designated recycling collection point for electrical devices. This ensures that the device is recycled professionally, and prevents the release of harmful substances * into the environment. Make sure that the device is in an uncharged state before disposal.

*Products labelled with Cd = cadmium, Hg = mercury, Pb = lead.

12. Specifications

EN

Model:	ND-9002
Adapter: Power input: Power output:	PS06C050K0500UU 100–240 V, ~50/60 Hz, 0.25A 5V - 500mAh
Protection class: - device: - USB adapter:	III II
Dimension (H × W × D) - device: - charging cradle:	185 × 48 × 41mm 54 × 77 × 77mm
Weight:	316g

13. Signs & symbols

	Declaration of Conformity: Products labelled with this symbol meet all applicable provisions of the European Economic Area.
IPX0	No protection against water.
	This symbol identifies electrical devices that belong to the protection class 3 (device).
	Do not dispose with household waste! The WEEE directive requires that you dispose of this electrical device (including all accessories and parts) at a designated recycling collection point for electrical devices.
	Devices marked with this symbol are only suitable for household use (dry indoor rooms).
	This symbol identifies electrical devices that belong to the protection class 2 (USB adapter).

14. Warranty

What Does This Warranty Cover?

Home Skinovations, its applicable affiliates (hereinafter: Home Skinovations) or the authorized dealer from whom you purchased your Silk'n product (hereinafter: the Device) shall repair or replace the Device, if any material defects caused by Home Skinovations are discovered in the Device during the warranty period, whether in the body of the Device or its operation.

Who is entitled to warranty by virtue of this warranty certificate?

Only you, the original buyer of the Device who purchased the Device from Home Skinovations or its authorized dealers, are entitled to receive the warranty specified in this warranty certificate. In order to be entitled to the warranty service, you need to provide a proof that you are the original buyer of the Device by presenting an original receipt for purchase of the Device. This warranty is not an international warranty. The warranty will only be honored in the country in which the authorized dealer or the Home Skinovations entity from whom you purchased the Device, are located. For a list of Home Skinovations' authorized dealers see: www.silkn.com/For_International_Visitors.

What period does this warranty certificate cover?

The warranty covers repair requests received within one year (1) from the Device's original purchase date. Repair requests after one year (1) from the date of purchase are not covered by this certificate of warranty. Repair requests regarding a Device that was repaired or replaced must be submitted no later than: (a) one year (1) from the original purchase date of the Device or (b) ninety (90) days from the date that the repaired or replaced Device was returned to you, the later to occur. If a sticker is attached to your Device stating a warranty period longer than that indicated above, the warranty period mentioned on the sticker will apply, and the Device will be repaired or replaced as part of the warranty until such period expires.

If you purchased the Device in an EU country, you have an additional year (according to the EU-Directive 1999/44) to claim for deficiencies or malfunctioning of your Device, which were present at the time of purchase, but were not apparent then, as determined by Home Skinovations or the authorized dealer. In order to make this extended warranty effective you must contact the authorized dealer from whom you purchased the Device or the local service center and notify them about the production fault immediately, and no later than 60 days after detecting the fault.

How to get service for the warranty?

For further instructions how to obtain warranty service and for contact details of the service centers close to you, visit: www.silkn.com/For_International_Visitors, or contact your local authorized dealer, or call our Help Desk at:

North America: Tel: 1-877-DO-SILKN (1-877-367-4556) / Email: contact@silkn.com

European countries: (Including all member states of the European Union, Switzerland, Norway and Russia) Tel: +31-0-180-330-550 / Email: silkn@inno-essentials.com

Japan: Tel: 0120-110-469 / Email: beauty0120110469@sensepil.jp

Israel: Tel: 1-700-700-309, Email: info@silkn.com

Other countries: Email: info@silkn.com

When asking for customer services you may be asked to fill in a form for returning the product and you shall receive a reference number. In case of a defect that is not covered by the warranty for whatever reason, Home Skinovations can charge you for the shipping and repair costs if you agreed to repair the Device at your expense.

What is not covered by this certificate of warranty?

This warranty does not cover any damage resulting from an accident, misuse (including using the Device other than as described in the use instructions accompanying the Device), defective use, lack of reasonable care, irregular and unreasonable use including without limitation commercial or professional use, adding parts not supplied together with the Device, losing parts, using unauthorized spare parts or exposing the Device to any electrical current other than the indicated current. The warranty will be cancelled if the Device was treated by anyone not authorized by Home Skinovations and Home Skinovations will not be liable for any damage resulting from repairs performed by anyone not authorized by Home Skinovations. The warranty does not include replacing used non-durable parts.

WARRANTY DISCLAIMER

THE EXPRESS WARRANTY SET FORTH IN THIS DOCUMENT IS IN LIEU OF ALL OTHER WARRANTIES INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY, AND FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS THE SOLE AND EXCLUSIVE REMEDY FOR ANY DEFECT OF THE DEVICE.

LIMITATIONS

IN NO EVENT SHALL HOME SKINOVATIONS, ITS AFFILIATES OR AUTHORIZED DEALERS, BE LIABLE FOR ANY INDIRECT, SPECIAL, PUNITIVE, EXEMPLARY, INCIDENTAL OR CONSEQUENTIAL DAMAGES (INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF USE OR CHARGES FOR SUBSTITUTE EQUIPMENT), REGARDLESS OF THE LEGAL THEORY ON WHICH THE CLAIM IS BASED EVEN IF THEY HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. IN THE EVENT THAT THE REMEDIES SET FORTH IN THIS WARRANTY ARE HELD TO FAIL OF THEIR ESSENTIAL PURPOSE OR IF HOME SKINOVATIONS, ITS AFFILIATES OR AUTHORIZED DEALERS ARE HELD TO BE LIABLE FOR DAMAGES ON THE BASIS OF ANY LEGAL THEORY, IN NO EVENT SHALL HOME SKINOVATIONS; ITS AFFILIATES OR AUTHORIZED DEALERS, LIABILITY FOR DAMAGES EXCEED THE WHOLESALE PRICE OF A SINGLE UNIT OF THE DEVICE AS OF THE DATE THAT THE CAUSE OF ACTION ACCRUES. YOU ASSUME ALL RISK AND LIABILITY FOR LOSS, DAMAGE OR INJURY TO YOU AND YOUR PROPERTY AND TO OTHERS AND THEIR PROPERTY ARISING OUT OF THE USE OR MISUSE OF, OR INABILITY TO USE, THE DEVICE EXCEPT TO THE EXTENT CAUSED DIRECTLY BY THE NEGLIGENCE OF HOME SKINOVATIONS.

THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM COUNTRY TO COUNTRY. SOME COUNTRIES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES OR IMPLIED WARRANTIES, SO THE ABOVE DISCLAIMERS AND LIMITATIONS MAY NOT APPLY TO YOU.

Legal Notice

Home Skinovations Ltd reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Home Skinovations Ltd is believed to be accurate and reliable at the time of publication. However, Home Skinovations Ltd assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Home Skinovations Ltd.

No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Home Skinovations Ltd. Data is subject to change without notification.

Home Skinovations Ltd has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Home Skinovations Ltd. Specifications are subject to change without notice.

SILK'N and the SILK'N logo are registered trademarks of Home Skinovations Ltd, Tavor Building, Shaar Yokneam P. O. Box 533, Yokneam 2069206, ISRAEL, www.silkn.com info@silkn.com.

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