## PRODUCT CARE

Hand wash with mild detergent and water below 30°C without wringing.

Dry in shade on a flat surface.

Do not dry-clean.

Do not iron.

Do not bleach.

Don't stretch the band for long duration.

Avoid contact with oils/balms/chemicals/solvents

Keep away from sharp objects that can puncture the product

Keep away from direct sunlight or high temperatures

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
## English

### Old-age, Rehabilitation and Recovery

Exercise	Back Pain	Hip Care	Shoulder	Arthritis	Osteoporosis	Diabetes
Side Bend		✓			✓	
Reverse Fly			✓	✓	✓	
Leg Press				✓		✓
Leg Curl	✓			✓	✓	

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1.

Side Bend [Target: oblique muscles in the abdomen]

- Feet shoulder-width apart.
- One end under the right foot, other end held at level of hip.
- Lean towards left against the resistance.
- Repeat, then switch sides.

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2.

Reverse Fly [Target: Trapezius muscles of Upper Back]

- Stand straight
- Hold the band with both hands, shoulder-width apart
- Open arms sideward against the resistance.
- Repeat

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### Biceps Curl [Target: Biceps]

- Standing on the center of the band, firmly hold both ends.
- With elbows close to your ribs, slowly bend arms against band resistance.
- Slowly return to start position and repeat.

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6.

### Pelvic Lift [Target: Glutes and Para-spinal muscles in lower back]

- Lie on back, place the band over hips and secure it firmly in your hands.
- Slowly lift hips to a comfortable height engaging glutes without momentum. Hold for 1-2 seconds.
- Return to start position

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## Tyroband

### English

Tyrobands are resistance bands which offer a calculated and consistent counter force for any kind of Isotonic exercise plan. Series of bands with increase in resistance force are available with color codification. Each one to suit the specific needs of each patient and a series of them to closely match the rehabilitation and strength-gain exercise needs.

#### CONSTRUCTION:

A ribbon made of high quality polymer/latex with consistent properties.

#### APPLICATIONS:

- Post injury/post operative rehabilitation
- Arthritic rehabilitation therapy
- Any disease related neural /muscular weakness
- Old-age weakness
- Joint function recovery
- Rehabilitation injuries involving joints
- Improve flexibility and range-of-motion of joints
- Increase muscle strength and endurance
- Gym and Sports

#### INSTRUCTIONS: [1]

- Use this product under strict guidance of a qualified doctor.
- Discontinue use and seek guidance of a qualified doctor, in case of
  - Impaired sensation
  - Pain increases or persists
- Consult your doctor in case of concurrent use of another medical device.
- Misuse or improper use of the device can lead to adverse effects and decreased effectiveness.
- Follow instructions mentioned on the box carefully.

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