

背面

(TechCare	(TechCare	ClechCae	(TechCare	(TechCare	(TechCare	(Tech Care	(Tech Care	
Certificate of Quality & Warranty Card Product Model Product Name Date of Purchase Purchase Location Retailer's Name Retailer Telephone Number Registerer Info Name Adress Phone Number Email You can register your product online at www.techcaremassager.com	2. Large contacting surface makes the belt easily to use. Functions 1. Eliminating fatigue of arms, waist abdomen and legs. 2. Is specifically designed to give optimum relief for the lower back. Safety Precautions 1. Do not use the belt if you have heart disease. 2. Do not apply the belt over the neck or mouth, as severe muscle spasms of the laryngeal and pharyngeal portion muscle s may occur. 3. Do not use the belt over head or chest. 4. Do not apply the belt over infected or irritated skin areas. 5. Do not use the belt during pregnancy. 6. A 30-minutes daily routine is recommended. 7. Keep the belt out of reach of children Cautions The belt should be cleaned regularly with warm water. Discoloration of the conductive area is, however, normal. Do not wash the belt in the washing machine nor use fabric conditioner to wash. Check the belt before every use.	ABDOMINAL & LOWER BACK BELT OPERATING INSTRUCTIONS 1. Connect the belt and place the belt where needed. 2. Press the two snap connectors at the end of the auxiliary wire respectively on either button found beneath the cloth in which the Techcare logo is printed. 3. Insert the other end of the auxiliary wire into the output jack found at the bottom of controller. 4. Turn on controller. MASSAGE SHOES Techcare Massage Shoes are high-Tech product. Regular and lasting use can relieve fatigue and enhance circulation. Features 1. Combination of soft and hard material makes the massage shoes flexible, comfortable and durable.	Functions 1. Patients with heart disease, hypertension, calf pain and women who are pregnant are prohibited from using the massage shoes, People with unexplained joint and muscle inflammation within 72 hours of using massage shoes should consult a doctor. 2. Not allowed for children usage 3. Do not operate this product in wet places such as bathing/showering. 4. Do not attempt to disassemble or repair. Cautions When massaging, please keep the points and the massage shoes clean. For most users the output strength of this machine is high enough, if the maximum strength level of this device does meet your requirements, you may moisten the massage points with a little clean water to enhance the effect. MASSAGE SHOES INSTRUCTIONS Snap on the wire to the shoes and the other end of the auxiliary wire to bottom of the controller.	Place bare feet into the massager shoes. Operate the controller to massage your feet(please refer to the Operating Instructions of the controller) 3. Cleaning and maintenance A. For the control unit: 1) To keep the controller clean, use a soft and dry cloth for dust or a soft damp cloth for any dirt and smudges. Do not use any cleaning solutions to clean the controller and its pads. 2) Do not use or store the device where there are magnetic fields or electric waves (near TV set or speakers). 3) Do not place the devices in areas of high temperature, high humidity, or under direct sunlight. 4) Keep the device out of reach of children. 5) All worn accessories should be disposed of according to your local regulations. * Refer to the user manual of the pads manufacturer (510K092546), or reference to the following suggestions: 1) Unplug the output cord from the output jack of the controller after each use. Cover both pads with the protective film before	storage. Never fold the electrode pads. 2) Between uses, store the electrodes in the reusable bag in a shady place. Stockroom temperature: +5°C~+27°C (41-80°F) and humidity of 30%~80%. No need to sterilize. 3) Never apply the pads to any other surface other than your skin. If the pads become soiled or dirty, the adhesive power may decrease. In this case, moisten the surface of the pads with water and wipe away the dirty portion. This will allow a temporary restoration of the adhesive power. However, too much water will result in loss of the adhesive power. 4) The life of the electrodes varies depending on skin conditions, storage, amount of use, type of stimulation, and stimulation site. Electrode life may be extended by carefully following this Instruction for Use. The expired electrodes are to be recycled and do not harm environment. Warning: The electrodes are intended for single patient use only! 4. Adverse Reactions 1) Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin; 2) Patients may experience headache and other painful sensations during or following the application of electrical stimulation. 3) Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.	Tips: Dear customers, in order to keep the stability of the machine, we suggest you charging it insevery 3 month, even you don't use it, thank you!!	Need Replacement Pads? FDA 510(k) Cleared Electrode Pads with Discounted Price www.amazon.com/techcaremassagers Techcare www.techcaremassager.com	Plus Combo Set SM TENS & PMS User's Manual Please read the user's manual carefully before using the product

两处虚线为易撕处

展开尺寸: 675*120MM正负1MM

折叠尺寸: 75*120MM

材质: 128G铜版纸 四色印刷 九折页