

背面

(TechCare	(TechCare	(TechCare	(Tech Care	(TechCare	(TechCare	
Certificate of Quality & Warranty Card	3. Cleaning and maintenance  A. For the control unit:  1) To keep the controller clean, use a soft and dry cloth for dust or a soft damp cloth for any dirt and smudges. Do not use any cleaning solutions to clean the controller and its pads.  2) Do not use or store the device where there are magnetic fields or electric waves (near TV set or speakers).  3) Do not place the devices in areas of high temperature, high humidity, or under direct sunlight.  4) Keep the device out of reach of children.  5) All worn accessories should be disposed of according to your local regulations.  * Refer to the user manual of the pads manufacturer (510K092546), or reference to the following suggestions:  1) Unplug the output cord from the output jack of the controller after each use. Cover both pads with the protective film before storage. Never fold the electrode pads.  2) Between uses, store the electrodes in the reusable bag in a shady place. Stockroom temperature: +5°C~+27°C (41-80°F) and humidity of 30%~80%. No need to sterilize.  3) Never apply the pads to any other surface other than your skin. If the pads become soiled or dirty, the adhesive power may decrease. In this case, moisten the surface of the pads with water and wipe away the dirty portion. This will allow a temporary restoration of the adhesive power. However, too much water will result in loss of the adhesive power.  4) The life of the electrodes varies depending on skin	conditions, storage, amount of use, type of stimulation, and stimulation site. Electrode life may be extended by carefully following this Instruction for Use. The expired electrodes are to be recycled and do not harm environment.  Warning: The electrodes are intended for single patient use only!  4. Warnings  1) The long-term effects of chronic electrical stimulation are unknown.  2) Do not use stimulation over the main arteries around your neck.  3) Do not use stimulation over the neck or mouth. Severe spasm of muscles may occur and the contractions may be strong enough to close the airway or cause difficulty in breathing.  4) Do not apply stimulation directly to the heart area across the chest or on the chest.  5) Do not use stimulation on your head or across your head.  6) Do not use stimulation over swollen, infected, or inflamed areas or skin eruptions, e.g., phlebitis, thrombophlebitis, varicose veins.  7) Do not use stimulation over, or in proximity to, cancerous lesions.  8) Do not use stimulation when charging the device.	5. Precautions  1) Safety of powered muscle stimulator for use during pregnancy has not been established.  2) Caution should be used for patients with suspected or diagnosed heart problems.  3) Caution should be used for patients with suspected or diagnosed epilepsy.  4) Caution should be used in the presence of the following:  a. When there is a tendency to hemorrhage following acute trauma or fracture;  b. Following recent surgical procedures when muscle contraction may disrupt the healing process;  c. Over the menstruating or pregnant uterus;  d. Over areas of the skin which lack normal sensation.  5) Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium or alternate electrode placement.  6) Electrode placement and stimulation settings should be based in the guidance of the prescribing practitioner.  7) Power muscle stimulators should be kept out of the reach of children.  8) Power muscle stimulators should be used only with the leads and electrodes recommended for use by the manufacturer.  9) Portable power muscle stimulators should not be used while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.	6. Precautions  1) Do not use this device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device, because this may cause electric shock, burns, electrical interference, or death.  2) Do not use this device on patients whose pain syndromes are undiagnosed.  3) Do not use this device during pregnancy.  4) Do not use this device on babies or infants who cannot express themselves.  7. Adverse Reactions  1) Patients may experience skin irritation and burns beneath the stimulation electrodes applied to the skin;  2) Patients may experience headache and other painful sensations during or following the application of electrical stimulation.  3) Patients should stop using the device and should consult with their physicians if they experience adverse reactions from the device.  Tips:  Dear customers, in order to keep the stability of the machine, we suggest you charging it in every 3 months, even you don't use it, thank you!!	Need Replacement Pads?  FDA 510(k) Cleared Electrode Pads with Discounted Price www.amazon.com/techcaremassagers  Techcare  www.techcaremassager.com	SE SM TENS & PMS User's Manual  Please read the user's manual carefully before using the product

两处虚线为易撕处

展开尺寸: 525\*120MM正负1MM

折叠尺寸: 75\*120MM

材质: 128G铜版纸 四色印刷 七折页

封面