

How to Use the Surround[®] Toothbrush

Latex-free!



1

Apply a small amount of toothpaste to the center of the toothbrush.



2

Place the Surround[®] Toothbrush over the teeth so that one brush head rests on the back (tongue side) of the teeth, the other brush head rests on the front (face side) of the teeth, and the center brush head rests on the biting surface.



3

Move the brush back and forth gently in short (tooth-wide) strokes. If necessary, tilt the brush towards the face side of the teeth and then the tongue side as you brush, to be sure the bristles reach down to the gums.

The Surround[®] Toothbrush has a unique brush head that surrounds the teeth to clean front, back and biting surfaces at the same time. The side bristles are positioned at the 45-degree angle recommended by the American Dental Association to clean both teeth and gums.

The full length of the Surround[®] Toothbrush is covered with a soft, rubber-like material to help protect the oral environment if the person moves unexpectedly during brushing. The brush head and neck are reinforced with a flexible material to allow the brush to bend, not break, if the user bites down on the toothbrush.