
Dual Digital Thermometer

FDTH-V0-13

easy@Home®
Digital Clinical Thermometer

-EN-

Overview

The Easy@home Dual Digital Thermometer (FDTH-V0-13) is intended for the measurement and monitoring of human body temperature by consumers at home. The Dual Digital Thermometer has two modes - rapid read thermometer or basal body temperature (BBT) thermometer. In the mode of rapid read thermometer, the oral body temperature can be measured in about 8 to 12 seconds. When you feel unwell, you may use this mode to help you measure your body temperature quickly. The BBT thermometer mode is design to help measure BBT for ovulation tracking purpose.

Note:

- This product is only for monitoring body temperature, not for diagnosis of disease.

What is Basal Body Temperature, and why do you need a digital basal thermometer that has BBT thermometer mode to measure it?

Basal body temperature (BBT) is your temperature when you are fully at rest. It is taken immediately upon awakening before any activity. Unlike traditional thermometers, basal body thermometers can measure with increased accuracy to easily detect your fertile window. The Digital Basal Thermometer is sensitive enough to measure subtle changes in body temperature with the high precision of 1/100° degree, and only takes about 60-120 seconds for the reading with proper use in the BBT thermometer mode.

Note: It is important to record your temperature every day to identify your BBT trend for ovulation prediction.

How does your Basal Body Temperature (BBT) help you pinpoint ovulation? What do you need to know?

Your BBT temps will likely surge right after ovulation, with a 0.20°C-0.60°C (0.40°F-1.00°F) or more increased shift.

What's included

- Thermometer * 1
- Manual * 1
- Storage case * 1

than 30 minutes, so that it reaches the equilibrium state before measuring.

3. Please do not measure under strong electromagnetic interference environment (working microwave oven, induction cooker, mobile phone nearby, etc.), this may cause incorrect measurement or the inability to measure.

4. Users taking measurements after taking medications that change their body temperature (e.g. aspirin, acetaminophen, ibuprofen), which may cause measurement bias.

5. The subjects may cause temperature fluctuations after intense exercise, crying, eating, etc. It is recommended to rest for 30 minutes before measurement.

6. When measuring, please keep quiet and steady, according to the measurement method shown in this manual, so as not to cause the measurement deviation caused by improper operation.

7. The slight fluctuation of human body temperature is normal, and the body temperature collected from different parts will also vary. It is recommended to measure and observe the body temperature in the same part of the mouth several times under similar conditions.

8. Different thermometer brands cannot be compared. It is recommended to use the same thermometer brand to measure temperature to better know the temperature fluctuations.

9. Keep new and used batteries out of the reach of children. Seek immediate medical attention if a battery is ingested, and follow any other consensus medical advice from a trusted medical advisor.

10. Choking Hazard-Small parts, not for children under 3 years or any individuals who tend to place inedible objects in their mouths.

Warnings or precautions

1. Do not forcibly collide, drop, step on, or shake the product, as it may cause product failure or cause inaccurate measurements.

2. Do not disassemble, repair, or modify the product by yourself.

3. Do not use in electromagnetic interference environment, which may affect the measurement results.

4. Please keep the product in accordance with the product transportation / storage conditions to avoid use problems or measurement abnormalities.

5. If thermometer is not used for a long time, the battery should be removed to prevent battery leakage or corrosion.

Standards

FDTH-V0-13 complies with the following standards

IEC 60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance

IEC 60601-1-2 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Colateral standard: Electromagnetic compatibility. Requirements and tests

IEC 60601-1-11 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Colateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process

EN ISO 15223-1 Medical device - symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements

ISO 80601-2-56 Medical electrical equipment - Part 2-56: Particular requirements for basic safety and essential performance of clinical thermometers for body temperature measurement

Warranty

This product is warranted by the manufacturer for one year from the date of retail purchase. It does not cover damages or wear resulting from an accident, misuse or abuse, commercial use, or an unauthorized adjustment or repair of the product.

Explanation of symbols

	Refer to instruction
	Keep out of reach of children
	Manufacturing date
	Temperature limit of -25°C - 55°C (-13°F - 131°F)
	Atmospheric pressure limitation of 70kPa - 106kPa
	CE marking of conformity
	Authorized representative in the European Community

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