# invertrac<sup>®</sup>

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## A - Femur Length Adjustment Dials

- five positions obtainable.
- start dial at 2 and adjust up for adults.
- dial settings must be identical.
- obtain comfort behind knee.
- move hips forward over pad allowing torso to hang straight.

### Please Note

Dial adjustment correlates with a higher number for a longer leg.

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## IMPORTANT

Check all nuts and bolts on a regular basis.

Total Inversion Traction Workout should take an average of only 2-3 minutes. Repeat session up to three times daily or as directed by your Specialist or Physician.

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## B - Height Adjustment – See STEP 1

### Please Note

During initial Invertrac sessions, user may experience discomfort due to congestion, however with repeated use this condition will diminish allowing maximum traction to occur.

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## C - Adjust thigh dimensions

- latch adjusts up or down
- obtain comfort behind knee without restriction or pressure

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**Medical Warning** - Consult Physician or Specialist prior to use of Invertrac if any medical conditions such as Hypertension, Glaucoma, or Osteoporosis is suspected.

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## Step 1

Step into unit and insure pelvic (hip) bones are a minimum of 3 cm above the top of thigh support/pad. Height adjustment may be required.

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## Step 2

Lean forward while holding handlegrips. To raise restraint bar, raise one foot.

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## Step 3

Grasp side-grips and tip forward. Continue to M-Beam until fully inverted. Relax.....

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## Step 4

This is the correct position for natural traction. Head should not be touching floor, place hands in above position. Adjust height if necessary.

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## Step 5

**Warning** - The following exercises should not be attempted without prior medical consent.

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Spinal Rotations - Keep back straight, turn slowly to left and hold. Return to center and repeat right.

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## Step 6

Extra-Gentle Extensions – hold M-Beam or sidegrips, slowly straighten arms, relax back muscles.

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## Step 7

Extension Exercise – arch from 45 - 90 degrees using back muscles, hold 5 seconds, relax 5 seconds, repeat.

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## Step 8

Gentle Flexion Exercise. Raise hands gently toward knees, hold and relax.

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## Step 9

To dismount hold both sidegrips evenly and come up slowly insuring the back is not twisted. Keep knees bent.

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## Step 10

Transfer hands to handgrips and return feet to floor.

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## Step 11

Inversion traction session completed. Step over restraint bar, relax for 2 minutes.

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