


RELAXONCHAIR

YUKON-4D

INSTALLATION GUIDE

&

USER MANUAL



Please read SAFETY INSTRUCTIONS carefully before using the chair.



CAUTION

Do not plug in the chair
until the end of the installation.


STEP 1: PREPARING ASSEMBLY

Checklist

Check all the boxes to make sure that you have everything you need. Make sure that the chair is reclined and is away from objects and obstructions for safe assembly.

Box 1 (Main Carton Box): Massage Chair Main Body & 2 Accessory Boxes

Massage chair main body



Accessories
from the accessory boxes



Headrest

Additional
Head cushion

Backrest



Controller
pocket



Electric
power cable




The armrest
exterior parts



Controller
& Cable

Box 2 (Armrest Carton Box): Left and right armrests




Left armrest

Right armrest

Box 3 (Footrest Carton Box): Footrest

Footrest




Preparing to Install the Footrest


* Locate the body of the chair on a soft surface in order to begin assembly. Identify the location of the footrest screws and attachments on the body of the chair. Using a Phillips-head screwdriver, remove the footrest screws and put aside for later.

STEP 2: INSTALLING THE FOOTREST


1. Place the footrest near the body of the chair, and unscrew the the footrest screws.




2. Plug in the electrical cord and air hose into their respective socket and nozzle.




3. Lift up the footrest attachment covers, then carefully lift up the footrest to place the footrest into the footrest attachments.



4. Close the footrest attachment covers, and screw the footrest screws in.



5. Remove the velcro covering from the footrest flap, and attach the velcro to the body of the chair.




STEP 3: INSTALLING THE ARMRESTS


Preparing to Install the Armrests

Locate the left-side air hoses on the body of the chair that will be connected to the left-side arm panel. Locate the right-side air hoses and electrical plug on the body of the chair that will be connected to the right-side arm panel.

Note: The hoses and electrical plug may be underneath the chair seat, so they may not be visible without taking them out.




Electrical Plug (*Right side only)




Air hoses
(Both left and right side)

Remove the rubber plugs located on the inside of the arm panels, and set aside for later.




The armrest exterior parts should snap into place.





Installing the Armrests

1. Carefully angle one of the arm panels so that the bracket on the body of the chair is angled with the brace on the inside portion of the arm panel. Then proceed to slide it halfway into the mounting bracket.



STEP 3: INSTALLING THE ARMRESTS



2. Holding the arm panel in place halfway, reach and connect the air hoses to the nozzles located on the arm panel.




(The grey hose connects to the right nozzle while the black hose connects to the left nozzle respectively.)

(The electrical plug connects to the socket on the right arm panel. There is no electrical plug and socket on the left side.)

3. Then align the panel with the body of the chair in order to slide the arm panel into the slots.




4. Tighten the screw fully in order to secure the armrest to the chair. Once fully screwed in, plug the hole with the rubber plug that was previously removed.



Rubber plug


5. Repeat Steps 1-4 for the another arm panel.



STEP 4: INSTALLING OTHER ACCESSORIES


Installing the backrest, additional head cushion and headrest

The backrest can be secured with zippers along the length of the chair. Then place the additional head cushion at the top of the chair, which is also fastened with zippers. The headrest is secured by sliding the tab in between the open handle at the top.



Installing the controller and controller pocket


Connect the remote controller pocket to the right armrest by the velcros. The remote control can be plugged into the right side of the chair, near the air hoses and electrical plug.



STEP 5: PLUG AND PLAY

Power cable can be plugged into the back of the chair and the power can be turned on by flipping the power switch on. You can now enjoy your RELAXONCHAIR Yukon-4D Massage Chair.


Note: Please refer to the user manual for instructions on how to operate the massage chair.



RELAXONCHAIR

YUKON-4D

USER MANUAL



Please read SAFETY INSTRUCTIONS carefully before using the chair.

YUKON - 4D

Function Introduction

1. With the addition of patented 4D massage design built into the roller mechanism, the massage chair is able to deliver a powerful, creative massage that utilizes techniques such as kneading, tapping, a combination of kneading and tapping, shiatsu, knocking, and rolling in order to imitate human hands along your body.
2. Heating has been implemented into the rolling mechanism with a heating ball, allowing the heat therapy to reach all the way from the neck to the buttocks.
3. Long L-track designed to apply massage techniques all along the length of the body from your head to your buttocks.
4. Space-saving function allows for the massage chair to be close to the wall, as the chair slides forward as the back reclines.
5. Zero-gravity functionality is available with just a-click-of-a-button. It relieves pressure off of your spine and other joints by achieving a flat posture in order to feel completely relaxed.
6. A simple hand controller with full-color LCD screen for displaying massage programs and status in real time.
7. Wireless Bluetooth connectivity allows you to play music or videos from your device, and lets you set the mood in order to match your style.
8. Equipped with omnibearing rolling and air compression technology, the foot massager is able to relax your feet in addition to a roller heating function.
9. The calves are kneaded with multiple air cells by using pulsing air pressure which can promote blood circulation and restore vitality to your legs.
10. Arms feel relaxed from the air compression massage, which can promote blood circulation.
11. Experience a full body compression massage with airbags located on the waist, hip, thigh, arms, calves, feet, back, and shoulders to relieve muscle tension.
12. The footrest can be extended up to 8.7 in, accommodating for various heights.

Thank you for choosing the Relaxonchair Yukon Massage Chair. Please read this manual very carefully before use of the product. Keep this manual for future reference. Please pay attention to the SAFETY INSTRUCTIONS.

Note: We reserve the right for the design to be changed. It is subject to change without notice. The pictures used in this manual are just for reference, and we ask you to please look at the actual product.

Table of Contents

PREPARATION

| | |
|--------------------------------|---|
| Safety Instructions | 1 |
| External Structure Parts | 5 |
| Internal Structure Parts | 6 |
| Mechanical Components | 7 |

INSTALLATION

| | |
|--------------------|---|
| Installation | 9 |
|--------------------|---|

MASSAGE OPERATION

| | |
|------------------------------------|----|
| Adjustments Before Massaging | 14 |
| Manual Control | 17 |
| Operation Instruction..... | 22 |
| Function Illustrations | 23 |
| Voice Recognition | 25 |
| After Massage | 26 |

MAINTENANCE



| | |
|--------------------------------|----|
| Cleaning and Maintenance | 27 |
|--------------------------------|----|

OTHERS





| | |
|----------------------|----|
| Malfunctions | 28 |
| Technical Data | 29 |

1 SAFETY INSTRUCTIONS



- Please read this manual carefully to ensure that the installation and operation of the massage chair is correct.
- Warning and Caution markings are provided for important safety information.





| | |
|---|---|
|  Warning | This sign indicates possible cause of death or severe injury to the body. |
|  Caution | This sign indicates possible injury or damage to things. |

SYMBOL MEANING






| | | | |
|---|---|---|---|
|  | This sign indicates forbidden behavior. |  | This sign indicates that disassembly is forbidden. |
|  | This sign indicates instructions for operating. |  | This sign indicates the correct operation of a cable. |





- The following should be kept for any further references.

|  Warning | | |
|---|---|--|
|  | <ul style="list-style-type: none"> • People with osteoporosis, vertebral fractures, and acute painful diseases such as sprains or muscle strains. • Heart disease patients (including those who have pacemakers or are using other medical electronics). • Cancer patients. • People with perceptual disorders due to high levels of peripheral circulatory disturbances such as diabetes. • People with abnormal or crooked spine. • People with Protruded intervertebral Disc. • Pregnant women and those who have just completed childbirth, as well as people with poor health. • People with trauma, such as thrombosis, acute varices, and various skin inflammations. • People who have been identified as having lost their ability to discern temperature. • People with high fever (acute inflammation: such as burnout, chills, blood pressure changes, physical weakness). • Elderly people over 80. • Persons with impaired consciousness, mentally handicapped persons, or those who are incapable of operating the unit are prohibited from using this unit. • People who can not take care of themselves, infants, and drunkards. • Do not use this unit when it is wet or when you are wet. • It must be taken care of, or not used by people who are not sensitive to heat or people who are allergic to heat. • Do not fall asleep in this unit. <p>People with above symptoms should not use massage chair, or use it under doctor's advice. User's discretion is advised if use the chair.</p> | <p>Not suitable users.</p> <p>Risk of causing accident or body discomfort.</p> |

|  Warning | | |
|--|--|-----------------------------------|
|  | <ul style="list-style-type: none"> ● Before sitting, please confirm that there are no foreign matter inclusions between the various parts of the machine. ● When using the backrest and footrest, be sure to check that there are no people or objects in the area (behind, under, or in front of the chair). ● In the event of a physical abnormality or abnormality in use, stop using it immediately and ask a doctor. ● Do not over massage the head. ● Be careful when you massage the area around your neck. ● Avoid excessive neck massage. ● Do not massage the abdomen or knees with this machine. ● If you want to leave your seat during the massage, you must stop the massage chair before you leave the seat. ● Do not use it with other instruments. ● Do not allow children or pets to play on or around the unit (especially in use); nor can they sit or stand on backrests ,armrests or legrest. ● When people are inconvenient for action or lack of communication skills, they must have strict and continuous care. ● Before moving the unit, the backrest and legrest must be angled back to their original position. | Accident /injury |
| | <ul style="list-style-type: none"> ● When assembling, turn off the switch and unplug it after the chair returns. ● Insert the plug of the chair into correct and suitable power socket. ● Make sure the plug is completely inserted into the socket to reduce short circuit or fire risk. ● Unplug the power when not in use,stop all massage functions and turn off on/off key before unplugging. | Fire and electric shock |
|  | <ul style="list-style-type: none"> ● Do not fold the massage chair cushions. ● The machine is forbidden for no-load running for a long time. ● An appliance should never be left unattended when plugged in. | Accident |
| | <ul style="list-style-type: none"> ● Do not use the chair if any of the cables and plugs are frayed. ● Do not use other electric power except rated voltage in case of fire, electric shock or any breakdown. ● Any broken, re-process, knotted, twisted power cord or pressed under heavy cargo is not allowed to use. ● Never operate this appliance if it has a damaged cord or plug, if it is not working. ● Keep the chair away from heat sources. | Fire and electric shock |
|  | <ul style="list-style-type: none"> ● If the chair works incorrectly or is broken, splashed with water or falls into water, please have it repaired by professional repairman or send the chair to assigned service agent. Self-modification, disassembly or maintenance is forbidden. ● If power cord or plug is damaged, please don't replace it yourself, please contact us or appointed dealer for replacement. ● To avoid electric shock or bodily harm, please don't open or disassemble any part of the chair, there is not an attached part of the chair for maintenance by customers. ● Check if the leather or cloth material on the back and other areas is damaged or not. If broken, even though it is small, please stop using the chair, unplug the power and have it repaired. | Disassembly and maintenance |


3 SAFETY INSTRUCTIONS

|  Caution | | |
|--|--|--|
|  | <ul style="list-style-type: none"> ● Please operate the chair in a flat or level place. | To avoid sudden falling. |
| | <ul style="list-style-type: none"> ● The maximum suggested operating time is 20mins. 5mins for each position beginning with the slight- step massage. A long continuous massage in the same position may cause an uncomfortable feeling. | Too long time massage is bad for your body. |
| | <ul style="list-style-type: none"> ● Before sitting on the chair to avoid injury please inspect everything to confirm that the massage heads are in the correct position. | To avoid accidental injury. |
| | <ul style="list-style-type: none"> ● Remote controller should be timely put into the remote controller bag. Do not arbitrarily put away, to avoid squeezing or drop damage. | |
| | <ul style="list-style-type: none"> ● Please make sure that the power switch is at the“OFF” position before inserting the plug into the socket. | |
| | <ul style="list-style-type: none"> ● Turn off all controls before removing plug from outlet. | To avoid the accident when turn on the power. |
| | <ul style="list-style-type: none"> ● Make sure the plug completely insert into the socket. | To reduce short circuits and any fire hazard. |
| | <ul style="list-style-type: none"> ● Checking whether any damage on wire,make sure no water or rubbish on wire. | To avoid electric shock in any leakage or breakdown. |
| <ul style="list-style-type: none"> ● Installing the plug in the correct way. The plug must be inserted into socket with plug, the plug couldn't be bend or cut off.  | | |
|  | <ul style="list-style-type: none"> ● Do not press with hands or step on the massage heads with feet. ● Do not put hands between the working massage heads. ● Do not put hands or feet in the gap of mechanical parts. ● Do not drop or embed anything into the slot of massage chair. ● Do not stand or sit or put anything on the backrest, armrest or leg-rest. ● Do not move the chair during operation. ● Do not stand on the chair during operation. ● Do not pull out the plug or turn off the power during operation. ● Do not put your hands or head between leg-rest or beneath. ● Do not use the massage chair while wearing hair ornaments. ● Do not put hard things in pocket when take buttocks and huckle massage. ● Do not sit or press on the backrest when the chair is in reclining position. ● Do not use on bare skin or thin clothes. Thin clothes might achieve a better massage effect, but can also cause irritation to the skin or burns. ● Do not use massager with loose clothing or jewelry. ● Keep long hair tied or tucked away from massager while in use. | To avoid accidents or injury |
|  | <ul style="list-style-type: none"> ● Make sure to install the plug in correct way (Incorrect installation way may cause electric shock. If you don't know how to install, please find a professional electrician.) ● Do not use plug adapter (If the enclosed plug doesn't match with socket, please don't change the plug by yourself and find a professional electrician to install a matched socket.) | To avoid electric shock |

|  Caution | | |
|--|--|--|
|  | <ul style="list-style-type: none"> ● Keep the chair away from hot, wet, sharp, corrosive, inflammable, and explosive environments. ● To use electric blanket or other heating products on the chair is forbidden. | To avoid fade, discolor on upholstery or hardened, deteriorated on PU. |
| | <ul style="list-style-type: none"> ● Keep the chair away from stove or other heating products. ● Do not put the chair directly under the sun or other high temperature places. ● Do not put a dusty space or corrosive gas environment. ● Do not put the chair in a small space or use in a bad air circulation place. | |
| | <ul style="list-style-type: none"> ● Don't use the chair outdoor. | |
|  | <ul style="list-style-type: none"> ● Keep Dry-Don't use or store the chair in the wet environment like nearing the swimming pool, or bathroom. ● Keep the chair and remote control dry all the time. ● Do not use it when user's body is wet. | To avoid leakage or electric shock. |
|  | <ul style="list-style-type: none"> ● Unplug the power before cleaning the chair. ● Unplug the power after using the chair. | To avoid electric shock or fire hazard. |
| | <ul style="list-style-type: none"> ● Unplug the power when any damage on the power. | |
| | <ul style="list-style-type: none"> ● Hold the plug to unplug the power but not try to pull the wire to unplug the power. | To avoid accidents. |
| | <ul style="list-style-type: none"> ● Please turn off the power switch and unplug the power after massage. | To avoid damage or accidents when power on again. |
| | <ul style="list-style-type: none"> ● Unplug the power immediately when electricity goes off. | |

GROUNDING INSTRUCTION

The chair is equipped with a cord with a built-in grounding conductor. The plug must be plugged into an appropriate socket in order to ensure that it is properly installed and matches the plug.





REPLACING THE FUSE

When the fuse has burned out, please turn off the power switch located at the back of the chair. Open the lid of the cover, which is covering the switch box. Replace the old fuse with the new one. Then place back on the cover.


5

EXTERNAL STRUCTURE PARTS




INTERNAL STRUCTURE PARTS

6



7 MECHANICAL COMPONENTS

MECHANICAL COMPONENTS DIAGRAM



MECHANICAL COMPONENTS

MECHANICAL COMPONENTS LIST

| NO | Name | Qty | NO | Name | Qty |
|----|------------------------------|-----|----|-----------------------------|-----|
| 1 | Remote control | 1 | 12 | Left armrest assembly | 1 |
| 2 | Right armrest assembly | 1 | 13 | Left armrest frame assembly | 1 |
| 3 | Massage machine | 1 | 14 | Footrest electric actuator | 1 |
| 4 | Power case assembly | 1 | 15 | Seat frame surface assembly | 1 |
| 5 | Left side fender | 1 | 16 | Seat bottom frame assembly | 1 |
| 6 | Right armrest frame assembly | 1 | 17 | Pillow | 1 |
| 7 | Under frame assembly | 1 | 18 | Shoulder pad | 1 |
| 8 | Seat electric actuator | 1 | 19 | Back and seat cushion | 1 |
| 9 | Backrest electric actuator | 1 | 20 | Backrest assembly | 1 |
| 10 | Right side fender | 1 | 21 | Shoulder device assembly | 1 |
| 11 | Foot and leg assembly | 1 | 22 | Back cover assembly | 1 |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |
| | | | | | |


9 INSTALLATION

STEP 1: PREPARING ASSEMBLY

- Check the Boxes

Check all the boxes to make sure that you have everything you need. Make sure that the chair is reclined and is away from objects and obstructions for safe assembly.

BOX 1 (MAIN CARTON BOX): MASSAGE CHAIR MAIN BODY & 2 ACCESSORY BOXES




Massages chair main body

Accessories from the accessory boxes


- Headrest
- Additional Head cushion
- Backrest
- Controller pocket
- Electric power cable
- The armrest exterior parts
- Controller & Cable

BOX 2 (ARMREST CARTON BOX): LEFT AND RIGHT ARMRESTS



Left armrest Right armrest

BOX 3 (FOOTREST CARTON BOX): FOOTREST




Footrest

- Check the environment

Check all the boxes to make sure that you have everything you need. Make sure that the chair is reclined and is away from objects and obstructions for safe assembly.

Due to the weight of the chair, we recommend using a blanket or soft cover in order to protect your floor from getting ruined.


| | |
|---|---|
|  Caution | <ul style="list-style-type: none"> • Don't use the chair in wet environments like near the bathroom or swimming pool to avoid leakage or electric shock. • Please use the chair in a flat area to prevent the chair from tipping over, falling, creating noises or other unexpected trouble. • Don't use the massage chair under direct sunlight or near heating places like the stove or fireplace in order to avoid the upholstery from deteriorating. |
|---|---|

STEP 2: INSTALLING THE FOOTREST


- **Preparing to Install the Footrest**

Locate the body of the chair on a soft surface in order to begin assembly. Identify the location of the footrest screws and attachments on the body of the chair. Using a Phillips-head screwdriver, remove the footrest screws and put aside for later.


1. Place the footrest near the body of the chair, and unscrew the the footrest




2. Plug in the electrical cord and air hose into their respective socket and




3. Lift up the footrest attachment covers, then carefully lift up the footrest to place the footrest into the footrest attachments.



4. Close the footrest attachment covers, and screw the footrest screws



5. Remove the velcro covering from the footrest flap, and attach the velcro to the body of the chair.




11 INSTALLATION

STEP 3: INSTALLING THE ARMRESTS


• Preparing to Install the Armrests

Locate the left-side air hoses on the body of the chair that will be connected to the left-side arm panel. Locate the right-side air hoses and electrical plug on the body of the chair that will be connected to the right-side arm panel.


Note: The hoses and electrical plug may be underneath the chair seat, so they may not be visible without taking them out.




Electrical Plug (*Right side only) Air hoses (Both left and right side)




1. Remove the rubber plugs located on the inside of the arm panels, and set aside for later.




2. The armrest exterior parts should snap into place.




3. Carefully angle one of the arm panels so that the bracket on the body of the chair is angled with the brace on the inside portion of the arm panel. Then proceed to slide it halfway into the mounting bracket.




4. Holding the arm panel in place halfway, reach and connect the air hoses to the nozzles located on the arm panel.




5. Then align the panel with the body of the chair in order to slide the arm panel into the slots.



6. Tighten the screw fully in order to secure the armrest to the chair. Once fully screwed in, plug the hole with the rubber plug that was previously removed.




7. Repeat Steps 1- 6 for the another arm panel.



STEP 4: INSTALLING OTHER ACCESSORIES

- **Installing the backrest, additional head cushion and headrest**


The backrest can be secured with zippers along the length of the chair. Then place the additional head cushion at the top of the chair, which is also fastened with zippers. The headrest is secured by sliding the tab in between the open handle at the top.



13 INSTALLATION

- **Installing the controller and controller pocket**


Connect the remote controller pocket to the right armrest by the velcros. The remote control can be plugged into the right side of the chair, near the air hoses and electrical plug.




STEP 5: PLUG AND PLAY

Power cable can be plugged into the back of the chair and the power can be turned on by flipping the power switch on. You can now enjoy your RELAXONCHAIR Yukon-4D Massage Chair.

Note: Please refer to the user manual for instructions on how to operate the massage chair.




CHECK ENVIRONMENT

| | |
|---|--|
|  Caution | <ul style="list-style-type: none"> • Don't operate the chair in wet environments in order to avoid leakage or electric shock. • Don't operate under direct sunlight or near heating sources in order to avoid upholstery deterioration. • Please operate the chair on a flat surface in order to avoid falling, injuries, unnecessary noise, etc. |
|---|--|


• Check Places in Front

Check the areas around the chair before operation in order to make sure that no people, pets, or objects are in danger of being injured/damaged.




• Check Power Cable

Ensure that the cable is not frayed, twisted, tied, or underneath a heavy object in order to have safe operation.




• Space Size for using the chair

The chair requires at least 4 inches from the back of the chair to the wall in order to accommodate for the full recline of the chair.



• Check Surroundings

Pay attention to any person or things within the vicinity of the chair and for ample space when the backrest and footrest reclines and extends.



The chair requires at least 24 inches in front of the leg in order to accommodate for the full recline of the chair.

15 ADJUSTMENTS BEFORE MASSAGING


HOW TO MOVE MESSAGE CHAIR

⚠ Caution

• Don't move the chair while it is in operation in order to prevent accidents or damages.

• Preparing to Move


- Bind the cables and put them in the seat to avoid damage when moving.
- Do not use the roller on the back of the chair if the flooring can be easily damaged or ruined.
- It may take 2 or more people to properly move the chair.
- Check to make sure that the power has been turned off, and that the power cable is away.



• How to Move

Press down on the head portion to create leverage on the rest of the massage chair. Using the rollers on the back of the chair, move the chair to your desired location slowly and steadily. Then proceed to set down the chair slowly.

The footrest may need to also be held up as shown in the picture in order to prevent scratches and scrapes on the leather/fabric.




ADJUSTMENTS BEFORE MASSAGING 16

MASSAGE ADJUSTMENT

• Massage Adjustment


Make sure the back rollers stop at the upper portion of the back, then proceed to sit on the chair slowly.

When the massage chair cannot find the shoulders during body scans, it will default to the middle of the back portion.



Properly Seated:

Please rest your back up against the seat as much as possible, and make sure the head rests against the pillow.





The remote control will beep once the body scan begins into function. the up-down key can be used to adjust the shoulder height.

Following these adjustments should provide a more pleasurable experience as comfortability influences the experience.


• Additional Option
















If the massage is not strong enough for you, the headrest and head cushion can be removed.

Due to the weight of the chair, it may ruin the flooring in a long-term setting. We recommend having a blanket or something soft underneath the chair to help protect the flooring.



17 MANUAL CONTROL



| | |
|---|--|
| <p>On/Off</p>  | <p>1. The ON/OFF button powers the chair ON and OFF. While the chair is in any position and OFF, the ON/OFF button turns the chair ON, and leaves the chair in the same position. The massage chair is then able to be operated. Once finished, while the chair is ON, the ON/OFF button will turn off the chair and return the chair to an upright position. Using the ON/OFF button in the middle of operation will stop everything and return the chair to an upright position.</p> |
| <p>Heating</p>  | <p>2. This is the HEATING button. Press the button to turn on heating for the back and legs.</p> |
| <p>Pause</p>  | <p>3. The PAUSE button will stop all operations of the massage chair, and other buttons will not work. Pressing the PAUSE button again will resume the operation of the chair once again.</p> |
| <p>Menu</p>  | <p>4. The MENU button will change the display from the main interface into the menu interface. Where you can navigate through all the possible functions and options the chair has to offer. Pressing the MENU button returns the display back to the main interface.</p> |
| <p>Experience</p>  | <p>5. The EXPERIENCE button will give you a quick relaxation program that will last 8 minutes.</p> |
| <p>Auto</p>  | <p>6. The AUTO button will start the first auto program built into the chair "Recovery". Pressing the AUTO button again will cycle through to the next auto program. Keep pressing the AUTO button to cycle through all the programs. The name of the program will be visible on the screen. There are a total of 7 auto programs with an 8th program designed to stop everything.</p> |
|  | <p>7. Press and holding down this button will raise the top of the massage chair while the footrest goes down to adjust the position of the chair.</p> |
|  | <p>8. Press and holding down this button will lift the footrest up to adjust the position of the footrest.</p> |
|  | <p>9. Press and holding down this button will lower the footrest down to adjust the position of the footrest.</p> |
|  | <p>10. Press and holding down this button will lower the top of the massage chair while the footrest lifts up to adjust the position of the chair.</p> |
| <p>Zero</p>  | <p>11. The ZERO Button will adjust the position of the chair to a zero-gravity position. Pressing the ZERO button again will change the position of the chair to the next zero-gravity position until it resets to the upright position. There are a total of 3 different zero-gravity positions.</p> |
|  | <p>12. There are 4 directional buttons used for navigating the MENU button. Press UP and DOWN to quickly browse through the options, while the LEFT and RIGHT buttons are used to go from and to the previous and next menu. In manual mode, the UP and DOWN button, on the main interface, can be used for moving the roller position up and down.</p> |
| <p>OK</p>  | <p>13. The OK button is used to confirm menu options in the MENU interface.</p> |
| <p>Time +</p>  <p>Time -</p>  | <p>14. The TIME- and TIME+ buttons are used to adjust the massage times. Pressing the TIME+ button will increase the time up to the closest time up in intervals of 10 minutes, 20 minutes, and 30 minutes, but the time cannot go above 30 minutes at one time. Pressing the TIME- button will decrease the time to the closest time down in intervals of 10 minutes, 20 minutes, and 30 minutes, but the time cannot be reduced below 10 minutes.</p> |

19 MANUAL CONTROL

LCD DISPLAY GRAPHIC DESCRIPTION

• Main Interface

Pressing the **on/off** button when the power switch is flipped on at the back of the massage chair will power on the chair and the display on the remote control will light up after 2-3 seconds. The main interface shows the massage position, program, time, current massage status, etc. All information is updating accordingly on the display of the remote control.

The main interface LCD display includes the following elements:

- Speed:** Level 1, 2, 3, 4, 5 Status Display (represented by 5 bars)
- Width:** Narrow, Mid, Wide (represented by 3 bars)
- Air Intensity:** Level 1, 2, 3, 4, 5 Status Display (represented by 5 bars)
- 4D:** Level 1, 2, 3, 4 (represented by 4 bars)
- Massage Position:** A blue silhouette of a person sitting in a chair, indicating the current massage posture.
- Heating:** A wavy line icon representing the heating function.
- MassagTime Left:** A clock icon and '20 min' indicating the remaining massage duration.
- Temperature Level of Massage Machine Heat:** A thermometer icon and a bar chart showing the heat level.
- Massage Function:**
 - Kneading (represented by a red hand icon)
 - Tapping (represented by a red hand icon with a small square)
 - Knocking (represented by a red hand icon with a small square)
 - Shiatsu (represented by a red hand icon with a small square)
- Zero Gravity:** A red wavy line icon with the text 'Zero Gravity'.
- Instructions:** 'Set function by pressing Menu key.'

• Menu Interface

The MENU button will switch the remote display form the main interface to the menu interface. Pressing it again returns you to the main interface.

The menu interface LCD display shows the following options:

- Auto** (selected)
- Manual
- Airbag
- Roller
- Heating
- Timer
- Settings
- Stop
- Recover
- Relief
- Stretch
- Soft
- Renew
- Calming
- Special

The **Special** option is expanded to show sub-options:

- Refresh
- Energy
- Revive
- Active
- Relax
- Healing


Instructions: 'Select by the navigation keys'

• Auto Mode:

Press the MENU button, select the auto program you enjoy by pressing the OK button. The list of programs are **Recovery, Relief, Stretch, Soft, Renew, Calming, and Special**. The Special program has 6 different programs under it called **Refresh, Energy, Revive, Active, Relax, and Healing**.


Manual Mode:

Press the MENU button, and select the manual option using the directional keys. There are 6 options: **Mode**, **Position**, **Width**, **Speed**, **4D**, and **Footrest**. You can choose the corresponding options with the OK button. Mode has 5 options: **Kneading**, **Tapping**, **Knocking**, **Kneading & Tapping**, and **Shiatsu**. There are 6 positions available: **Full body**, **Shoulder**, **Back**, **Waist**, **Seat**, and **Fixed Point**. The massage width and speed is able to be changed. The 4D function can be adjusted to different levels of intensity. The footrest could be extended outwards and retracted inwards with the **In** and **Out** options.




Airbag Massage:

On the menu interface, the airbag option allows you to choose the position and intensity of compression massages. The **Full Body**, **Shoulder**, **Hand**, **Seat**, **Leg**, and **Seat**. The intensity can be adjusted between 5 different levels.



Foot Roller:


Select the **Roller** option with the navigation key. Choose between 3 options: **Quick**, **Soft**, or **Stop**.



21 MANUAL CONTROL


Heating:

Press the MENU button, and navigate towards the **heating** option. The **Back** and **Foot** options can be toggled ON/OFF, which heats up the back and foot portions. If the heating is too hot or not hot enough, the TEMP+ and TEMP- buttons can be used to control the temperature of the heating function.




Timer:

Navigate the MENU to the **Timer** option. The time can be adjusted to intervals of 10 minutes. You can choose between 10 minutes, 20 minutes, and 30 minutes. After the time is up, the chair will automatically turn off.



Setting:


After pressing the MENU button, you have the ability to set the language to either Chinese or English.



* When the chair is powered on, you can start the auto program. The auto program starts with a body scan, and will be displayed on the screen. After the body scan is finished, a beeping sound will beep 5 times over 5 seconds. At about this point in the auto program, the program will be at the shoulders. At this point, the shoulder height is adjustable up and down with the navigation keys. The default massage time is 20 minutes, and will turn off the chair once the massage is over.

Stop:


Press MENU key, select stop through the navigation key, it will stop massage, and get back to main interface.



ON/OFF

• Step 1. Power On

1. Put one end of the power cable to the massage chair.
2. Put the power cable into the socket.
3. Proceed to turn on the power switch located at the back of the massage chair.
4. Press the “ON/OFF” button




• Step 2. Massage



1. Operate according to the controller operation instructions. Please refer to the remote control page.
2. The massage chair will return to the upright position once the massage ends.



• Step 3. Turn off the power

1. Press ON/OFF button while the chair is on to return it to the upright position.
2. Turn off the power switch at the back of the base of the chair.
3. Unplug the power cable from the socket.



| | | | |
|---|--|--|--|
|  Warning | <ol style="list-style-type: none"> 1. Before use, check to make sure that PU material is not damaged or broken. If it is broken or damaged, please stop use immediately, unplug the chair, and ask for repairs. 2. Only plug the chair into a properly working power socket. |  Attention | <ol style="list-style-type: none"> 1. Before sitting down, make sure that nothing is being pinched between the areas of the armrest, footrest, backrest, and seat of the chair. 2. The footrest is equipped with a safety protection system, so that if there is anything between the leg and foot portion of the footrest, it will stop retraction automatically. 3. Make sure that the chair is in the upright position before using the chair. 4. Do not stand on the chair while using it. |
|---|--|--|--|

23 FUNCTION ILLUSTRATION


• Footrest Up

Press and hold this button on the controller to lift up the footrest. Releasing the button stops moving the footrest, and locks the position of the footrest.



• Footrest Down

Press and hold this button on the controller to lower the footrest. Releasing the button stops moving the footrest, and locks the position of the footrest.



• Footrest Up & Backrest Down

Press and hold this button on the controller to lift the footrest and lower the backrest. Releasing the button stops moving the footrest and the backrest, and locks the position of the footrest and backrest.



• Footrest Down & Backrest Up

Press and hold this button on the controller to lower the footrest and lift the backrest. Releasing the button stops moving the footrest and backrest, and locks the position of the footrest and backrest.




• Zero Gravity

Press the ZERO button on the remote control to adjust the position of the footrest and backrest into a zero gravity position. This massage position is used to relieve pressure on the spine and joints. There are 3 different stages, and can be cycled until the chair returns back to the upright position.



*Tips: Zero-gravity, according to the NASA space capsule relax theory, is to make the legs and heart be on the same elevation, which disperses the weight and pressure from the body and legs to the rest of the body, where each part is equally affected by the effects of gravity. This achieves a comfortable, relaxing massage posture.




Warning

1. Please make sure that there is nothing within the vicinity of the chair that can be damaged or injured.
2. When adjusting the chair, be careful of the gap between the armrest and the seat of the chair.
3. While the footrest is moving, please do not stand or put pressure on the footrest, nor should you leave from the seat until it has completely stopped in order to prevent causing damage to the internal mechanisms of the chair.
4. Please Refer to the control instructions for more details.


- **Extending and Retracting the Footrest**

The footrest can extend and retract within 8.7 inches automatically.




- **Shoulder Massage**

Multi-layered airbags are located on both shoulders to assist with massages. When 4D massage rollers push up against your back, the shoulder airbags will grip your shoulders and hold you down in order to provide deeper massages.




- **New Curve Rail**

This massage chair has adopted 4D technology, allowing the massages to reach multi-angle integrated massage for the back and buttocks. It has the S-track built in, better fitting the contour of the human back. It also has the components of the L-track, letting the mechanism reach all the way down to the buttocks and glutes. This SL-track, with the addition of 4D technology, gives your body the best massage possible, hitting all of the right notes.




- **Arm Massage**

The massage chair uses a capsule type armrest with a built-in air pressure massage function. Put the arms in between the armrest's airbags. After filling up with air, the arms and hands are kneaded by rhythmic air compressions for comfortable massages.




- **Warm Foot Therapy**

An innovative two-way rotary foot kneading technique with the addition of foot heating allows the footrest to give masterful, relaxing foot massages.



- **Music Therapy**

This massage chair takes multiple speakers and the mechanism in order to rhythmically massage you to the beat of your music. The music allows you to feel comfortable and lets you enjoy a massage with accompanying music.



25 VOICE RECOGNITION

4D VOICE RECOGNITION COMMAND

- To Turn On Voice Command

Say 'Turn on speech recognition mode'.

* If your voice command is recognized, the speaker will respond after 1 second.

- Commands

You may now give your chair any of the commands below within 30 seconds. Make sure to speak clearly and slowly.

- Chinese Massage
- Pain Relief
- Body Stretch
- Experience
- Relaxation
- Zero Gravity Start
- Massage Chair Shutdown

- To Turn Off Voice Command

At any time after 30 seconds of your last command, the Voice Command function will turn off automatically.


To turn off manually, simply say 'Voice off' to turn voice recognition off.

* If you want to activate a different command, it must be awakened again by saying the activation phrase, 'Turn on speech recognition mode'.


| Voice recognition command | |
|---|--|
|  Turn on speech recognition mode | |
|  Chinese Massage |  Pain Relief |
|  Body Stretch |  Experience |
|  Relaxation |  Zero Gravity Start |
|  Voice Off |  Massage chair shutdown |

RESTORATION OF MASSAGE CHAIR


1. Make sure that there is nothing within the vicinity of the chair's backrest and footrest.




2. The backrest and footrest can be restored to the upright position by pressing the ON/OFF button on the chair to turn off the power.




3. After the chair returns to the upright position, store the remote control in the remote control pocket on the right armrest.



4. Switch of the power switch located at the back of the base of the chair.



5. Unplug the power cable from the electrical socket by tugging on the base of the plug and not the cable.



6. Clean the chair if necessary.



Attention

- After using, press the ON/OFF button to end all operations on the chair, and please wait for the chair to fully stop moving before getting out of the chair.
- Make sure to unplug the electricity. Avoid moisture, dirt, etc. in case of causing damage.

27 CLEAN AND MAINTENANCE

SYNTHETIC & GENUINE LEATHER (SITTING AREA)

- **General Cleaning**

Clean with tidy, soft, and dry cloth slightly.

Note: Do not use thinner, gasoline, alcohol or other chemicals.

- **If the leather is dirty**

1. Dip the soft cloth into a 3-5% neutral detergent solution, wring out most of the liquid, then clean the stain on the surface very gently.
2. Wash the cloth, squeeze it out, then wipe with just water.
3. Clean the area with a soft dry cloth.
4. Allow for the area to dry naturally.

CLOTH MATERIAL

1. Dip cloth into weak neutral detergent, squeeze out most of the liquid, and gently wipe down the area.

Note: Do not use thinner, gasoline, alcohol, or other chemicals.

2. Brush down the area with a weak neutral detergent.

Note: excessive cleaning may cause damage to the material.

3. Wipe with water-dipped but squeezed cloth.
4. Allow to dry naturally.

PLASTIC PARTS

1. Dip cloth into weak neutral detergent, squeeze out most of the liquid, and wipe the area.
2. Dip cloth into water, squeeze out most of the liquid, and wipe down the area.
3. Allow to dry naturally.

CLEAN THE CONTROLLER

1. Wipe with soft wet cloth that has most of its liquid squeezed out.
2. Allow to dry

ABOUT CLEANING

Before cleaning, unplug the electricity. Do not unplug something electrical with wet hands. Make sure that your hands are dry when plugging the chair back into the socket.

Note: Failure to follow the cleaning instructions may get you shocked or burned.

FAQ AND TROUBLESHOOTING

| Trouble | Possible reason | Troubleshooting |
|---|---|--|
| The massage chair can't work after power is on. | Check if the remote controller is on or not. | Open control switch |
| | Not choosing a massage function. | Choose massage function |
| | The plug and socket is in poor contact with the chair. | Make sure plug and socket are connected. |
| | Wire or plug has damaged. | Please ask manufacturer or appointed or authorized maintenance man for checking and replacement. |
| | The fuse burned out. | Replace fuse of same specification. |
| | Something wrong with inner circuit | Please ask appointed or authorized maintenance man for checking and repairing. |
| There is a beeping sound inside the massager while working. | It is sound of air pump, motor and mechanical parts. | Normal sound, no need to pay attention to it. |
| Different height of the left roller and right roller. | Massage rollers work with oscillations. | Normal, no need to pay attention to it. |
| Sound suddenly becomes much louder while working. | Check if it is due to long time working. | Turn off the power to have the massager rest for over half an hour before next turn working. |
| | Inner parts worn away due to age. It's normal. | Please ask appointed or authorized maintenance man for checking and repairing. |
| Massage function suddenly stops while working. | Suddenly power-down | Check the electricity power again |
| | The session time is up. | Turn off the power to have the massager rest for over half an hour before using again. |
| Neither the backrest nor the legrest can be lay down. | Check if there are any obstacles; the chair is overloads; for the sake of safety, sometimes the massage machine may automatically stop working. | Remove the obstacles; turn off the power to have the massager rest for over half an hour before next turn working. |
| The chair can't return to its original position. | Excessive operation; setting time is up. | Please ask manufacturer or appointed or authorized maintenance man for checking and replacement. |
| Electric wire or plug is heating up. | It may work overload | Turn off the power to have the massager rest for over half an hour before next round working. |
| | Other abnormal case | Stop using, and ask appointed maintenance man for help. |



Note: If the troubles are not caused by above reasons, please turn off the power switch and cut off the power, then contact the manufacturer or service agent to repair your chair.

29 TECHNICAL DATA

TECHNICAL DATA

| | | | |
|---------------|----------------|-------|---|
| Name | Massage Chair | Power | 260W |
| Model | Yukon - 4D | NW/GW | Main body: 236/283 lbs Armrest: 38/48 lbs Legrest & footrest: 70/78 lbs |
| Input voltage | □ AC120V, 60Hz | | Packing size |
| Safety design | Class I | | |