

# User Guide

1

Thank you for purchasing your Bristl! We'll help you maximize your oral hygiene. Bristl comes with cutting-edge wavelength light therapy technology backed by dentistry clinical studies to help prevent and treat gingivitis, strengthen gums and make teeth whiter. To fully benefit from the product, please read the important safeguards and operation instructions.

Please read all safety instructions before use.

Important notices:

The product is not intended for use by persons with an oral sensory disorder or reduced oral motor skills. The product is a personal oral health care device, and it is not intended for sharing among multiple users.

Consult your dental professional before using the product if:

- you have had an oral surgery within last 60 days
- you are currently receiving gum treatment or
- you have one or more dental implants

Do not attempt to disassemble the product as it is not intended for user repair and contains factory serviceable only parts.

Keep the product and loose brushes out of reach of children.



2

## CONTENTS

1. Brush cap to cover the brush for traveling and storing.
2. Replaceable brush with microfiber bristles.
3. Built-in brush head mount with sonic vibration and light emitters.
4. Power on/off button with power indicator and mode operation.
5. Main body with ergonomic grip design.
6. Micro USB charging cable.
7. Wall charging block with USB port

3

## CHARGING YOUR BRISTL

1. Connect the given micro USB cable to the charging block, and plug into standard outlet.
  2. Gently lift rubber USB charging cover to expose charging port. Insert cable.
- Power button will blink when charging.
  - A full charge lasts 3 weeks – assuming 2 minutes of brushing 3 times per day.
  - When approaching low battery, light will alternate between colors.

## Table of Contents

1. Notice
2. Contents
3. Charging
4. Operation Instructions
5. Removal/Cleaning of Brush Head
6. Brushing with Bristl
7. Dangers and Warnings
8. Brush Head Replacement
9. Auto-Timer
10. Warranty
11. Storage
12. Specifications

4

## OPERATING YOUR BRISTL

1. Press the power button once to turn on. Default mode is regular vibration with both red and blue light activated.
2. One quick press toggles vibration modes between regular, gentle and no vibration mode.
3. Press for two seconds and release to toggle between colors.
4. Press and hold for three seconds to power off.

## Bristl's SONIC VIBRATION MODES

Regular mode (default): For deep cleaning – 24,000 RPM  
Gentle mode: For sensitive area cleaning – 12,000 RPM

5

## CLEANING BRUSH HEAD

Detach brush head from the brush (see below for instructions). Rinse the brush head with clean water. Use a cotton swab to clean the inside if necessary. Clean once a week for optimal results.

- Make sure USB port cover is closed when cleaning.
- Do not use soap water or cleaning solution.
- Do not pick inside the brush head with sharp objects.
- Do not immerse the entire product in water.

### Detaching the brush head:

- Hold neck of brush in one hand.
- Use two fingers on your other hand to pull and slightly twist brush head until it releases.

### Attaching the brush head:

- Slide one end of the brush head base into the slot on the brush. Align the brush and push to lock.
- You should hear a click when the brush head is properly locked.

6

## BRUSHING WITH BRISTL

1. Wet the bristles and apply toothpaste.
2. Press the power button to turn on the brush.
3. Select preferred vibration strength.
4. Apply gentle pressure bristles and move the brush head back and forth across teeth and gumline.
5. While brushing, select Red Light to soothe and heal your gums.
6. While brushing, use Blue Light to whiten teeth and kill bacteria.
7. Toggle between colors to achieve desired outcomes for specific areas. Even select both colors (purple) to soothe, heal, kill bacteria, and whiten in the same area.

7

## Dangers:

- To avoid the risk of electrocution:
- Always charge the product at a safe distance from water.
- Never charge the product with a damaged cable.
- If the product falls in water while charging, unplug it from the wall immediately. In this case, DO NOT reach for the product directly.

## Warnings:

- Do not immerse the entire product in water.
- Do not charge the product while cleaning the product.
- Do not charge the product next to a heat source (e.g. stove, oven, etc.).
- Do not use the product during a shower or bath.
- Make sure the rubber USB port cover is closed when using.
- Dry the product thoroughly before charging.
- Do not point the emitting light directly into your eyes.
- Ensure the brush is securely attached to the brush head before use.
- Do not swallow the brush head.
- Do not use non-genuine attachments that do not come with the product.
- Discontinue the use of product if it appears damaged (brush head, handle, power button, USB port).

8

Bristl's microfiber brush is the only part you replace. Replace brush head every 2 months for optimum results. Use only Bristl's genuine replacement brush. Bristl's microfiber brushes are available online at [www.bristlscience.com](http://www.bristlscience.com)

9

Bristl will turn off automatically after 2 minutes. Power on again if using for more than 2 minutes

10

## WARRANTY

The product is under a 90 day manufacturer's warranty from the date of product purchase. The manufacturer will service or replace factory defects within the warranty period. This warranty will not cover damages caused by disassembling the product, damages caused by using a third party attachment, or damages caused by user negligence.

11

## STORAGE

For storing the product for an extensive period, clean the inside of the brush head cover provided. Place head cover on brush head and store in a cool, dry place. Do not store the product close to a heat source. Avoid direct contact with chemicals, heating or cooling units, and any other health hazards.

12

1. Vibration: 12,000 RPM- 24,000 RPM
2. Weight: 2.65 Ounces
3. Dimensions: 26mm (W) x 37mm (L) x 219mm (H)

For more information please visit [www.bristlscience.com](http://www.bristlscience.com)