

# BIOS living

Analog Glass Scale

Instruction Manual

SC424



## Introduction

Thank you for purchasing the BIOS Living Analog Glass Scale. This scale is equipped with an easy to read dial and has a capacity of 260 lbs / 120 kg.

## Operation

1. Place the scale on a hard, flat surface. Always weigh yourself on the same scale placed on the same surface.  
**NOTE:** If the scale is placed on a soft surface, such as a carpet, or on an uneven floor the results may be inaccurate.
2. Step onto the scale slowly and gently (do not jump onto the scale). Stand vertically without swaying and place your feet evenly apart - the dial will display your weight.

## Advice for Use and Care

- Avoid using platform when wet - platform will be slippery when wet.  
Keep scale dry.
- Avoid excessive impact or vibration to the unit.
- Clean the scale with a damp cloth and avoid water getting inside scale.
- Store scale in horizontal position.

## Five Year Warranty

If this product proves to be defective in material or workmanship within five years of purchase, please return it to the address below. It will be repaired or replaced without charge upon receipt of the unit prepaid with \$5.00 to cover handling, packaging and return postage. Please include proof of purchase, your full name, address, daytime phone number or email address.

Thermor Ltd.  
16975 Leslie Street  
Newmarket, ON L3Y 9A1  
[www.biosmedical.com](http://www.biosmedical.com)