

BIOS living

SC415 Personal Digital Scale Instruction Manual



Introduction

Thank you for purchasing the SC415 BIOS Living Personal Digital Scale. This Personal Digital Scale was developed with state-of-the-art technology and digital electronics; this precision device has a capacity of 400 lbs / 181 kg, and it equipped with an auto power off function.

Features

- Large 2" Blue LED readout
- Displays weight in lbs / kg
- Clean white platform
- Auto on
- Auto zero
- Graduation: 0.2 lbs / 0.1 kg
- 4 x AAA batteries included
- Low battery indicator
- Auto off after 8 seconds
- Accurate to 400 lbs / 181 kg

Starting your Scale

1. Remove the battery cover from the back of the scale. Insert four (4) AAA size batteries into the battery compartment according to the polarity indicated, close the battery cover.
2. Before starting the scale, change the unity of measurement. On the back of the scale there is a small switch labeled "unit", push the switch to select the desired weight.

Operation

1. Place the scale on a hard, flat surface (avoid carpet or other soft surfaces).
NOTE: If the scale is placed on a soft surface, such as a carpet, the results may be inaccurate or an error message may occur.
2. Turn ON the scale by lightly tapping on the right lower corner of the scale platform. Wait a few seconds until the display prompts to "0.0" on the screen, the scale is now ready to be used.
NOTE: The display will show "Err" if you step on the scale before it shows "0". If this occurs repeat step 2.
3. Step onto the scale slowly and gently (do not jump onto the scale). Stand vertically without swaying and place your feet evenly apart.
4. When your weight is displayed on the screen it will start to flash. The flashing will continue for 8 seconds and then the scale will turn off.

Troubleshooting

Error warnings will appear on the screen, as shown, if the following occur:

Error Displayed	Reason for Error
Lo	Indicates low battery
ERR	Indicates error during the reading
----	If the scale is overloaded - Over capacity (400 lbs / 181 kg) remove weight immediately from scale

Five Year Warranty

If this product proves to be defective in material or workmanship within five years of purchase, please return it to the address below. It will be repaired or replaced without charge upon receipt of the unit prepaid with \$5.00 to cover handling, packaging and return postage. Please include proof of purchase, your full name, address, daytime phone number or email address.

Thermor Ltd.
16975 Leslie Street,
Newmarket, ON L3Y 9A1
www.biosmedical.com