BIOSliving

Analog Glass Scale Instruction Manual SC424



Introduction

Thank you for purchasing the BIOS Living Analog Glass Scale. This scale is equipped with an easy to read dial and has a capacity of $260 \, \text{lbs} / 120 \, \text{kg}$.

Operation

- I. Place the scale on a hard, flat surface. Always weigh yourself on the same scale placed on the same surface. NOTE: If the scale is placed on a soft surface, such as a carpet, or on an uneven floor the results may be in
- Step onto the scale slowly and gently (do not jump onto the scale). Stand vertically without swaying and place your feet evenly apart - the dial will display your weight.
- 3. Recalibration:
 Adjust the pointer on the dial with your finger
 (as shown in **Figure A**) to '0' for recalibration.

Figure A

Advice for Use and Care

- Avoid using device when wet platform will be slippery when wet. Keep scale dry.
- · Avoid excessive impact or vibration to the unit.
- Clean the scale with a damp cloth and avoid water getting inside scale.
- Store scale in horizontal position.

Five Year Warranty

The Analog Glass Scale has a 5 year warranty to be free of manufacturing defects in materials and workmanship under normal applications for 5 years of the original owner. If this product becomes inoperable due to defect and requires repair, return the product with all component pieces and proof of purchase to the address listed below. This warranty does not cover any shipping/transport costs. This warranty does not apply if the product is subject to misuse, neglect, rough handling or damage.

Ship the unit prepaid and insured (at owners option) to:

Thermor Ltd.

Attn: Repair Department 16975 Leslie Street Newmarket, ON L3Y 9A1

www.biosmedical.com Email: support@biosmedical.com



BIOS I Medical

MANUFACTURED BY / FABRIQUÉ PAR:
THERMOR LTD.

16975 LESUE STREET
NEWMARKET, ON LBY 9A1

MADE IN CHINA / FABRIQUÉ EN CHINE
WWW BIOSAFIDICAL COM