

User Manual

Vibrant™ Back

Y3MV



Thank you for choosing Vibrant™.

For customer support or to learn more about our products, please contact us at:

Email: support@vibrant.me Neurovative Technologies Inc.
 Toll-free: 1-855-584-2726 656 Colby Drive, Unit 201
 Local 519-888-8184 Waterloo, ON, N2V 1A2

Before using the Vibrant™ Back, please read this manual thoroughly and retain for future reference.

About Your Device:

Vibrant™ Back is non-invasive solution for pain relief. Benefits may include:

- Pain reduction
- Increased range of motion
- Less reliance on medication
- Improved overall quality of life

This package contains:



● Charging cable




● USB adapter

First-time use:



Vibrant™ Back comes ready to use for your first treatment, remove your device from its package and follow the instructions below.

Using your Vibrant™ Back:



Turn the device ON by pushing the finger sized button located at the bottom of your device. An initial vibration indicates that your device is ON and ready. After a short pause your 16 minute 3s treatment will start.

Laying down or on a comfortable chair place your Vibrant device on the are where you feel pain.




During the 16 minutes treatment, your device will alternate between 4 pre-defined patterns. Alternatively, the intensity and duty cycle settings can be adjusted through the Vibrant™ App in your smart phone or tablets. Please note this functionality is only available in Vibrant BLE versions.

At the end of the treatment, the vibration stops automatically, and your device goes into sleep mode.

To extend the battery life turn the device OFF by pressing the ON/OFF button.

Charging instructions:




- Connect magnetic connector to the charging connector at the bottom of the device.

Insert charger cable to the AC adaptor and plug the adaptor into an electrical outlet.



Once successfully connected:



Red light ON: device is charging

- Green or Blue Light ON: device is charged

Please note that charging may take up to 90 min to complete if the battery is very low.

Using the Mobile Application:

Download the free mobile application to your smartphone or tablet and unleash the power of Vibrant™! Search for "Vibrant" in the Google Play Store, Apple App Store or Blackberry App World to download the App for free. Please note: Vibrant™ products do NOT require the app to work effectively, however we strongly recommend using the app for customized treatment. The Vibrant™ app allows you to create different profiles and personalize your treatment by adjusting the vibration patterns, speed and intensity. Please also note that **Vibrant™ app only works with Vibrant™ BLE version products.**

Connecting your device to the Application:

Please ensure that Bluetooth is enabled on your mobile device or tablet. After downloading and opening the Vibrant™ application on your mobile device, follow the instructions for setting up a user profile. Turn your Vibrant™ device ON. Once your device is ON, select "Begin Treatment" on the App and your Vibrant™ device will automatically connect to your mobile device.

For support please contact us at:

Tel: 1-855-584-2726

Email: support@vibrant.me

