

frasaco

Strip Crowns



Frasaco strip crowns are elastic, thin and easy-to-cut celluloid crowns. In the transparent execution, they are light-transmissive for polymerization light. They are used as moulds for creating temporary restorations on preparations and as matrices for restoration of corners and edges. They are placed without problem using any conventional material for temporary restorations. A strip crown in the required size and shape is adjusted to the cervical area and the gingival margin with crown scissors. Before placement, ensure that the strip crown and preparation are clean and dry. Fill the strip crown with temporary filling material and place it. Strip crowns can be used as matrices

when restoring corners and edges – also when using composites. When doing so, cut out that section of the strip crown corresponding to the part of the tooth being restored, allowing for the strip crown to overlap the margins of the cavity. Fill the matrix with the definite filling material and press it into place. Remove the matrix after the filling material has been set. Should strip crowns remain on a tooth for a few days, please make sure that their gingival edge is no deeper than 1 mm inside the sulcus to avoid possible traumatisation of the marginal periodontium. Under no circumstances should strip crowns remain in the mouth for longer than 30 days.

1. Adult Strip Crowns Anterior/Posterior



Crown Forms Strip Frasaco Adult Anterior Kit-174

- Anterior assortment with 174 transparent strip crowns.
- 1 x Plastic box with transparent cover.
- 3 x Strip crowns of the same size in 58 compartments.

Jaw	Crown No.
Upper Jaw	111 – 116, 120 – 124, 131, 133, 141 – 143, 161, 162
Upper Jaw	211 – 216, 220 – 225, 231, 233, 241 – 243, 261, 262
Lower Jaw	311, 312, 321, 322, 330, 331, 341, 342, 360 – 362
Lower Jaw	411, 412, 421, 422, 430, 431, 441, 442, 460 – 462

Code: FR-STFK174
Kit-174

\$749.⁹⁵



Crown Forms Strip Frasaco Adult Anterior And Posterior Kit-58

- Compact assortment (anterior and posterior) 58 transparent strip crowns.
- 1 x Plastic box with transparent cover.
- 1 x Strip crown of one size in 58 compartments.

Jaw	Crown No.
Upper Jaw	110 – 115, 120 – 125, 130 – 135
Upper Jaw	210 – 215, 220 – 225, 230 – 235
Lower Jaw	411 – 413, 420 – 422, 430 – 434
Lower Jaw	311 – 313, 320 – 322, 330 – 334

Code: FR-STFBK58
Kit-58

\$349.⁹⁵



Crown Forms Strip Frasaco Adult Refill

- Contains 5 x transparent strip crowns of the same size, refill package.

Code: FR-ST5K
Vial-5

\$28.⁹⁵



Explanation of the description of Frasaco strip crowns:

1st digit = quadrant | 2nd digit = position of the tooth | 3rd digit = size backward-sorted



2. Pedodontal Strip Crowns

Placing prefabricated frasco pedodontal strip crowns is an excellent technique for retaining badly damaged deciduous teeth. Formerly, deciduous teeth with carious defects were extracted once these exceeded a certain stage. Nowadays, they are left in the mouth as a spacer as it is known that the deciduous teeth are important for the unimpeded development of children's jaws.

The progress made in conservative dentistry has changed the way in which this field is approached. Reference should also be made to aesthetic and phonetic aspects. We have kept pace with these changes by supplying our pedodontal strip crowns in several sizes. Allergies, which can be caused by other crown materials, are virtually eliminated.



20%
OFF!



Crown Forms Strip Frasco Pedodontal Kit - 84


- Pedodontal Strip Crown Assortment with 84 x Transparent strip crowns for upper anteriors.
- 1 x Plastic box with transparent cover.
- 3 x Strip crowns of the same size in 28 compartments.

Jaw	Crown No.
Upper Jaw	511 – 514
Upper Jaw	521 – 524
Upper Jaw	531 – 534
Upper Jaw	611 – 614
Upper Jaw	621 – 624
Upper Jaw	631 – 634

Code: FR-STKK84
Kit-84

\$399.⁹⁵

BUY 3
GET 1
BONUS!



Crown Forms Strip Frasco Pedodontal Refill

- 1 x Plastic box with transparent cover.
- Contains 5 x transparent strip crowns of the same size, refill package.

Code: FR-STK5K
Vial-5

\$29.⁹⁵