



## **Raised Toilet Seat**



Thank you for purchasing the BIOS Diagnostics™ Raised Toilet Seat, please read all instructions carefully before installing your raised toilet seat.

## **CAUTION**

- Ensure that the seat is correctly and securely locked in place before using.
- The adjustment knob, located at the front of the seat, must be tight at all times to secure the seat to your toilet.
- This product is designed for use with most common types of toilets. It will not fit all toilet types. If you have a style of toilet that this product will not fit properly or cannot be tightened securely - discontinue use.
- The padded arms are for the purpose of providing assistance while getting on or off the raised toilet seat – DO NOT use the padded arms to support full body weight during transfer, as this may result in injury.

## **Installation**

- Before each use adjust the padded arms to your specific needs. When the desired position is reached, lock the arms in place by aligning buttons, until they 'click' into their holes.
- Check that the arms are locked in place by gently pulling on them.
- Lift your toilet lid and seat assembly and place the rear section of the raised toilet seat

inside the inner edge of the toilet bowl (some toilets may require that the seat assembly be removed to accommodate a proper fit).

- Ensure that the front adjustment bracket is sitting inside the toilet bowl clearing the toilet bowl rim.
- With the raised toilet seat resting correctly on the toilet rim, straighten and center the product over the toilet bowl.
- Turn the adjusting knob clockwise to tighten the raised toilet seat securely to the toilet bowl.

Hand tighten only – use no tools and be careful not to over-tighten to prevent stripping of the adjustment bracket.

- Your weight should be centered over the toilet seat, not too far forward or too far to one side, as unit may tip.
- Ensure that the raised toilet seat is correctly installed, fits properly and is stable before using.

## **Lifetime Warranty:**

This product is warranted to be free from defects in material or workmanship for the life of the original owner. It will be repaired or replaced without charge upon receipt of the unit prepaid with \$5.00 to cover handling, packaging and return postage to the address below. Please include your full name, address and daytime phone number or email address.

Toll Free: 800-387-8520



16975 Leslie Street,  
Newmarket, ON L3Y 9A1  
[www.biosexactly.com](http://www.biosexactly.com)