



PATIENT INSTRUCTIONS: OPS COLD/COMPRESSION WRAPS

Cold & Compression in a self contained wrap. Place the wrap in the freezer to cool, then apply.

It's that simple! A handy pump is included with each wrap to add compression.

Just press a button to release excess pressure. Turn a valve and squeeze the pump to add pressure.

If you have any difficulty with the use of this product, please call a J.I.M. representative at 262-547-4276 or 866-JIM-4CPM. A J.I.M. Medical Representative is on call 24 hours a day, 7 days a week.