

Bledsoe Philippon

Post-op Hip Brace



**Bledsoe
Brace
Systems™**

The Bledsoe Philippon Post-op Hip Brace was developed in collaboration with Dr. Marc Philippon, MD, Director of Sports Medicine/ Hip Disorders at the University of Pittsburgh.



- **The Philippon Post-op Hip Brace:** is designed to allow controlled hip flexion in 15° increments from -30° to 105° following hip arthroscopy, or it can be locked in any position from -30° to 60°.



- **This easy to fit brace can be applied in less than 3 minutes:** saving valuable time in the operating room.

- **It features easily formable aluminum shells:** covered with a double-thick ultra-breathable polyurethane foam padding with a comfortable loop pile surface.

- **Three simple straps:** hold it securely affixed to the patient's hip and leg.

- **It is available:** in left and right in both regular and short sizes.

- **The regular size:** will fit most patients 5'4" (163 cm) and above.

- **The short size:** will fit patients less than 5'4" (163 cm) tall.

- **The brace circumference:** is universal.



New lightweight and easy to use Post-op Hip Brace

Bledsoe Philippon Post-op Hip Brace

Restrictions on use

This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient or for extended use beyond 180 days. This device is generally not intended for patients in excess of 250 lbs. (114 kg). This device is not intended to prevent injuries or reduce or eliminate the risk of injury.

Warning

Do not leave this device in the trunk of a car on a hot day or in any place where the temperatures may exceed 150°F (65°C) for any extended period of time. If any additional pain or symptoms occur while using this device, seek medical attention.



**Bledsoe
Brace
Systems**

by Medical Technology, Inc.

2601 Pinewood Drive
Grand Prairie, TX 75051 USA

1-888-BLED SOE
(1-888-253-3763)

1-800-527-3666

972-647-0884

Local Fax:

972-660-5495

International Fax:

972-606-0649



www.bledsoebrace.com

Indications

This device is indicated for post-operative or post-injury use for limited range-of-motion control of hip flexion/extension; as well as limited abduction and adduction of the hip joint following total hip revisions, arthroscopic hip repairs or other hip joint surgeries, injuries or problems that can benefit from range-of-motion control.

Contraindications

Contraindicated for controlling fractures of the femur or pelvis.



Ordering Information

Philippon Post-op Hip Brace

HCFA Code: L1686

Size	Leg	Part Number		
		Brace	Brace w/Positioning Kit	Brace w/Positioning Kit and Post Pad
Regular 5'4" (163 cm) & Taller	Left	AT090105	AT094105	AT093105
	Right	AT090205	AT094205	AT093205
Short Shorter than 5'4" (163 cm)	Left	AT090115		
	Right	AT090215		

Accessories

- KAF Positioning Kit
- Positioning with Post Pad Kit
- Perineal Post Pad Kit
- Strap Sets

Call for more information

Bledsoe limited warranty

This device is warranted to be free from defects in material and workmanship for a period of 90 days for wraps, pads, straps and normal wear components. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.



For more detailed information on this brace including fitting and adjustment instructions, visit our Web site: www.bledsoebrace.com

For product information or questions pertaining to sales or service, please contact the sales representative in your area or Bledsoe Brace Systems.



CP010208 Rev C 5/04