

THE INNOVATIVE TELESCOPING ELBOW BRACE

The Easy Brace

The Telescoping Elbow Brace from Bledsoe makes everything easier

- | | |
|------------------|---------------------------------|
| Easy to apply | No tools – not even scissors! |
| Easy to fit | Push-button telescoping |
| Easy to position | -10° to 110° range of motion |
| Easy to lock | Locks every 10° |
| Easy to form | Aluminum plates form to the arm |



Less Inventory

*Length adjusts to fit any patient,
Simply designate right or left*



Versatile Locking

*Locks the elbow in any position,
-10° to 100°*

*No other elbow brace offers
this wide a range*



Better Fit

*Aluminum frame easily molds
to your patient for a better, more
comfortable fit*



Optional Hand Position Control

*Maintain neutral, supination or
pronation positions of the hand
and wrist with this optional
attachment*




**Engineered for Simplicity.
Engineered by Bledsoe.**

- Dislocations
- Arthroscopic Procedures
- Stable Fracture
- Bicep Tendon Repair
- Tommy John Surgery
- Nerve Entrapment
- Internal Fixation
- Radial Head FX
- FX Olecranon

Better fit. Fewer SKUs. Easier to apply. Greater comfort.

With the Bledsoe Telescoping Elbow Brace, the joint heals in a protected range of motion. Rehab time and pain are reduced!
Protects the elbow AND promotes healing.



Arm plates can be shaped by hand

Quick push button length adjustment
- No tools required

Replicates anatomical motion of the arm

Ski-buckles to remove brace easily

Easily set wide range of motion from -10° to 110°
- Lateral hinge means no rubbing against body or clothing

Surged Pad Edges
- No pad trimming necessary

Soft, comfortable strapping

Optional Hand Attachments adjusts for neutral, supination or pronation of hand

Telescoping Elbow Brace ordering information:

PART NUMBERS

AE023100 Left One-Size-Fits-All
AE023200 Right One-Size-Fits-All

OPTIONS

Supination/Pronation Attachment

AE012123B-- Left AE012223B-- Right

Neutral Hand Position Attachment

AE012111B-- Small, Left AE012211B-- Small, Right
AE012113B-- Large, Left AE012213B-- Large, Right

INDICATIONS

Indicated for immediate application following injury or operative procedures to the elbow, or for stable fractures of the distal humerus or proximal to middle radius or ulna.



7" - 9.5" (18 - 24 cm) Small 8.5" - 12.5" (21.5 - 32 cm) Large

Sizing is for Neutral Hand Position Attachment only. Measure from ulnar styloid to medial humeral epicondyle (A to B).



1.888.BLEDSOE • LOCAL PHONE 972.647.0884 • INTERNATIONAL FAX 972.606.0649

www.BledsoeBrace.com

