


# Breast Pump Assembly Diagram



Wash and sanitize all parts before use.

\*Do not sterilize tubing.



*motif*  
LUNA

## Double Pumping



Manufactured For: Motif Medical  
8 The Green, Suite #5601 Dover, DE 19901  
Made in Korea



**Caution:** Please refer to the Motif Luna® Operation Manual for a complete overview of this product. You can download the operation manual at [motifmedical.com](http://motifmedical.com).

## Intended Use

The Motif Luna® Double Electric Breast Pump is intended for use by lactating women to express and collect breast milk. The Luna® can be used as a single or double pump and is intended for personal use only. \*Sharing a breast pump can present a health hazard.

## Product Configuration

The Motif Luna® Double Electric Breast Pump utilizes a vacuum to express breast milk. The motor is classified as a medium vacuum with an intermittent suction device.

## Safety Information

The following symbols are important for user safety. Failure to follow unit instructions can lead to personal injury or damage to the breast pump. When used in conjunction with the following words, the symbols indicate:



**Warning:** A hazard of serious injury or death.



**Caution:** A hazard of minor injury or damage to the system.

## Need More Help?

Visit us at [motifmedical.com](http://motifmedical.com) to chat with one of our friendly customer service representatives.

844-272-8390

[info@motifmedical.com](mailto:info@motifmedical.com)

