



## Discover **chromotek**<sup>®</sup> Nano-Traps: affinity beads for immunoprecipitation


Nano-Traps are the benchmark in immunoprecipitation and allow fast and reliable one-step pulldowns of even low expression proteins. They consist of Nanobodies coupled to beads and are ready to use.



**Order your free Sample**

## chromotek<sup>®</sup> Nano-Traps

- Absence of heavy & light antibody chains results in cleaner SDS-PAGE with no extra bands
- Able to withstand stringent washing conditions to remove background proteins
- Recombinantly expressed for consistent quality
- Short incubation (30-60 minutes)
- Suitable for IP, Co-IP, MS, & CHIP/RIP



▲ Immunoprecipitation of a protein of interest by Nano-Traps (left) compared with conventional antibodies coupled to Protein A/G beads (right) analyzed by SDS-PAGE. M: Marker, B: Bound.