

# ViroVac™

Surgical Smoke Plume Evacuator



## Compact Design

Allows for easy placement on a shelf or counter and accommodates stackability.

## Quiet, Yet Powerful Operation

Assists in providing a safe operating environment without disruption.

## 35 Hour Maximum Filter Life

Delivers unsurpassed value and reduces per procedure costs.

## 3-Port Filter Design

Puts the power of hose and accessory options in your hand and eliminates adapters.

## Covered Filter Ports

Filter ports are covered for additional safety.

## Simplicity

Reduces training needs and enhances utilization.

## Technical Information

### Noise Levels\*

@ maximum power < 55 dBA

### Dimensions (H x W x D)

6 in x 11 in x 15.5 in  
15.2 cm x 27.9 cm x 39.4 cm

### Flow Rates\*

25 cfm (708 lpm) using 7/8 in (22 mm) tubing

### Weight

With filter installed: 12.0 lbs (5.0 kg)  
Without filter installed: 9.5 lbs (4.0 kg)

### Warranty

One (1) Year Full Warranty

## Product Ordering Information

### Catalog Number:

VV120 100/120 Volt Systems

VV220 220/240 Volt Systems

(includes one filter and one footswitch)

VV120 NF 100/120 Volt Systems

VV220 NF 220/240 Volt Systems

(includes one footswitch)



### Buffalo Filter, LLC

5900 Genesee Street  
Lancaster, NY 14086  
U.S.A.

Phone: 800.343.2324 (US)

Phone: +1.716.835.7000

Fax: 716.835.3414

Email: [smoke-info@buffalofilter.com](mailto:smoke-info@buffalofilter.com)

Web: [www.buffalofilter.com](http://www.buffalofilter.com)

\*Nominal values

Made in USA


PN 903209 REVC

# ViroVac™

## Surgical Smoke Plume Evacuator

ViroVac™ is designed to be simple to use, flexible in its options for tubing and accessories, and sets new precedents for quiet operation. Providing effective surgical smoke plume evacuation that encompasses 4 layers of filtration and up to 35 hours of filter life, ViroVac™ will deliver the performance necessary to effectively manage surgical smoke plume with unsurpassed value.

EZLink, an automatic activation device (sold separately), can be used to remotely turn the evacuator on and off, saving filter life.



Pneumatic footswitch (included)



EZLink (sold separately; part: EZLINK01).

Buffalo Filter® is a registered trademark of Buffalo Filter, LLC in the United States and other countries. All other trademarks used herein are the property of their respective owners. ©2014, Buffalo Filter. All rights reserved.